

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-005945

Location: TOXL

- | | Check When Done |
|---|-----------------|
| A. Is the breath tube heated? | <u>✓</u> |
| B. Is the diagnostic check complete? | <u>✓</u> |
| C. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| D. Print test. (Attach test record.) | <u>✓</u> |
| E. Low AC. Use <0.03 AC in ACA mode. (Attach test record.) Sim SN <u>DR3803</u> Lot No. <u>13070</u> Rep AC <u>0.020</u> | <u>✓</u> |
| F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.) Sim SN <u>DR5187</u> Lot No. <u>13080</u> Rep AC <u>0.300</u> | <u>✓</u> |
| G. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode. (Attach test record.) Sim SN <u>DR3375</u> Lot No. <u>1CS2</u> Rep AC <u>0.10 Act + 0.05% acetone</u> | <u>✓</u> |
| H. RFI check. Run CMS mode. Key radio on first room air. (Attach test record.) | <u>✓</u> |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode. (Attach test record.) | |

| | | |
|---------------------------|---------------------|-----------------------------|
| Lot No. <u>16512080A2</u> | Cyl No. <u>22</u> | Exp. Date <u>08/01/2014</u> |
| Test 1 <u>0.080</u> | Test 4 <u>0.080</u> | Test 7 <u>0.079</u> |
| Test 2 <u>0.080</u> | Test 5 <u>0.080</u> | Test 8 <u>0.079</u> |
| Test 3 <u>0.079</u> | Test 6 <u>0.079</u> | Test 9 <u>0.080</u> |
| Average <u>0.080</u> | | |

J. Tank pressure, Level 3 Func D Sub G P: 848 psi Reg 850 psi

K. Remarks/Maintenance record: Calibration, Annual inspection

04 Apr 2014
Date

04 APR 2014
Date

Deb Kashur
Field Inspector's Signature

Charles E. Edm
Reviewed by

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005945
 04/03/2014 15:55:57

Auto Calibration

pg 1 of 2

| | <<<<< | 3um | >>>>> | <<<<< | 9um | >>>>> |
|-------------|--------|-------------|-----------|-------------|--------------|---------------|
| ----- | | | | | | |
| Solution = | 0.000 | g/210L | or 0.0000 | mg/l, | Samples = 4, | Discarded = 1 |
| Sample | % Abs | (% Abs Ref) | % Abs | (% Abs Ref) | | |
| Sample #1 | 0.0920 | (0.0090) | 0.2210 | (0.0050) | | |
| Sample #2 | 0.0790 | (0.0380) | 0.2060 | (0.0150) | | |
| Sample #3 | 0.0920 | (0.0560) | 0.2060 | (0.0010) | | |
| Sample #4 | 0.0830 | (0.0420) | 0.2270 | (-0.0200) | | |
| Avg % Abs | 0.0847 | (0.0453) | 0.2130 | (-0.0013) | | |
| STD DEV | 0.0067 | (0.0095) | 0.0121 | (0.0176) | | |
| REL STD DEV | 7.864 | (20.849) | 5.692 | (1321.221) | | |

0.000 AC
 Milli-Q H₂O
 G10329

| | | | | | | |
|-------------|--------|-------------|-----------|-------------|--------------|---------------|
| ----- | | | | | | |
| Solution = | 0.050 | g/210L | or 0.2381 | mg/l, | Samples = 4, | Discarded = 1 |
| Sample | % Abs | (% Abs Ref) | % Abs | (% Abs Ref) | | |
| Sample #1 | 0.9530 | (-0.0160) | 1.8840 | (-0.0290) | | |
| Sample #2 | 0.9510 | (-0.0170) | 1.8880 | (-0.0280) | | |
| Sample #3 | 0.9360 | (0.0090) | 1.8670 | (-0.0070) | | |
| Sample #4 | 0.9410 | (0.0040) | 1.8950 | (-0.0250) | | |
| Avg % Abs | 0.9427 | (-0.0013) | 1.8833 | (-0.0200) | | |
| STD DEV | 0.0076 | (0.0138) | 0.0146 | (0.0114) | | |
| REL STD DEV | 0.810 | (1034.710) | 0.774 | (56.789) | | |

0.050 AC
 lot: 13170
 Exp: 06/19/15
 DR5144

| | | | | | | |
|-------------|--------|-------------|-----------|-------------|--------------|---------------|
| ----- | | | | | | |
| Solution = | 0.080 | g/210L | or 0.3810 | mg/l, | Samples = 4, | Discarded = 1 |
| Sample | % Abs | (% Abs Ref) | % Abs | (% Abs Ref) | | |
| Sample #1 | 1.5220 | (-0.0220) | 3.0290 | (-0.0060) | | |
| Sample #2 | 1.5040 | (-0.0040) | 3.0050 | (0.0070) | | |
| Sample #3 | 1.4740 | (0.0150) | 2.9550 | (0.0300) | | |
| Sample #4 | 1.4690 | (0.0230) | 2.9390 | (0.0440) | | |
| Avg % Abs | 1.4823 | (0.0113) | 2.9663 | (0.0270) | | |
| STD DEV | 0.0189 | (0.0139) | 0.0344 | (0.0187) | | |
| REL STD DEV | 1.277 | (122.368) | 1.161 | (69.191) | | |

0.080 AC
 lot: 13190
 Exp: 07/16/15
 DR5132

| | | | | | | |
|-------------|--------|-------------|-----------|-------------|--------------|---------------|
| ----- | | | | | | |
| Solution = | 0.200 | g/210L | or 0.9524 | mg/l, | Samples = 4, | Discarded = 1 |
| Sample | % Abs | (% Abs Ref) | % Abs | (% Abs Ref) | | |
| Sample #1 | 3.4710 | (0.0040) | 6.7920 | (0.0020) | | |
| Sample #2 | 3.4910 | (0.0230) | 6.7820 | (0.0340) | | |
| Sample #3 | 3.4790 | (0.0240) | 6.7910 | (0.0100) | | |
| Sample #4 | 3.4680 | (0.0330) | 6.7580 | (0.0280) | | |
| Avg % Abs | 3.4793 | (0.0267) | 6.7770 | (0.0240) | | |
| STD DEV | 0.0115 | (0.0055) | 0.0171 | (0.0125) | | |
| REL STD DEV | 0.331 | (20.653) | 0.252 | (52.042) | | |

0.200 AC
 lot: 12120
 Exp: 08/10/14
 DR5113

| | | | | | | |
|-------------|--------|-------------|-----------|-------------|--------------|---------------|
| ----- | | | | | | |
| Solution = | 0.400 | g/210L | or 1.9048 | mg/l, | Samples = 4, | Discarded = 1 |
| Sample | % Abs | (% Abs Ref) | % Abs | (% Abs Ref) | | |
| Sample #1 | 6.7160 | (-0.0150) | 13.0720 | (0.0030) | | |
| Sample #2 | 6.6750 | (0.0270) | 12.9670 | (0.0580) | | |
| Sample #3 | 6.6730 | (0.0360) | 12.9990 | (0.0590) | | |
| Sample #4 | 6.6980 | (0.0270) | 12.9730 | (0.0430) | | |
| Avg % Abs | 6.6820 | (0.0300) | 12.9797 | (0.0533) | | |
| STD DEV | 0.0139 | (0.0052) | 0.0170 | (0.0090) | | |
| REL STD DEV | 0.208 | (17.321) | 0.131 | (16.805) | | |

0.400 AC
 lot: 12105
 Exp: 05/23/14
 DR 5185

Deb Kashur
 calibration

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005945
 04/03/2014 15:55:57

Auto Calibration

pg 2 of 2

<<<<< 3um >>>>>

<<<<< 9um >>>>>

 Zero Order Coef -225.51
 First Order Coef 2706.67
 Second Order Coef 26.62

 -310.30
 1394.52
 7.50

| Act (g/210L) | Fit (g/210L) | Residual (g/210L) |
|-----------------|-----------------|----------------------|
| 0.000 | 0.000 | -0.0001 |
| 0.050 | 0.049 | 0.0007 |
| 0.080 | 0.081 | -0.0007 |
| 0.200 | 0.200 | 0.0002 |
| 0.400 | 0.400 | -0.0000 |

| Act (g/210L) | Fit (g/210L) | Residual (g/210L) |
|-----------------|-----------------|----------------------|
| 0.000 | -0.000 | 0.0003 |
| 0.050 | 0.049 | 0.0008 |
| 0.080 | 0.082 | -0.0017 |
| 0.200 | 0.199 | 0.0008 |
| 0.400 | 0.400 | -0.0001 |

<<<<< 3um >>>>>

<<<<< 9um >>>>>

Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1

| Sample | 3um | 9um |
|-----------------------|-----------|-----------|
| Sample #1 | 3126.00 | 3146.00 |
| Sample #2 | 3080.00 | 3171.00 |
| Sample #3 | 3099.00 | 3161.00 |
| Sample #4 | 3079.00 | 3189.00 |
| Avg | 3086.0000 | 3173.6667 |
| STD DEV | 11.2694 | 14.1892 |
| REL STD DEV | 0.365 | 0.447 |
| H2O adjust (mg/l*10k) | 723 | 636 |

0.080 AC
 lot: 16S12080A2
 Exp: 8/01/2014
 cyl: 22

Atmospheric Pressure = 1011

*****CALIBRATION SUCCESSFUL*****

*Deb Kashur
 Calibration*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005945
Location = TOXL 8164.13.00 06/09
04/04/2014 08:20

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 16512080A2
Standard Cyl #? 22
Standard Expiration? 08/01/2014
Oper No? 777777

Flow Cal. Date: 06/20/2013
Slope 637
Intercept -566407

IR Calibration Date: 04/03/2014
3um 9um

0th Coef(*100): -22551 -31030
1st Coef(*100): 270667 139452
2nd Coef(*100): 2662 750
H2O adj(mg/l*10k): 723 636

***** Printer Test End *****

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Print Test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005945
Location = TOXL 8164.13.00 06/09
04/04/2014 08:14

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 08:15 |
| 02 Std. Sol. | 0.019 | 08:15 |
| 03 Room Air | 0.000 | 08:16 |
| 04 Std. Sol. | 0.020 | 08:17 |
| 05 Room Air | 0.000 | 08:17 |
| 06 Std. Sol. | 0.019 | 08:18 |
| 07 Room Air | 0.000 | 08:18 |

08 Sim Temp = 34.0°C

Simul Ser No = DR3803

Std Sol No = 13070

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Low AC 0.020 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005945
Location = TOXL 8164.13.00 06/09
04/03/2014 16:48

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 16:49 |
| 02 Std. Sol. | 0.304 | 16:51 |
| 03 Room Air | 0.000 | 16:51 |
| 04 Std. Sol. | 0.303 | 16:52 |
| 05 Room Air | 0.000 | 16:52 |
| 06 Std. Sol. | 0.302 | 16:53 |
| 07 Room Air | 0.000 | 16:54 |

08 Sim Temp = 34.0°C

Simul Ser No = DR5187
Std Sol No = 13080
County = 08 Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:

High AC 0.300 AC


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005945
Location = TOXL 8164.13.00 06/09
04/04/2014 08:21

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 08:22 |
| 02 *Subject Test | INT* | 08:22 |
| 03 Room Air | 0.000 | 08:23 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 16512080A2
Cyl No = 22
Expiration Date = 08/01/2014
County = 08 Oper No. = 777777



Operator Signature
DEB KASHUR

Remarks:

Interferent test 0.100AC + 0.05% acetone

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005945
Location = TOXL 8164.13.00 06/09
04/04/2014 08:30

| Test | AC | Time |
|---------------|-------|-------|
| 01 Diagnostic | OK | 08:30 |
| 02 Room Air | RFI* | 08:31 |
| 03 Room Air | 0.000 | 08:31 |

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 16512080A2
Cyl No = 22
Expiration Date = 08/01/2014
County = 08 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
DEB KASHUR

Remarks:

RFI test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005945
Location = TOXL 8164.13.00 06/09
04/04/2014 08:33

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 08:34 |
| 02 Std. Gas | 0.080 | 08:34 |
| 03 Room Air | 0.000 | 08:34 |
| 04 Std. Gas | 0.080 | 08:35 |
| 05 Room Air | 0.000 | 08:35 |
| 06 Std. Gas | 0.079 | 08:36 |
| 07 Room Air | 0.000 | 08:36 |

Lot No = 16512080A2
Cyl No = 22
Exp Date = 08/01/2014
County = 08 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005945
Location = TOXL 8164.13.00 06/09
04/04/2014 08:42

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 08:42 |
| 02 Std. Gas | 0.080 | 08:43 |
| 03 Room Air | 0.000 | 08:43 |
| 04 Std. Gas | 0.080 | 08:43 |
| 05 Room Air | 0.000 | 08:44 |
| 06 Std. Gas | 0.079 | 08:44 |
| 07 Room Air | 0.000 | 08:45 |

Lot No = 16512080A2
Cyl No = 22
Exp Date = 08/01/2014
County = 08 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005945
Location = TOXL 8164.13.00 06/09
04/04/2014 08:46

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 08:46 |
| 02 Std. Gas | 0.079 | 08:47 |
| 03 Room Air | 0.000 | 08:47 |
| 04 Std. Gas | 0.079 | 08:48 |
| 05 Room Air | 0.000 | 08:48 |
| 06 Std. Gas | 0.080 | 08:49 |
| 07 Room Air | 0.000 | 08:49 |

Lot No = 16512080A2

Cyl No = 22

Exp Date = 08/01/2014

County = 08

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080AC

Form 106-I8000