

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004193

Location: TOXL

- | | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5132</u> Lot No. <u>12210</u> Rep AC <u>0.020</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR3841</u> Lot No. <u>13620</u> Rep AC <u>1/15/15</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR3846</u> Lot No. <u>—</u> Rep AC <u>0.100 AC + acetone</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>23411080A1</u>	Cyl No. <u>15</u>	Exp. Date <u>10/01/2013</u>
Test 1 <u>0.080</u>	Test 4 <u>0.079</u>	Test 7 <u>0.080</u>
Test 2 <u>0.080</u>	Test 5 <u>0.080</u>	Test 8 <u>0.079</u>
Test 3 <u>0.079</u>	Test 6 <u>0.079</u>	Test 9 <u>0.080</u>
Average <u>0.080</u>		

L. Tank pressure, Level 3 Func D Sub G P: 577 psi Reg 600 psi

M. Remarks/Maintenance record: Annual inspection, replaced external battery, recalibrated optical bench. ^(3/18/13)
08 Apr 2013 Date Deb Kashner Field Inspector's Signature
08 APR 13 Date Charles E. Edr Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = TOXL 8164.13.00 06/09
04/08/2013 09:24

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:25
02 Std. Sol.	0.021	09:26
03 Room Air	0.000	09:26
04 Std. Sol.	0.021	09:27
05 Room Air	0.000	09:28
06 Std. Sol.	0.021	09:28
07 Room Air	0.000	09:29

08 Sim Temp = 34.0°C

Simul Ser No = DR5132

Std Sol No = 12210

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Low AC ^{DK} ~~0.021~~ 0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = TOXL 8164.13.00 06/09
04/08/2013 09:30

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:31
02 Std. Sol.	0.299	09:32
03 Room Air	0.000	09:32
04 Std. Sol.	0.302	09:33
05 Room Air	0.000	09:33
06 Std. Sol.	0.303	09:34
07 Room Air	0.000	09:35

08 Sim Temp = 34.0°C

Simul Ser No = DR3841

Std Sol No = 13020

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Linearity test 0.300 AC

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = TOXL 8164.13.00 06/09
04/08/2013 09:36

Test	AC	Time
01 Room Air	0.000	09:37
02 *Subject Test	INT*	09:37
03 Room Air	0.000	09:38

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 23411080A1
Cyl No = 15
Expiration Date = 10/01/2013
County = 08 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Interferent check 0.180AC + acetone

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = TOXL 8164.13.00 06/09
04/08/2013 09:38

Test	AC	Time
01 Diagnostic	OK	09:38
02 Room Air	RFI*	09:39
03 Room Air	0.000	09:39

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 23411080A1
Cyl No = 15
Expiration Date = 10/01/2013
County = 08 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

RFI check

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = TOXL 8164.13.00 06/09
04/08/2013 09:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:40
02 Std. Gas	0.080	09:40
03 Room Air	0.000	09:41
04 Std. Gas	0.080	09:41
05 Room Air	0.000	09:42
06 Std. Gas	0.079	09:42
07 Room Air	0.000	09:42

Lot No = 23411080A1
Cyl No = 15
Exp Date = 10/01/2013
County = 08 Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = TOXL 8164.13.00 06/09
04/08/2013 09:43

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:43
02 Std. Gas	0.079	09:44
03 Room Air	0.000	09:44
04 Std. Gas	0.080	09:44
05 Room Air	0.000	09:45
06 Std. Gas	0.079	09:45
07 Room Air	0.000	09:46

Lot No = 23411080A1
Cyl No = 15
Exp Date = 10/01/2013
County = 08 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = TOXL 8164.13.00 06/09
04/08/2013 09:46

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:47
02 Std. Gas	0.080	09:47
03 Room Air	0.000	09:48
04 Std. Gas	0.079	09:48
05 Room Air	0.000	09:49
06 Std. Gas	0.080	09:49
07 Room Air	0.000	09:49

Lot No = 23411080A1
Cyl No = 15
Exp Date = 10/01/2013
County = 08 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080AC

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-004193
 04/08/2013 08:38:52

Auto Calibration

pg 1 of 2

<<<<<			3um	>>>>>			<<<<<			9um	>>>>>		
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1													
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		
Sample #1	0.0110	(0.0040)		0.0250	(0.0000)		0.0250	(0.0000)		0.0250	(0.0000)		
Sample #2	0.0430	(0.0100)		0.0650	(-0.0110)		0.0650	(-0.0110)		0.0650	(-0.0110)		
Sample #3	0.0340	(0.0420)		0.0460	(0.0200)		0.0460	(0.0200)		0.0460	(0.0200)		
Sample #4	0.0430	(0.0490)		0.0370	(0.0210)		0.0370	(0.0210)		0.0370	(0.0210)		
Avg % Abs	0.0400	(0.0337)		0.0493	(0.0100)		0.0493	(0.0100)		0.0493	(0.0100)		
STD DEV	0.0052	(0.0208)		0.0143	(0.0182)		0.0143	(0.0182)		0.0143	(0.0182)		
REL STD DEV	12.990	(61.760)		28.975	(181.934)		28.975	(181.934)		28.975	(181.934)		
Solution = 0.020 g/210L or 0.0952 mg/l, Samples = 4, Discarded = 1													
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		
Sample #1	0.4310	(0.0070)		0.8760	(-0.0120)		0.8760	(-0.0120)		0.8760	(-0.0120)		
Sample #2	0.4510	(0.0120)		0.8950	(-0.0180)		0.8950	(-0.0180)		0.8950	(-0.0180)		
Sample #3	0.4230	(0.0380)		0.8830	(-0.0030)		0.8830	(-0.0030)		0.8830	(-0.0030)		
Sample #4	0.4230	(0.0380)		0.8830	(0.0040)		0.8830	(0.0040)		0.8830	(0.0040)		
Avg % Abs	0.4323	(0.0293)		0.8870	(-0.0057)		0.8870	(-0.0057)		0.8870	(-0.0057)		
STD DEV	0.0162	(0.0150)		0.0069	(0.0112)		0.0069	(0.0112)		0.0069	(0.0112)		
REL STD DEV	3.739	(51.174)		0.781	(198.350)		0.781	(198.350)		0.781	(198.350)		
Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1													
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		
Sample #1	0.7760	(0.0080)		1.5800	(0.0000)		1.5800	(0.0000)		1.5800	(0.0000)		
Sample #2	0.7890	(0.0000)		1.5530	(0.0210)		1.5530	(0.0210)		1.5530	(0.0210)		
Sample #3	0.8000	(0.0150)		1.5930	(0.0020)		1.5930	(0.0020)		1.5930	(0.0020)		
Sample #4	0.7920	(0.0150)		1.5790	(0.0010)		1.5790	(0.0010)		1.5790	(0.0010)		
Avg % Abs	0.7937	(0.0100)		1.5750	(0.0080)		1.5750	(0.0080)		1.5750	(0.0080)		
STD DEV	0.0057	(0.0087)		0.0203	(0.0113)		0.0203	(0.0113)		0.0203	(0.0113)		
REL STD DEV	0.716	(86.603)		1.289	(140.868)		1.289	(140.868)		1.289	(140.868)		
Solution = 0.200 g/210L or 0.9524 mg/l, Samples = 4, Discarded = 1													
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		
Sample #1	3.4970	(0.0100)		6.8730	(-0.0020)		6.8730	(-0.0020)		6.8730	(-0.0020)		
Sample #2	3.5270	(0.0160)		6.8900	(0.0280)		6.8900	(0.0280)		6.8900	(0.0280)		
Sample #3	3.5130	(0.0330)		6.9220	(0.0210)		6.9220	(0.0210)		6.9220	(0.0210)		
Sample #4	3.5260	(0.0350)		6.9390	(0.0350)		6.9390	(0.0350)		6.9390	(0.0350)		
Avg % Abs	3.5220	(0.0280)		6.9170	(0.0280)		6.9170	(0.0280)		6.9170	(0.0280)		
STD DEV	0.0078	(0.0104)		0.0249	(0.0070)		0.0249	(0.0070)		0.0249	(0.0070)		
REL STD DEV	0.222	(37.287)		0.360	(25.000)		0.360	(25.000)		0.360	(25.000)		
Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1													
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		
Sample #1	5.1760	(0.0020)		10.0260	(-0.0020)		10.0260	(-0.0020)		10.0260	(-0.0020)		
Sample #2	5.2110	(0.0200)		10.0810	(0.0210)		10.0810	(0.0210)		10.0810	(0.0210)		
Sample #3	5.1750	(0.0330)		10.0650	(0.0370)		10.0650	(0.0370)		10.0650	(0.0370)		
Sample #4	5.1780	(0.0300)		10.0710	(0.0330)		10.0710	(0.0330)		10.0710	(0.0330)		
Avg % Abs	5.1880	(0.0277)		10.0723	(0.0303)		10.0723	(0.0303)		10.0723	(0.0303)		
STD DEV	0.0200	(0.0068)		0.0081	(0.0083)		0.0081	(0.0083)		0.0081	(0.0083)		
REL STD DEV	0.385	(24.603)		0.080	(27.451)		0.080	(27.451)		0.080	(27.451)		

G10329
0.000 AC

DR5132
Lot: 12210
Exp: 11/2/14
0.020 AC

DR5113
Lot: 11150
Exp: 8/9/13

DR5144
Lot: 11120
Exp: 7/19/13

DR3803
Lot: 12050
Exp: 3/21/14

Calibration
Deb Kashur

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-004193
 04/08/2013 08:38:52

Auto Calibration

pg 2 of 2

```

<<<<<      3um      >>>>>
-----
Zero Order Coef   -161.95
First Order Coef  2638.40
Second Order Coef 28.79
  
```

```

<<<<<      9um      >>>>>
-----
                -139.56
                1296.13
                13.69
  
```

```

-----
Act      Fit      Residual
(g/210L) (g/210L)   (g/210L)
0.000    -0.001    0.0012
0.020     0.021   -0.0007
0.040     0.041   -0.0010
0.200     0.199    0.0008
0.300     0.300   -0.0003
-----
  
```

```

-----
Act      Fit      Residual
(g/210L) (g/210L)   (g/210L)
0.000    -0.002    0.0016
0.020     0.021   -0.0014
0.040     0.041   -0.0007
0.200     0.199    0.0009
0.300     0.300   -0.0004
-----
  
```

```

<<<<<      3um      >>>>>
-----
  
```

```

<<<<<      9um      >>>>>
-----
  
```

Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1

Sample

Sample #1	3336.00	3303.00
Sample #2	3306.00	3311.00
Sample #3	3238.00	3291.00
Sample #4	3268.00	3269.00
Avg	3270.6667	3290.3333
STD DEV	34.0783	21.0079
REL STD DEV	1.042	0.638
H2O adjust (mg/l*10k)	539	519

Atmospheric Pressure = 1012

*****CALIBRATION SUCCESSFUL*****

*Calibration
 Deb Kashner*