

INTOXILYZER 8000 INSPECTION

Serial No.: 80-003069

Location: TOX LAB

- |   | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes?  | <u>✓</u>        |
| B. Is a three-pronged grounded outlet used?   | <u>✓</u>        |
| C. Is the breath tube heated?   | <u>✓</u>        |
| D. Is the diagnostic check complete?  | <u>✓</u>        |
| E. Is the time, date, and year correct (re-set if necessary)?   | <u>✓</u>        |
| F. Print test. (Attach test record.)  | <u>✓</u>        |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5190</u> Lot No. <u>—</u> Rep AC <u>0.019 AC</u>  | <u>✓</u>        |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5113</u> Lot No. <u>—</u> Rep AC <u>0.285 AC</u>   | <u>✓</u>        |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR5135</u> Lot No. <u>—</u> Rep AC <u>0.10 AC + 0.05% Acetone</u> | <u>✓</u>        |
| J. RFI check. Run CMS mode. Key radio on first room air.<br>(Attach test record.)   | <u>✓</u>        |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.)   |                 |

Lot No. <u>659358</u>	Cyl No. <u>29</u>	Exp. Date <u>09/03/2011</u>
Test 1 <u>0.081</u>	Test 4 <u>0.080</u>	Test 7 <u>0.081</u>
Test 2 <u>0.081</u>	Test 5 <u>0.081</u>	Test 8 <u>0.081</u>
Test 3 <u>0.081</u>	Test 6 <u>0.081</u>	Test 9 <u>0.081</u>
Average <u>0.081</u>		

L. Tank pressure, Level 3 Func D Sub G P: 350 psi Reg 360 psi

M. Remarks/Maintenance record: Annual Inspection

<u>16 April 2010</u> Date	<u>Mary Pearson</u> Field Inspector's Signature
<u>19 Apr 2010</u> Date	<u>Deb Kashner</u> Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-003069  
Location = TOXL                    8164.13.00 06/09  
04/16/2010                            11:32

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Wet  
Standard Value?            0.000  
Standard Solution #?        00000  
Sim Serial #?              000000  
Oper No?                    107501

Flow Cal. Date:            04/08/2008  
Slope                      691  
Intercept                  -599264

IR Calibration Date:        06/04/2008  
                                    3um                      9um

-----  
0th Coef(\*100):            -14879                    -13696  
1st Coef(\*100):            280678                    136270  
2nd Coef(\*100):            177                        615  
H2O adj(mg/l\*10k):        391                        185

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
-----  
Operator Signature  
MARGY PEARSON

Remarks:

*Paint*

Form 106-I8000

80-003069  
6/8/05 KM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003069  
Location = TOXL      8164.13.00 06/09  
04/16/2010      11:44

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:45
02 Std. Sol.	0.018	11:46
03 Room Air	0.000	11:46
04 Std. Sol.	0.019	11:47
05 Room Air	0.000	11:48
06 Std. Sol.	0.019	11:48
07 Room Air	0.000	11:49


08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 00000

County = 08

Oper No. = 107501



Operator Signature  
MARGY PEARSON

Remarks:

*Low AC 0.019*

Form 106-I8000

*80-003069  
6/15 PM*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003069  
Location = TOXL      8164.13.00 06/09  
04/16/2010      11:50

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:51
02 Std. Sol.	0.283	11:52
03 Room Air	0.000	11:53
04 Std. Sol.	0.284	11:53
05 Room Air	0.000	11:54
06 Std. Sol.	0.284	11:55
07 Room Air	0.000	11:55


08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 00000

County = 08

Oper No. = 107501



Operator Signature  
MARGY PEARSON

Remarks:

Linearity 0.285 AC

Form 106-I8000

80-003069  
7/15 FN

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003069  
Location = TOXL      8164.13.00 06/09  
04/16/2010      11:25

Test	AC	Time
01 Diagnostic	OK	11:26
02 Room Air	0.000	11:26
03 *Subject Test	INT*	11:29
04 Room Air	0.000	11:29

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS821456  
Simul Ser No = 000000  
Std Sol No = 00000  
County = 08      Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

0.10% AC + Acetone 0.05%

Form 106-I8000

80-003069  
7/75 KM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-003069  
Location = TOXL                    8164.13.00 06/09  
04/16/2010                            11:30

Test	AC	Time
01 Diagnostic	OK	11:31
02 Room Air	RFI*	11:31
03 Room Air	0.000	11:32

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS821456  
Simul Ser No = 000000  
Std Sol No = 00000  
County = 08                            Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

RFI

Form 106-I8000

80-003069  
7/25 PM


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003069  
Location = TOXL      8164.13.00 06/09  
04/16/2010      11:58

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:58
02 Std. Gas	0.081	11:59
03 Room Air	0.000	11:59
04 Std. Gas	0.081	12:00
05 Room Air	0.000	12:00
06 Std. Gas	0.081	12:00
07 Room Air	0.000	12:01

Lot No = 659358  
Cyl No = 29  
Exp Date = 09/03/2011  
County = 08      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

*Cal Check*

Form 106-I8000

*80-003069  
7/25 FH*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003069  
Location = TOXL      8164.13.00 06/09  
04/16/2010      12:02

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:02
02 Std. Gas	0.080	12:02
03 Room Air	0.000	12:03
04 Std. Gas	0.081	12:03
05 Room Air	0.000	12:04
06 Std. Gas	0.081	12:04
07 Room Air	0.000	12:05

Lot No = 659358  
Cyl No = 29  
Exp Date = 09/03/2011  
County = 08      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

*Cal Check*

Form 106-I8000

*80-003069  
7/15/10*




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003069  
Location = TOXL      8164.13.00 06/09  
04/16/2010      12:05

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:06
02 Std. Gas	0.081	12:06
03 Room Air	0.000	12:06
04 Std. Gas	0.081	12:07
05 Room Air	0.000	12:07
06 Std. Gas	0.081	12:08
07 Room Air	0.000	12:08

Lot No = 659358  
Cyl No = 29  
Exp Date = 09/03/2011  
County = 08      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

*Cal Check*

Form 106-18000

80-003069  
75/75-FAH