



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Tray N Nielsen</i>		
Location <i>PD3 1230 255250</i>	Intoxilyzer® Serial Number <i>80-006495</i>	
Gas Lot Number <i>1981708 OAH</i>	Gas Cylinder Number <i>14</i>	Gas Expiration Date <i>9-5-19</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. ~~0.00~~ AC *TM* 0.080
 - B. 0.080 AC
 - C. 0.080 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i>	Date <i>6-11-18</i>
Reviewed By (Crime Laboratory Use Only) <i>Rebecca Huegel-Morris</i>	Date <i>11 Jun 18</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.14.00 09/16
06/11/2018 16:07

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:08
02 Std. Gas	0.080	16:08
03 Room Air	0.000	16:09
04 Std. Gas	0.080	16:09
05 Room Air	0.000	16:10
06 Std. Gas	0.080	16:10
07 Room Air	0.000	16:11

Lot No = 19817080A4
Cyl No = 14
Exp Date = 09/05/2019
County = 09 Oper No. = 110104



Operator Signature
TROY NIELSEN

Remarks: