



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Pamela B. Pehl</i>		
Location <i>Stanton, ND</i>	Intoxilyzer® Serial Number <i>80-006685</i>	
Gas Lot Number <i>17618080 A3</i>	Gas Cylinder Number <i>069</i>	Gas Expiration Date <i>8/5/2020</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

- A. 0.080 AC
- B. 0.081 AC
- C. 0.080 AC

- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>P. Pehl</i>	Date <i>2/6/2019</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Yungu-nimmo</i>	Date <i>2/7/2019</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006685
Location = STAT 8164.14.00 09/16
02/06/2019 17:00

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:01
02 Std. Gas	0.080	17:01
03 Room Air	0.000	17:02
04 Std. Gas	0.081	17:02
05 Room Air	0.000	17:03
06 Std. Gas	0.080	17:03
07 Room Air	0.000	17:04

Lot No = 17618080A3
Cyl No = 69
Exp Date = 08/05/2020
County = 29 Oper No. = 120906



Operator Signature
PAMELA PEHL

Remarks:

Form 106-I8000