

**INTOXILYZER® 8000 INITIAL INSPECTION**

Intoxilyzer® 8000 Serial Number: 80-006490

**A. GENERAL Setup and Checks:**

1.  Power On
  - 5 min. countdown during initial warm-up
  - 1 min. countdown to begin test after standby mode
2.  Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:
  - Level 1 (Operator)    Level 2 (Field Inspector)    Level 3 (Lab Personnel)
4.  Printer Setup (Level 3,S,P)
5.  Forms Load (CMI, Inc., CoBrA Software)

**B. LEVEL 1 Setup and Checks:**

1.  Internal Printer Test and Paper Feed (Level 1,P). Attach record.
2.  External Printer Test (Level 1,P). Attach record.
3. Barcode Reader (Level 1,B):
  - Operator ID    ND License    R key: Refusal. Attach record.

**C. LEVEL 2 Setup and Checks:**

1.  Continuous Air Blank (Level 2,A)
2.  Diagnostic (Level 2,D). Attach record.
3.  Prelim Data Entry (Level 2,E)
4. Ser Number/Version Info (Level 2,V):
  - Inst Ser#: 80-006490   DSP Ver: 1362.14
  - Core Ver: 8164.14   CPLD Ver: 7

**D. LEVEL 3 Setup and Checks:**

1. DVM Monitor (Level 3,D,D): 3: 12506   s: 1   c: 47.06  
9: 13281   s: 2   B: 45 0  
(RFI indicator with keyed radio):
2. Function (F) Keys: :  F5 – Tone                       F9 - Fan  
 F6 – Green LED                     F10 – Pump  
 F7 – Red LED                         F11 – Solenoid
3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 958   CF: 1.05
4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 193   F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T): C: 47.0 A: 47.0 I: 0.0  
B: 45.0 V: 0.0
6. Tank Monitor (Level 3,D,G): P (display): 328 psi P (gauge): 350 psi
7.  General Setup (Level 3,S,G)
8.  Configure Start Test (Level 3,S,T)
9.  Network (Level 3,S,N)

**E. TESTS (attach all test records):**

**Configure wet bath simulator standards before running tests (Level 1,S):**

1. Zero AC (Level 1,C):  
Sim ser #: MP3002 Lot #: MILLI-Q H<sub>2</sub>O Exp date: NA  
Known value: 0.000 AC  
Results: 0.000 AC 0.000 AC 0.000 AC
2. Low AC (Level 1,C):  
Sim ser #: MP3063 Lot #: 20130 Exp date: 3.10.22  
Known value: 0.030 AC  
Results: 0.029 AC 0.029 AC 0.029 AC
3. High AC (Level 1,C):  
Sim ser #: MP3069 Lot #: 201902D Exp date: 2.20.21  
Known value: 0.300 AC  
Results: 0.295 AC 0.295 AC 0.295 AC

**Configure dry gas standard before running tests (Level 1,S):**

4. Ethanol + 0.05% Acetone (Level 1,B):  
Sim ser #: DR3451 Lot #: ICS6 Exp date: NA  
Known value: 0.10 AC + 0.05% ACETONE  
 Display reads "Interferent Detect"
5. Acetone (Level 1,B):  
Sim ser #: DR5188 Lot #: ACE1 Exp date: NA  
Known value: 0.05% ACETONE  
 Display reads "Interferent Detect"
6. Methanol (Level 1,B):  
Sim ser #: DR3455 Lot #: MET2 Exp date: NA  
Known value: 0.10% METHANOL  
 Display reads "Interferent Detect"
7. Isopropanol (Level 1,B):  
Sim ser #: DR5189 Lot #: IPA2 Exp date: NA  
Known value: 0.10% ISOPROPNOL  
 Display reads "Interferent Detect"

- 8.  Ambient Fail (Level 1,B): Display and test record read correct.
- 9.  Invalid Sample (Level 1,B): Display and test record read correct.
- 10.  Sequence Aborted (Level 1,B): Display and test record read correct.
- 11.  ABA (Level 1,B): Test record reads correct.
- 12.  Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- 13.  RFI Detected (B,C, or CMS mode): Display and test record read correct.
- 14.  Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 15.  Deficient Samples (CMS mode): Display and test record read correct.
- 16.  Improper Sample (CMS mode): Display and test record read correct.
- 17.  Difference Too Great (CMS mode): Display and test record read correct.
- 18.  CMS (CMS mode): Test record reads correct.
- 19. Stability Test (Level 3,T,S):

Ave: 0.0811 Std Dev: 0.0004 Rel Std Dev: 0.56%

Maintenance/Repairs/Notes: INITIAL INSPECTION DUE TO OPENING UP  
INSTRUMENT DURING REPAIR

---

---

---

---

---

---

Instrument is acceptable to be used in the field.

Charles E. Ehr  
Field Inspector's Signature

7/30/20  
Date

NA  
Reviewed by

NA  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-006490  
Location = TOXL 8164.14.00 09/16  
07/30/2020 14:02

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
abcdefghijkl  
ABCDEFGHIJKL

16 char per line  
abcdefghijklmnop  
ABCDEFGHIJKLMN

21 chars per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNQRSTU

24 characters per line  
abcdefghijklmnopqrstuwx  
ABCDEFGHIJKLMNQRSTUWX

32 characters per line  
abcdefghijklmnopqrstuwx123456  
ABCDEFGHIJKLMNQRSTUWXYZ123456


42 characters per line  
abcdefghijklmnopqrstuwx1234567890abcdef  
ABCDEFGHIJKLMNQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup  
Data Entry Mode: Enabled  
Start Test Sequence: DABACABA  
Display Prelim Rslt? Yes  
Display Third Digit? Yes  
Inhib Printer(Y/N)? No  
Display Volume? No  
Disable On Menfull? Yes  
# of Print Copies? 1  
Select Std (D/W/I)? Dry  
Standard Value? 0.080  
Standard Lot #? 24119080A1  
Standard Cyl #? 9  
Standard Expiration? 11/05/2021  
Oper No? 666666

Flow Cal. Date: 08/18/2015  
Slope 675  
Intercept -635799

IR Calibration Date: 08/18/2015  
3um 9um  
-----  
0th Coef(\*100): -18073 -15729  
1st Coef(\*100): 265337 133887  
2nd Coef(\*100): 3965 1560  
H2O adj(ng/l\*10K): 767 457

\*\*\*\*\* Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks: INTERNAL  
PRINT  
TEST

Form 106-18000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                 8164.14.00 09/16  
07/30/2020                        14:03

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup


Data Entry Mode:                Enabled  
Start Test Sequence:            DABACABA  
Display Prelim Rslt?            Yes  
Display Third Digit?            Yes  
Inhib Printer(Y/N)?            No  
Display Volume?                No  
Disable On Memfull?            Yes  
# of Print Copies?              1  
Select Std (D/W/I)?            Dry  
Standard Value?                0.080  
Standard Lot #?                 24119080A1  
Standard Cyl #?                 9  
Standard Expiration?          11/05/2021  
Oper No?                        666666

Flow Cal. Date:                 08/18/2015  
Slope                            675  
Intercept                       -635799

IR Calibration Date:            08/18/2015  
                                    3um                            9um

-----  
0th Coef(\*100):                -18073                        -15729  
1st Coef(\*100):                265337                        133887  
2nd Coef(\*100):                3965                            1560  
H2O adj(mg/l\*10k):            767                             457

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:    EXTERNAL PRINT

TEST


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006490  
Location = TOXL      8164.14.00 09/16  
07/30/2020      14:05

Test	AC	Time
01 Room Air	0.000	14:06
02 *Subject Test	REF*	14:06
03 Room Air	0.000	14:07

\*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = R KEY REFUSAL  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: "R" KEY REFUSAL  
BARCODE READER

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                    8164.14.00 09/16  
07/30/2020                            14:08

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

DIAGNOSTIC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006490  
Location = TOXL      8164.14.00 09/16  
07/30/2020      14:17

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:17
02 Std. Sol.	0.000	14:18
03 Room Air	0.000	14:19
04 Std. Sol.	0.000	14:19
05 Room Air	0.000	14:20
06 Std. Sol.	0.000	14:21
07 Room Air	0.000	14:21

08 Sim Temp = 34.0°C

Simul Ser No = MP3002  
Std Sol No = MILLIQ H2O  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:      ZERO AC  
                         0.000 AC



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                    8164.14.00 09/16  
07/30/2020                            14:23

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:23
02 Std. Sol.	0.029	14:24
03 Room Air	0.000	14:25
04 Std. Sol.	0.029	14:26
05 Room Air	0.000	14:26
06 Std. Sol.	0.029	14:27
07 Room Air	0.000	14:27

08 Sim Temp = 34.0°C

Simul Ser No = MP3063  
Std Sol No = 20130  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:    *Low AC*  
                  *0.030 AC*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006490  
Location = TOXL      8164.14.00 09/16  
07/30/2020      14:29

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:29
02 Std. Sol.	0.295	14:30
03 Room Air	0.000	14:30
04 Std. Sol.	0.295	14:31
05 Room Air	0.000	14:32
06 Std. Sol.	0.295	14:33
07 Room Air	0.000	14:33

08 Sim Temp = 34.0°C

Simul Ser No = MP3069

Std Sol No = 201902D

County = 08

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                    8164.14.00 09/16  
07/30/2020                            14:35

Test	AC	Time
01 Room Air	0.000	14:36
02 *Subject Test	INT*	14:36
03 Room Air	0.000	14:37

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = INTERFERENT  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Interferent Check*  
*0.10AC + 0.05% ACETONE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                 8164.14.00 09/16  
07/30/2020                        14:38

Test	AC	Time
01 Room Air	0.000	14:39
02 *Subject Test	INT*	14:40
03 Room Air	0.000	14:40

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = INTERFERENT  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                    8164.14.00 09/16  
07/30/2020                            14:41

Test	AC	Time
01 Room Air	0.000	14:42
02 *Subject Test	INT*	14:42
03 Room Air	0.000	14:43

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = INTERFERENT  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Interferent Check*  
*0.10% METHANOL*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                    8164.14.00 09/16  
07/30/2020                            14:44

Test	AC	Time
01 Room Air	0.000	14:45
02 *Subject Test	INT*	14:45
03 Room Air	0.000	14:46

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = INTERFERENT  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10% ISOPROPAOL*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006490  
Location = TOXL      8164.14.00 09/16  
07/30/2020      14:48

Test	AC	Time
01 Room Air	AMB*	14:48
02 Room Air	0.000	14:49

\*Invalid Test  
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = AMBIENT FAIL  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*AMBIENT FAIL*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                8164.14.00 09/16  
07/30/2020                      14:49

Test	AC	Time
01 Room Air	0.000	14:50
02 Invalid Sample	X.XXX	14:51
03 Room Air	0.000	14:51

\*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                      Weight = 150  
Test = DUI                              Cit = INVALID SAMPLE  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08                              Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *INVALID SAMPLE*



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006490  
Location = TOXL      8164.14.00 09/16  
07/30/2020      14:52

Test	AC	Time
01 Room Air	0.000	14:53
02 *Subject Test	ABT*	14:53
03 Room Air	0.000	14:54

\*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = SEQUENCE ABORT  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *SEQUENCE ABORTED*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                    8164.14.00 09/16  
07/30/2020                            14:54

Test	AC	Time
01 Room Air	0.000	14:55
02 Subject Test 1	0.021	14:55
03 Room Air	0.000	14:56
04 Reported AC	0.021	14:55

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = ABA MODE  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*ABA MODE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                 8164.14.00 09/16  
07/30/2020                        14:57

Test	AC	Time
01 Diagnostic	OK	14:58
02 Room Air	0.000	14:58
03 Subject Test 1	0.007	14:58
04 Room Air	0.000	15:00
05 Std. Gas	0.088*	15:02
06 Room Air	0.000	15:02

\*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = CAL CHK OUT TOL  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*CALIBRATION CHECK  
OUT OF TOLERANCE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                    8164.14.00 09/16  
07/30/2020                            15:04

Test	AC	Time
01 Diagnostic	OK	15:05
02 Room Air	0.000	15:05
03 *Subject Test	RFI*	15:05
04 Room Air	0.000	15:06

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

RFI CHECK

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                    8164.14.00 09/16  
07/30/2020                            15:06

Test	AC	Time
01 Diagnostic	OK	15:07
02 Room Air	0.000	15:08
03 *Subject Test	0.018*	15:11
04 Room Air	0.000	15:13
05 Std. Gas	0.081	15:14
06 Room Air	0.000	15:16
07 *Subject Test	0.019*	15:19
08 Room Air	0.000	15:20

No RFI Detected

\*Deficient Sample - Value Printed was  
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = DEFICIENT SAMP  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08                            Oper No. = 666666

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the  
instructions displayed by the Intoxilyzer  
in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: DEFICIENT SAMPLES

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                    8164.14.00 09/16  
07/30/2020                            15:21

Test	AC	Time
01 Diagnostic	OK	15:22
02 Room Air	0.000	15:22
03 Reference	IPS*	15:22
04 Room Air	0.000	15:23

\*Improper Sample

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = IMPROPER SAMP  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*IMPROPER  
SAMPLE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                 8164.14.00 09/16  
07/30/2020                        15:24

Test	AC	Time
01 Diagnostic	OK	15:24
02 Room Air	0.000	15:25
03 Subject Test 1	0.019*	15:25
04 Room Air	0.000	15:27
05 Std. Gas	0.081	15:28
06 Room Air	0.000	15:30
07 Subject Test 2	0.146*	15:31
08 Room Air	0.000	15:32

\*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = DIFF TOO GREAT  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*DIFFERENCE*

*TOO GREAT*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                 8164.14.00 09/16  
07/30/2020                        15:35

Test	AC	Time
01 Diagnostic	OK	15:36
02 Room Air	0.000	15:37
03 Subject Test 1	0.020	15:37
04 Room Air	0.000	15:39
05 Std. Gas	0.080	15:40
06 Room Air	0.000	15:42
07 Subject Test 2	0.019	15:42
08 Room Air	0.000	15:43
09 Reported AC	0.019	15:42

Difference OK  
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = CMS MODE  
Dr. Lic. = ND/DIS921456  
Lot No = 24119080A1  
Cyl No = 9  
Expiration Date = 11/05/2021  
County = 08                          Oper No. = 666666

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*CMS MODE*



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006490  
Location = TOXL                 8164.14.00 09/16  
07/30/2020                        15:44

Test	AC	Time
01 Room Air	0.000	15:45
02 Std. Gas	0.081	15:45
03 Room Air	0.000	15:46
04 Std. Gas	0.080	15:46
05 Room Air	0.000	15:47
06 Std. Gas	0.081	15:47
07 Room Air	0.000	15:47
08 Std. Gas	0.081	15:48
09 Room Air	0.000	15:48
10 Std. Gas	0.082	15:49
11 Room Air	0.000	15:49
12 Std. Gas	0.081	15:50
13 Room Air	0.000	15:50
14 Std. Gas	0.081	15:51
15 Room Air	0.000	15:51
16 Std. Gas	0.081	15:51
17 Room Air	0.000	15:52
18 Std. Gas	0.081	15:52
19 Room Air	0.000	15:53
20 Std. Gas	0.081	15:53
21 Room Air	0.000	15:54
22 Std. Gas	0.081	15:54
23 Room Air	0.000	15:55
24 Std. Gas	0.081	15:55
25 Room Air	0.000	15:56
26 Std. Gas	0.081	15:56
27 Room Air	0.000	15:57
28 Std. Gas	0.081	15:57
29 Room Air	0.000	15:57
30 Std. Gas	0.081	15:58
31 Room Air	0.000	15:58
32 Std. Gas	0.081	15:59
33 Room Air	0.000	15:59
34 Std. Gas	0.082	16:00
35 Room Air	0.000	16:00
36 Std. Gas	0.081	16:01
37 Room Air	0.000	16:01
38 Std. Gas	0.081	16:02
39 Room Air	0.000	16:02
40 Std. Gas	0.082	16:02
41 Room Air	0.000	16:03
42 Cal Check Stats:		
43 Average	0.0811	
44 Std Dev	0.0004	
45 Rel Std Dev	0.56%	



Operator Signature  
CHARLES EDER

Remarks: 0.080 Ac  
LOT: 24119080A1  
CYL: 9  
Form 106-I8000  
Exp: 11.5.21

STABILITY TEST