



# INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

SFN 59281 (5-2011)

Serial Number <u>4951 80-004951</u> <u>DK</u>	Instrument Location <u>Mott ND</u>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

## Check When Done:

- ☒ 1. Surge Protector Installed/Properly Grounded.
- ☒ 2. Telephone Line Connected to Intoxilyzer® 8000.
- ☒ 3. Breath Tube Heated.
- ☒ 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- ☒ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- ☒ 7. Repair and/or Maintenance Performed (if Any): \_\_\_\_\_
- ☒ 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <u>[Signature]</u>	Date <u>5/4/12</u>
Reviewed By (Crime Laboratory Use Only) <u>Deb Resher</u>	Date <u>14 May 2012</u>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004951  
Location = MOTT      8164.13.00 06/09  
05/04/2012      15:22

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopQRSTUVWXYZ!@#%&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopQRSTUVWXYZ!@#%&\*()\_+?

Current Instrument Setup

Data Entry Mode:      Enabled  
Start Test Sequence:      DABACABA  
Display Prelim Rslt?      Yes  
Display Third Digit?      Yes  
Inhib Printer(Y/N)?      No  
Display Volume?      No  
Disable On Memfull?      Yes  
# of Print Copies?      1  
Select Std (D/W/I)?      Dry  
Standard Value?      0.080  
Standard Lot #?      23411080A1  
Standard Cyl #?      60  
Standard Expiration?      10/01/2013  
Oper No?      011119

Flow Cal. Date:      06/27/2011  
Slope      704  
Intercept      -660586

IR Calibration Date:      07/15/2011  
   3um      9um

0th Coef(*100):	-10545	-9200
1st Coef(*100):	258523	131513
2nd Coef(*100):	2314	1450
H2O adj(mg/l*10k):	588	659

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
Operator Signature  
JOSH MONSON

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004951  
Location = MOTT      8164.13.00 06/09  
05/04/2012      15:22

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:23
02 Std. Gas	0.082	15:23
03 Room Air	0.000	15:24
04 Std. Gas	0.081	15:24
05 Room Air	0.000	15:24
06 Std. Gas	0.080	15:25
07 Room Air	0.000	15:25

Lot No = 23411080A1

Cyl No = 60

Exp Date = 10/01/2013

County = 21

Oper No. = 011119

  
Operator Signature  
JOSH MONSON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004951  
Location = MOTT      8164.13.00 06/09  
05/04/2012      15:28

Test	AC	Time
01 Diagnostic	OK	15:30
02 Room Air	0.000	15:30
03 *Subject Test	RFI*	15:31
04 Room Air	0.000	15:31

\*Invalid Test  
Inhibited - RFI

Sub Name = MONSON, JOSHUA RUSSELL  
Sub DOB = 07/02/1982  
Sub Sex = Male      Weight = 240  
Test = DUI      Cit = NA  
Dr. Lic. = ND/MON825965  
Lot No = 23411080A1  
Cyl No = 60  
Expiration Date = 10/01/2013  
County = 21      Oper No. = 011119

I followed the Approved Method and the  
instructions displayed by the Intoxilyzer  
in conducting this test.

  
Operator Signature  
JOSH MONSON

Remarks: