



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
 SFN 59281 (06/2018)

Serial Number 80-006493	Instrument Location RZCHLAND COUNTY JAZZ
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): N/A
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature 	Date 06/11/21
-------------------------------	-------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By 	Certified Date 6-17-21
---------------------------	----------------------------------

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006493
Location = WAHP 8164.14.00 09/16
06/11/2021 21:18

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 02621080A1
Standard Cyl #? 33
Standard Expiration? 03/05/2023
Oper No? 132570

Flow Cal. Date: 05/23/2018
Slope 640
Intercept -624033

IR Calibration Date: 05/23/2018
 3um 9um

0th Coef(*100): -19857 -14695
1st Coef(*100): 263752 135087
2nd Coef(*100): 3253 1458
H2O adj(mg/l*10k): 643 483

***** Printer Test End *****



Operator Signature
CHRISTOPHER HUARD

Remarks: **PRZNT TEST FOR ZNSTALL**


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006493
Location = WAHP 8164.14.00 09/16
06/11/2021 21:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	21:21
02 Std. Gas	0.080	21:21
03 Room Air	0.000	21:22
04 Std. Gas	0.081	21:22
05 Room Air	0.000	21:23
06 Std. Gas	0.080	21:23
07 Room Air	0.000	21:24

Lot No = 02621080A1
Cyl No = 33
Exp Date = 03/05/2023
County = 39 Oper No. = 132570


Operator Signature
CHRISTOPHER HUARD

Remarks: **INSTALL**

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006493
Location = WAHP 8164.14.00 09/16
06/11/2021 21:28

Test	AC	Time
01 Diagnostic	OK	21:31
02 Room Air	0.000	21:32
03 *Subject Test	0.000*	21:35
04 Room Air	0.000	21:37
05 Std. Gas	0.081	21:38
06 Room Air	0.000	21:39
07 *Subject Test	0.000*	21:42
08 Room Air	0.000	21:43


No RFI Detected

*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = TEST, SUBJECT NA
Sub DOB = 01/01/1980
Sub Sex = Male Weight = 200
Test = OTH Cit = NA
Dr. Lic. = ND/TEST
Lot No = 02621080A1
Cyl No = 33
Expiration Date = 03/05/2023
County = 39 Oper No. = 132570

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
CHRISTOPHER HUARD

Remarks: *FAILED RFI TEST*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006493
Location = WAHP 8164.14.00 09/16
06/11/2021 21:44

Test	AC	Time
01 Diagnostic	OK	21:46
02 Room Air	0.000	21:46
03 *Subject Test	RFI*	21:46
04 Room Air	0.000	21:47

*Invalid Test
Inhibited - RFI

Sub Name = RFI, TEST NA
Sub DOB = 01/01/1980
Sub Sex = Male Weight = 200
Test = DUI Cit = NA
Dr. Lic. = ND/TEST
Lot No = 02621080A1
Cyl No = 33
Expiration Date = 03/05/2023
County = 39 Oper No. = 132570

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.


Operator Signature
CHRISTOPHER HUARD

Remarks: *RFI TEST*

Form 106-I8000