

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-007098

Location: LOGAN COUNTY
S.O. (NAPD)

A. General Setup and Checks:

1. Diagnostic Tests Pass and Instrument in Ready Mode
2. Breath tube heated
3. Date, time and location code (Level 2,E). Re-set if necessary.
4. Print test (Level 1,P). Sign and attach test record.
5. Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
 - a. Display: 996 psi Regulator: 1000psi

B. Tests (Sign and attach test records):

- Configure simulator for the following tests (Level 1,S).
1. Wet Calibration Check
 - a. Low AC (Level 1,C): Known Value ≤ 0.03 AC: 0.020 AC
Sim. Ser #: DR3375 Lot #: 17160 Exp. Date: 6.8.19
 - b. High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC
Sim. Ser #: DR7111 Lot #: 17350 Exp. Date: 10.11.19

Configure dry gas standard for the following tests (Level 1,S).

2. Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone
Sim. Ser #: DR7089 Lot #: ICS 5 Exp. Date: NA
 Display reads "Interferent Detect"
3. RFI Check (CMS Mode) Display reads "RFI Detect"
4. Dry Calibration Check (Level 1,C): Known Value 0.080 AC
Gas Cyl Lot #: 13518080 AG Cyl #: 24 Exp. Date: 8/5/2020

| | | |
|----------------------|---------------------|---------------------|
| Test 1 <u>0.081</u> | Test 4 <u>0.079</u> | Test 7 <u>0.080</u> |
| Test 2 <u>0.081</u> | Test 5 <u>0.080</u> | Test 8 <u>0.080</u> |
| Test 3 <u>0.080</u> | Test 6 <u>0.080</u> | Test 9 <u>0.080</u> |
| Average <u>0.080</u> | | |

C. Remarks/Maintenance: ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

[Signature]
Breath Analyst Signature

4.29.2019
Date

NA

NA

Reviewed by _____

Date _____

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007098
Location = NAPO 8164.16.00 09/18
04/29/2019 15:05

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 13518080A6
Standard Cyl #? 024
Standard Expiration? 08/05/2020
Oper No? 666666

Flow Cal. Date: 09/26/2018
Slope 684
Intercept -588446

IR Calibration Date: 09/26/2018
 3um 9um

0th Coef(*100): -40948 -15551
1st Coef(*100): 237555 134140
2nd Coef(*100): 2571 1433
H2O adj(mg/l*10k): 1097 457

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks: PRINT TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007098
Location = NAPO 8164.16.00 09/18
04/29/2019 15:28

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 15:29 |
| 02 Std. Sol. | 0.019 | 15:29 |
| 03 Room Air | 0.000 | 15:30 |
| 04 Std. Sol. | 0.019 | 15:31 |
| 05 Room Air | 0.000 | 15:31 |
| 06 Std. Sol. | 0.020 | 15:32 |
| 07 Room Air | 0.000 | 15:32 |


08 Sim Temp = 34.0°C

Simul Ser No = DR3375

Std Sol No = 17160

County = 24

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Low AC
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007098
Location = NAPO 8164.16.00 09/18
04/29/2019 15:34

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 15:34 |
| 02 Std. Sol. | 0.295 | 15:35 |
| 03 Room Air | 0.000 | 15:35 |
| 04 Std. Sol. | 0.297 | 15:36 |
| 05 Room Air | 0.000 | 15:37 |
| 06 Std. Sol. | 0.297 | 15:37 |
| 07 Room Air | 0.000 | 15:38 |


08 Sim Temp = 34.0°C

Simul Ser No = DR7111

Std Sol No = 17350

County = 24

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

HIGH AC
0.300 AC
Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007098
Location = NAPO 8164.16.00 09/18
04/29/2019 15:07

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 15:08 |
| 02 *Subject Test | INT* | 15:08 |
| 03 Room Air | 0.000 | 15:09 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 13518080A6
Cyl No = 024
Expiration Date = 08/05/2020
County = 24 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent Check
0.10 AC + 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007098
Location = NAPO 8164.16.00 09/18
04/29/2019 15:10

| Test | AC | Time |
|------------------|-------|-------|
| 01 Diagnostic | OK | 15:11 |
| 02 Room Air | 0.000 | 15:12 |
| 03 *Subject Test | RFI* | 15:12 |
| 04 Room Air | 0.000 | 15:12 |

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS821456
Lot No = 13518080A6
Cyl No = 024
Expiration Date = 08/05/2020
County = 24 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *RFI CHECK*

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007098
Location = NAPO 8164.16.00 09/18
04/29/2019 15:14

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 15:14 |
| 02 Std. Gas | 0.081 | 15:15 |
| 03 Room Air | 0.000 | 15:15 |
| 04 Std. Gas | 0.081 | 15:16 |
| 05 Room Air | 0.000 | 15:16 |
| 06 Std. Gas | 0.080 | 15:16 |
| 07 Room Air | 0.000 | 15:17 |

Lot No = 13518080A6
Cyl No = 024
Exp Date = 08/05/2020
County = 24 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007098
Location = NAPO 8164.16.00 09/18
04/29/2019 15:40

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 15:40 |
| 02 Std. Gas | 0.079 | 15:41 |
| 03 Room Air | 0.000 | 15:41 |
| 04 Std. Gas | 0.080 | 15:42 |
| 05 Room Air | 0.000 | 15:42 |
| 06 Std. Gas | 0.080 | 15:42 |
| 07 Room Air | 0.000 | 15:43 |

Lot No = 13518080A6
Cyl No = 024
Exp Date = 08/05/2020
County = 24 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007098
Location = NAPO 8164.16.00 09/18
04/29/2019 15:43

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 15:44 |
| 02 Std. Gas | 0.080 | 15:44 |
| 03 Room Air | 0.000 | 15:45 |
| 04 Std. Gas | 0.080 | 15:45 |
| 05 Room Air | 0.000 | 15:45 |
| 06 Std. Gas | 0.080 | 15:46 |
| 07 Room Air | 0.000 | 15:46 |

Lot No = 13518080A6
Cyl No = 024
Exp Date = 08/05/2020
County = 24 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK

0.080 AC

Form 106-I8000