

### INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number <b>20-08709</b> 6	Instrument Location Devils	416e
Reason for Install/Repair		
Install After Receiving From Crime Laboratory   Install	After Location Change	
Other (Specify)		
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3: Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Lo	cation: Level 2 Function F)	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S	5).	
6. Run Tests:		
X A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode	e or Level 1, Function B or C; Key Radio Du	uring Test).
7. Repair and/or Maintenance Performed (if any):	A	The state of the s
8. Complete the Top Portion of the Intoxilyzer® Record (SF	N50496, Form 120-G) and Place it by the I	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 12)	0-G) at the Intoxilyzer® Location at the Age	ency.
10. Send the Following to the Crime Laboratory: Completed 104-G), Print Test, ACA Test, and RFI Test.	Intoxilyzer® 8000 Installation and Repair C	Checkout (SFN59281, Form
Field Inspector/Signature		Inote of
rieid inspector ografia		Date 12/2/18
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved the date the Field Inspector performed the installation. This record county of Burleigh, North Dakota, is certified to be a true and correct	on file at the Office of Attorney General, Cri	
Reviewed/Certified By		Certified Date
- Roberta Lunger nummo		12/21/18
Ü		. 1
Note: Two print tests were recid.		
and another from 12/21/2018	at 10:12. Ryw 121	26/18

intoxilyzer test Record and Enecklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007096
Location = DEVL 8164.16.00 09/18
12/21/2018 10:09

\*\*\*\*\*\* Printer Test \*\*\*\*\*\*

# 12 char/line abcdefghijkl ABCDEFGHIJKL

## 16 char per line abcdefghijklmnop ABCDEFGHIJKLMNOP

## 21 chars per line abcdefghijklmnopqrstu ABCDEFGHIJKLMNOPQRSTU

24 characters per line abcdefghijklmnopgrstuvwx ABCDEFGHIJKLMNOPQRSTUVWX

32 characters per line abcdefghijklmnopgrstuvwxyz123456 ABCDEFGHIJKLMNOPQRSTUVWXYZ123456

42 characters per line abcdefghijklmnopgrstuvwyz123456789@abcoef ABCDEFGHIJKLMNOPURSTUVWXYZ123456789GABCDEF Current Instrument Setup Enabled Data Entry Mode: DABACABA Start Test Sequence: Yes Display Prelim Rsit? Yes Display Third Digit? No Inhib Printer(Y/N)? NG Display Volume? Yes Disable On Memfull? 1 # of Print Copies? Ory Select Std (D/W/I)? 0.080 Standard Value? 34917080Ĥ3 Standard Lot #? Standard Cyl #? 02/05/2020 Standard Expiration? 069108 Oper No? 10/02/2018 Flow Cal. Date: Slope -692605 Intercept 10/02/2018 IR Calibration Date: 3um Oth Coef(\*100): -33544 138054 1st Coef(\*100): 261464 2nd Coef(\*100): 3020 H2O adj(mg/l\*10k): 573

\*\*\*\*\*\*\*\*\* Test End \*\*\*\*\*\*\*

Operator Signature ROBERT KENNEDY

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007096
Location = DEVL 8164.16.00 09/18
12/21/2018 10:12

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34917080A3 Standard Cyl #? Standard Expiration? 02/05/2020 Oper No? 069108

Flow Cal. Date: 10/02/2018 Slope 695 Intercept -692605

IR Calibration Date: 10/02/2018
3um 9um

Oth Coef(\*100): -33544 -18475
1st Coef(\*100): 261464 138054
2nd Coef(\*100): 3020 1612
H2O adj(mg/l\*10k): 573 386

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature ROBERT KENNEDY

Remarks:

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = DEVL
12/21/2018

Alcohol Analyzer SN 80-007096 8164.16.00 09/18 10:22

#### DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:23
02 Std. Gas	0.080	10:23
03 Room Air	0.000	10:24
04 Std. Gas	0.080	10:24
05 Room Air	0.000	10:25
06 Std. Gas	0.079	10:25
07 Room Air	0.000	10:26

Lot No = 34917080A3

Cyl No = 3

Exp Date = 02/05/2020

County = 36

Oper No. = 069108

Operator Signature ROBERT KENNEDY

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = DEVL
12/21/2018

Alcohol Analyzer SN 80-007096 8164.16.00 09/18 10:31

## DRY CAL CHECK

Te	est	AC	Time
	Room Room	 RFI* 0.000	10:32 10:32

\*Invalid Test Inhibited - RFI

Lot No = 34917080A3

Cyl No = 3

Exp Date = 02/05/2020

County = 36

Oper No. = 069108

Operator Signature ROBERT KENNEDY

Remarks:

Form 106-I8000