



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59281 (06/2018)

| | |
|---|------------------------------------|
| Serial Number 80-007094 | Instrument Location CARS |
| Reason for Install/Repair | |
| <input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____ | |

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): **NONE**
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

| | |
|--|-----------------------|
| Field Inspector Signature <i>Serg Osh</i> | Date 1-3-19 |
|--|-----------------------|

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

| | |
|--|----------------------------------|
| Reviewed/Certified By <i>Roberta Mueger-nimmo</i> | Certified Date 1-17-19 |
|--|----------------------------------|

Note: Reason for Install/Repair - other. New Intoxilyzer location. Page 5 is a Diagnostic test sent with install paperwork. Page 6 is an additional RFI test ran on 1/16/19 by Terry Dillon. RSN 1/17/19.

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007094
Location = CARS 8164.16.00 09/18
01/03/2019 15:46

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 13518080A6
Standard Cyl #? 19
Standard Expiration? 08/05/2020
Oper No? 132318

Flow Cal. Date: 09/26/2018
Slope 672
Intercept -607823

IR Calibration Date: 09/26/2018
 3um 9um

0th Coef(*100): -28849 -15361
1st Coef(*100): 264571 134954
2nd Coef(*100): 3025 1343
H2O adj(mg/l*10k): 851 379

***** Printer Test End *****



Operator Signature
TERRY DILLON

Remarks:

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007094
Location = CARS 8164.16.00 09/18
01/03/2019 15:39

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 15:41 |
| 02 Std. Gas | 0.079 | 15:41 |
| 03 Room Air | 0.000 | 15:41 |
| 04 Std. Gas | 0.079 | 15:42 |
| 05 Room Air | 0.000 | 15:42 |
| 06 Std. Gas | 0.079 | 15:43 |
| 07 Room Air | 0.000 | 15:43 |

Lot No = 13518080A6
Cyl No = 19
Exp Date = 08/05/2020
County = 19 Oper No. = 132318


Operator Signature
TERRY DILLON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007094
Location = CARS 8164.16.00 09/18
01/03/2019 15:50

| Test | AC | Time |
|---------------|-------|-------|
| 01 Diagnostic | OK | 15:53 |
| 02 Room Air | RFI* | 15:53 |
| 03 Room Air | PUR** | 15:54 |

*Invalid Test
Inhibited - RFI
**Purge Fail

Sub Name = DILLO, JOHN DIL591111
Sub DOB = 05/14/1960
Sub Sex = Female Weight = NA
Test = DUI Cit = NA
Dr. Lic. = ND/DIL601960
Lot No = 13518080A6
Cyl No = 19
Expiration Date = 08/05/2020
County = 19 Oper No. = 132318

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
TERRY DILLON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007094
Location = CARS 8164.16.00 09/18
01/03/2019 16:22

DIAGNOSTIC

| | |
|---------------------------|------|
| Voltage/Current Test | Pass |
| RAM Test | Pass |
| EEPROM Checksum Test | Pass |
| Real Time Clock Test | Pass |
| DSP Test | Pass |
| Analytical Stability Test | Pass |
| Modem Test | Pass |
| Temp Regulation Test | Pass |

Oper No. = 132318



Operator Signature
TERRY TRREY

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007094
Location = CARS 8164.16.00 09/18
01/16/2019 13:59

| Test | AC | Time |
|---------------|-------|-------|
| 01 Diagnostic | OK | 14:03 |
| 02 Room Air | RFI* | 14:03 |
| 03 Room Air | 0.000 | 14:03 |

*Invalid Test
Inhibited - RFI

Sub Name = GILLON, NONE JOHN
Sub DOB = 05/14/1959
Sub Sex = Male Weight = NA
Test = DUI Cit = NA
Dr. Lic. = ND/GIL590000
Lot No = 13518080A6
Cyl No = 19
Expiration Date = 08/05/2020
County = 19 Oper No. = 132318

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
TERRY DILLON

Remarks:

Form 106-I8000