Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-007095 Location: ADAMS County A. General Setup and Checks: 5.6. (HETT)	11
 A. General Setup and Checks: Diagnostic Tests Pass and Instrument in Ready Mode Breath tube heated Date, time and location code (Level 2,E). Re-set if necessary. Print test (Level 1,P). Sign and attach test record. Tank monitor (Level 3,D,G). Display and Regulator within 50 psi. Display: Psi Regulator: 	
B. Tests (Sign and attach test records): Configure simulator for the following tests (Level 1,S). Wet Calibration Check a. ☐ Low AC (Level 1,C): Known Value ≤ 0.03 AC: Sim. Ser #: D 3375 Lot #: 1160 Exp. Date: 6.8.49 b. ☐ High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC Sim. Ser #: D 37111 Lot #: 17350 Exp. Date: 10.11.19 Configure dry gas standard for the following tests (Level 1,S). Interferent Check (Level 1,B): Known Value: 040 AC + 0.05% Acetone Sim. Ser #: D 7089 Lot #: 1755 Exp. Date: NA Display reads "Interferent Detect" A. ☐ Pry Calibration Check (Level 1,C): Known Value 0.080 AC Gas Cyl Lot #: 341708043 Cyl #: 144 Exp. Date: 2.5.20 Test 1 0.082 Test 4 0.082 Test 7 0.081 Test 2 0.082 Test 4 0.082 Test 9 0.082 Average 0.080	
C. Remarks/Maintenance: Awww. Juspecnow	
Instrument is acceptable to be used in the field. Breath Analyst Signature Reviewed by Date	

02/18

Issuing Authority - CEE

SN 80-007095

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007095
Location = HETT 8164.16.00 09/18
04/12/2019 11:57

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 34917080A3 Standard Expiration? 02/05/2020 Oper No? 666666

Flow Cal. Date: 09/27/2018 Slope 683 Intercept -637452

IR Calibration Date: 09/27/2018
3um 9um
----0th Coef(*100): -59629 -23097

1st Coef(*100): -59629 -23097 1st Coef(*100): 252379 137307 2nd Coef(*100): 2489 1549 H2O adj(mg/l*10k): 1414 594

******* Printer Test End *******

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = HETT
04/12/2019

Alcohol Analyzer SN 80-007095 8164.16.00 09/18 12:13

WET CAL CHECK

Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air 04 Std. Sol. 05 Room Air 06 Std. Sol. 07 Room Air	0.000 0.018 0.000 0.018 0.000 0.018 0.000	12:14 12:15 12:15 12:16 12:17 12:17

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = DR3375 Std Sol No = 17160 County = 01

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

0W 1

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = HETT
04/12/2019

Alcohol Analyzer SN 80-007095 8164.16.00 09/18 12:31

WET CAL CHECK

Test		AC	Time
01 Room 02 Std. 03 Room 04 Std. 05 Room 06 Std.	Sol. Air Sol. Air	0.000 0.300 0.000 0.301 0.000 0.301	12:32 12:32 12:33 12:34 12:34
07 Room		0.000	12:35

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = DR7111 Std Sol No = 17350

County = 01

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

O. 300 A

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = HETT 04/12/2019

Alcohol Analyzer SN 80-007095 8164.16.00 09/18 12:03

Te	est	AC ,	Time
02	Room Air	0.000	12:04
	Subject Test	INT	12:04
	Room Air	0.000	12:05

^{*}Invalid Test Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 01/01/1982

Sub Sex = Female

Weight = 150

Test = DUI

Cit = INTERFERENT

Dr. Lic. = ND/DIS821456

Lot No = 34917080A3

Cyl No = 14

Expiration Date = 02/05/2020

County = 01

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Interferent Check 0.10AC + 0.05% ACETONE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007095
Location = HETT 8164.16.00 09/18
04/12/2019 12:06

Т	est	AC	Time
02 03	Diagnostic Room Air *Subject Test Room Air	OK 0.000 RFI* 0.000	12:07 12:07 12:08 12:08

*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 01/01/1982

Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS821456

Lot No = 34917080A3

Cyl No = 14

Expiration Date = 02/05/2020

County = 01

Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = HETT 04/12/2019

Alcohol Analyzer SN 80-007095 8164.16.00 09/18 11:59

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:59
02 Std. Gas	0.082	12:00
03 Room Air	0.000 -	12:00
04 Std. Gas	0.082	12:01
05 Room Air	0.000	12:01
06 Std. Gas	0.081	12:01
07 Room Air	0.000	12:02

Lot No = 34917080A3

Cyl No = 14

Exp Date = 02/05/2020

County = 01

Oper No. = 666666

Operator Signature CHARLES EDER CAUBRATION CHECK

Remarks:

Form 106-I8000

0.080AC

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = HETT 04/12/2019 Alcohol Analyzer SN 80-007095 8164.16.00 09/18 12:37

DRY CAL CHECK

AC	Time
0.000	12:38
0.082	12:38
0.000	12:39
0.082	12:39
0.000	12:39
0.082	12:40
0.000	12:40
	0.000 0.082 0.000 0.082 0.000 0.082

Lot No = 34917080A3

Cyl No = 14

Exp Date = 02/05/2020

County = 01

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = HETT
04/12/2019

Alcohol Analyzer SN 80-007095 8164.16.00 09/18 12:40

DRY CAL CHECK

Test	AC .	Time
01 Room Air	0.000	12:41
02 Std. Gas	0.081	12:41
03 Room Air	0.000	12:42
04 Std. Gas	0.081	12:42
05 Room Air	0.000	12:43
06 Std. Gas	0.082	12:43
07 Room Air	0.000	12:44

Lot No = 34917080A3

Cyl No = 14

Exp Date = 02/05/2020

County = 01

Oper No. = 666666

Operator Signature
CHARLES EDER

CHARLES EDER

CHECKE

Remarks:

Form 106-I8000

0.08040