

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-007094

Location: GRANT COUNTY
S.O. (CARS)

A. General Setup and Checks:

1. Diagnostic Tests Pass and Instrument in Ready Mode
2. Breath tube heated
3. Date, time and location code (Level 2,E). Re-set if necessary.
4. Print test (Level 1,P). Sign and attach test record.
5. Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
 - a. Display: 1002 psi Regulator: 1000 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check

- a. Low AC (Level 1,C): Known Value ≤ 0.03 AC: 0.020 AC
Sim. Ser #: DR3375 Lot #: 17160 Exp. Date: 6.8.19
- b. High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC
Sim. Ser #: DR7111 Lot #: 17350 Exp. Date: 10.11.19

Configure dry gas standard for the following tests (Level 1,S).

2. Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone
Sim. Ser #: DR7089 Lot #: ICS 5 Exp. Date: NA
 Display reads "Interferent Detect"
3. RFI Check (CMS Mode) Display reads "RFI Detect"
4. Dry Calibration Check (Level 1,C): Known Value 0.080 AC
Gas Cyl Lot #: 13518080 A6 Cyl #: 19 Exp. Date: 8.5.20

| | | | | | |
|---------|--------------|--------|--------------|--------|--------------|
| Test 1 | <u>0.080</u> | Test 4 | <u>0.080</u> | Test 7 | <u>0.080</u> |
| Test 2 | <u>0.080</u> | Test 5 | <u>0.079</u> | Test 8 | <u>0.079</u> |
| Test 3 | <u>0.079</u> | Test 6 | <u>0.079</u> | Test 9 | <u>0.079</u> |
| Average | <u>0.079</u> | | | | |

C. Remarks/Maintenance: ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

Charles E. Ede
Breath Analyst Signature

4.17.2019
Date

NA

NA

Reviewed by _____

Date _____

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007094
Location = CARS 8164.16.00 09/18
04/17/2019 11:58

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 13518080A6
Standard Cyl #? 19
Standard Expiration? 08/05/2020
Oper No? 666666

Flow Cal. Date: 09/26/2018
Slope 672
Intercept -607823

IR Calibration Date: 09/26/2018
 3um 9um

0th Coef(*100): -28849 -15361
1st Coef(*100): 264571 134954
2nd Coef(*100): 3025 1343
H2O adj(mg/l*10k): 851 379

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

PRINT TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007094
Location = CARS 8164.16.00 09/18
04/17/2019 12:13

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 12:14 |
| 02 Std. Sol. | 0.019 | 12:15 |
| 03 Room Air | 0.000 | 12:15 |
| 04 Std. Sol. | 0.019 | 12:16 |
| 05 Room Air | 0.000 | 12:17 |
| 06 Std. Sol. | 0.019 | 12:17 |
| 07 Room Air | 0.000 | 12:18 |

08 Sim Temp = 34.0°C

Simul Ser No = DR3375

Std Sol No = 17160

County = 19

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Low AC
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007094
Location = CARS 8164.16.00 09/18
04/17/2019 12:24

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 12:24 |
| 02 Std. Sol. | 0.297 | 12:25 |
| 03 Room Air | 0.000 | 12:26 |
| 04 Std. Sol. | 0.297 | 12:26 |
| 05 Room Air | 0.000 | 12:27 |
| 06 Std. Sol. | 0.298 | 12:28 |
| 07 Room Air | 0.000 | 12:28 |


08 Sim Temp = 34.0°C

Simul Ser No = DR7111

Std Sol No = 17350

County = 19

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

HIGH AC
0.300 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007094
Location = CARS 8164.16.00 09/18
04/17/2019 12:04

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 12:04 |
| 02 *Subject Test | INT* | 12:05 |
| 03 Room Air | 0.000 | 12:06 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 13518080A6
Cyl No = 19
Expiration Date = 08/05/2020
County = 19 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Interferent Check*
0.10AC + 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007094
Location = CARS 8164.16.00 09/18
04/17/2019 12:06

| Test | AC | Time |
|------------------|-------|-------|
| 01 Diagnostic | OK | 12:07 |
| 02 Room Air | 0.000 | 12:07 |
| 03 *Subject Test | RFI* | 12:08 |
| 04 Room Air | 0.000 | 12:08 |

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS821456
Lot No = 13518080A6
Cyl No = 19
Expiration Date = 08/05/2020
County = 19 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *RFI CHECK*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007094
Location = CARS 8164.16.00 09/18
04/17/2019 12:00

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 12:01 |
| 02 Std. Gas | 0.080 | 12:01 |
| 03 Room Air | 0.000 | 12:02 |
| 04 Std. Gas | 0.080 | 12:02 |
| 05 Room Air | 0.000 | 12:02 |
| 06 Std. Gas | 0.079 | 12:03 |
| 07 Room Air | 0.000 | 12:03 |

Lot No = 13518080A6
Cyl No = 19
Exp Date = 08/05/2020
County = 19 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007094
Location = CARS 8164.16.00 09/18
04/17/2019 12:08

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 12:09 |
| 02 Std. Gas | 0.080 | 12:09 |
| 03 Room Air | 0.000 | 12:10 |
| 04 Std. Gas | 0.079 | 12:10 |
| 05 Room Air | 0.000 | 12:11 |
| 06 Std. Gas | 0.079 | 12:11 |
| 07 Room Air | 0.000 | 12:12 |

Lot No = 13518080A6
Cyl No = 19
Exp Date = 08/05/2020
County = 19 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007094
Location = CARS 8164.16.00 09/18
04/17/2019 12:30

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 12:31 |
| 02 Std. Gas | 0.080 | 12:31 |
| 03 Room Air | 0.000 | 12:32 |
| 04 Std. Gas | 0.079 | 12:32 |
| 05 Room Air | 0.000 | 12:33 |
| 06 Std. Gas | 0.079 | 12:33 |
| 07 Room Air | 0.000 | 12:34 |

Lot No = 13518080A6
Cyl No = 19
Exp Date = 08/05/2020
County = 19 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080AC

Form 106-I8000