

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-007091

Location: HETTINGER Co.  
S.O. (MOTT)

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 1005 psi Regulator: 1000 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check

- a.  Low AC (Level 1,C): Known Value  $\leq$  0.03 AC: 0.020 AC  
Sim. Ser #: DR3375 Lot #: 17160 Exp. Date: 6.8.19
- b.  High AC (Level 1,C): Known Value  $\geq$  0.25 AC: 0.300 AC  
Sim. Ser #: DR7111 Lot #: 17350 Exp. Date: 10.11.19

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR7089 Lot #: ICS 5 Exp. Date: NA  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 13518080A6 Cyl #: 1R Exp. Date: 8.5.20

Test 1	<u>0.078</u>	Test 4	<u>0.079</u>	Test 7	<u>0.078</u>
Test 2	<u>0.078</u>	Test 5	<u>0.079</u>	Test 8	<u>0.078</u>
Test 3	<u>0.078</u>	Test 6	<u>0.078</u>	Test 9	<u>0.078</u>
Average	<u>0.078</u>				

C. Remarks/Maintenance: ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

Charles E. Edr  
Breath Analyst Signature

4.17.2019  
Date

NA  
Reviewed by

NA  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-007091  
Location = MOTT                    8164.16.00 09/18  
04/17/2019                            13:31

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            13518080A6  
Standard Cyl #?            018  
Standard Expiration?       08/05/2020  
Oper No?                    666666

Flow Cal. Date:            09/25/2018  
Slope                        681  
Intercept                   -651030

IR Calibration Date:       09/25/2018  
                                    3um                      9um

-----  
0th Coef(\*100):        -39121                    -20906  
1st Coef(\*100):        240667                    132345  
2nd Coef(\*100):        3960                      1516  
H2O adj(mg/l\*10k):    899                        428

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

*PRINT TEST*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-007091  
Location = MOTT 8164.16.00 09/18  
04/17/2019 13:47

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:48
02 Std. Sol.	0.018	13:48
03 Room Air	0.000	13:49
04 Std. Sol.	0.018	13:50
05 Room Air	0.000	13:50
06 Std. Sol.	0.018	13:51
07 Room Air	0.000	13:51


08 Sim Temp = 34.0°C

Simul Ser No = DR3375

Std Sol No = 17160

County = 21

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC  
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-007091  
Location = MOTT      8164.16.00 09/18  
04/17/2019      13:57

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:57
02 Std. Sol.	0.298	13:58
03 Room Air	0.000	13:59
04 Std. Sol.	0.299	13:59
05 Room Air	0.000	14:00
06 Std. Sol.	0.299	14:01
07 Room Air	0.000	14:01

08 Sim Temp = 34.0°C

Simul Ser No = DR7111

Std Sol No = 17350

County = 21

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-007091  
Location = MOTT                                      8164.16.00 09/18  
04/17/2019    13:39

Test	AC	Time
01 Room Air	0.000	13:40
02 *Subject Test	INT*	13:40
03 Room Air	0.000	13:41

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                                      Weight = 150  
Test = DUI    Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 13518080A6  
Cyl No = 018  
Expiration Date = 08/05/2020  
County = 21    Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*

*0.10AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-007091  
Location = MOTT      8164.16.00 09/18  
04/17/2019      13:36

Test	AC	Time
01 Diagnostic	OK	13:37
02 Room Air	0.000	13:38
03 *Subject Test	RFI*	13:38
04 Room Air	0.000	13:38

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 13518080A6  
Cyl No = 018  
Expiration Date = 08/05/2020  
County = 21      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

RFI CHECK

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-007091  
Location = MOTT      8164.16.00 09/18  
04/17/2019      13:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:27
02 Std. Gas	0.078	13:28
03 Room Air	0.000	13:28
04 Std. Gas	0.078	13:29
05 Room Air	0.000	13:29
06 Std. Gas	0.078	13:30
07 Room Air	0.000	13:30

Lot No = 13518080A6  
Cyl No = 018  
Exp Date = 08/05/2020  
County = 21      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-007091  
Location = MOTT      8164.16.00 09/18  
04/17/2019      13:42

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:43
02 Std. Gas	0.079	13:43
03 Room Air	0.000	13:44
04 Std. Gas	0.079	13:44
05 Room Air	0.000	13:44
06 Std. Gas	0.078	13:45
07 Room Air	0.000	13:45

Lot No = 13518080A6  
Cyl No = 018  
Exp Date = 08/05/2020  
County = 21      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-007091  
Location = MOTT      8164.16.00 09/18  
04/17/2019      14:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:05
02 Std. Gas	0.078	14:05
03 Room Air	0.000	14:05
04 Std. Gas	0.078	14:06
05 Room Air	0.000	14:06
06 Std. Gas	0.078	14:07
07 Room Air	0.000	14:07

Lot No = 13518080A6  
Cyl No = 018  
Exp Date = 08/05/2020  
County = 21      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000