

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION Toxicology Section/Breath Alcohol Program SFN 59281 (10-2016)

Serial Number		Instrument Location				
80	-007092	UND-PD				
Reason for In	stall/Repair					
Install After Receiving From Crime Laboratory Install After Location Change						
(V) Other (Specify) After yearly check by C. Eder						
Check When	Done:					
VZm1.	Surge Protector Installed/Properly Grounded.					
WZM 2.	Telephone Line Connected to Intoxilyzer® 8000.					
VIRM 3.	Breath Tube Heated.					
Vizn 4.	Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
V ZM 5.	Scan/Enter Gas Cylinder Information (Level 1, Function S).					
Unm 6.	Run Tests:					
• •	A. Print Test (Level 1, Function P).71-W					
	B. ACA Test (Level 1, Function C). Thm					
(C. RFI Test (CMS Mode or Level 1, Fund	ction B or C; Key Radio During	Test).72M			
□2m 7.	Repair and/or Maintenance Performed (if	Any): None				
Fr. 1 8.	Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.					
Vicu 9.	File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.					
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkou (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.						
Field Inspector	Signature Mudani		Date 5-14-2019			
Crime Laborate	ry Use Only					
alcohol concent Attorney Genera	has been reviewed and the instrument is a ration from the date the Field Inspector per al, Crime Laboratory Division, in the County uments received.	formed the installation. This re	cord on file at the Office of ertified to be a true and correct			
Reviewed/Cert	ified by Charles Ed	4	Certified Date 5-20-2019			

Form 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007092
Location = UNDP 8164.16.00 09/18
05/14/2019 12:14

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #?
Standard Cyl #? 13518080A6 17 Standard Expiration? 08/05/2020 Oper No? 119409

Flow Cal. Date: 09/26/2018 Slope 670 Intercept -596123

IR Calibration Date: 09/26/2018

3um 9um

Oth Coef(*100): -27453 -17291

1st Coef(*100): 248629 134601

2nd Coef(*100): 2770 1393

H2O adj(mg/l*10k): 872 414

****** Printer Test End *******

Operator Signature TRACY MEIDINGER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-007092 Location = UNDP 8164.16.00 09/18 05/14/2019

12:15

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:15
02 Std. Gas	0.081	12:16
03 Room Air	0.000	12:16
04 Std. Gas	0.081	12:17
05 Room Air	0.000	12:17
06 Std. Gas	0.080	12:18
07 Room Air	0.000	12:18

Lot No = 13518080A6

Cyl No = 17

Exp Date = 08/05/2020

County = 18

Oper No. = 119409

Remarks:

Form 106-I8000

Operator Signature TRACY MEIDINGER

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-007092
Location = UNDP	8164.16.00 09/18
05/14/2019	12:19

Test	AC	Time
01 Room Air 02 Subject Test 1 03 Room Air 04 Room Air	0.000 0.000 RFI* 0.000	12:21 12:22 12:22 12:23

^{*}Invalid Test Inhibited - RFI

Sub Name = JANE, DOE NONE

Sub DOB = 01/01/1999

Sub Sex = Male Weight = 190 Test = DUI Cit = 111111111

Dr. Lic. = ND/111111111

Lot No = 13518080A6

Cyl No = 17

Expiration Date = 08/05/2020

> Operator Signature TRACY MEIDINGER

Remarks:

Form 106-I8000