

INTOXILYZER® 8000 CALIBRATION ADJUSTMENT

Intoxilyzer® 8000 Serial Number: 80-007093 Location: TOXL

*No Flow Sensor
 CALIBRATION.
 FLOW METERS
 HAVE BEEN SENT
 IN FOR CALIBRATION
 CHECK
 CEE
 12/16/19*

- A. *NA CEE 12/16/19* Flow Sensor Calibration and Verification Check (Level 3,M,C,F)
1. Replaced o-rings if damaged
 2. Flow Meter Serial Number: _____
 3. Air Supplied to Intoxilyzer® 8000 at:
 - a. 5 L/min 15 L/min 30 L/min
 4. Flow Rate Calibration Printout Attached
 - a. Correlation ≥ 0.99000
 5. Flow Sensor Calibration Verification (Level 3,D,F)
 - a. 10 L/min: 0. ___ L/S X 60 Sec/min = _____ L/min
 - b. 20 L/min: 0. ___ L/S X 60 Sec/min = _____ L/min
 - c. Flow Rates within ± 1 L/min of Expected Value

- B. Gas Tank Sensor Check (Level 3,D,G)
1. Display: 755 psi Regulator: 750 psi
 2. Display and Regulator within 50 psi
 3. Completed tare of tank sensor if needed (Level 3,M,C,G)

- C. Optical Bench Calibration and Verification Check (Level 3,M,C,O)
1. Autocalibration Printout Attached
 - a. Max Power Res Value ≥ 10
 - b. Auto Range Res Value ≥ 4
 2. Simulator Solutions for Optical Bench Calibration Adjustment
 - a. Set # Solutions to Run at 5

Soln.	g/210 L	Lot No.	Exp. Date	Simulator SN
1	0.000 <i>(STD) ACTUAL</i>	NA - MilliQ H ₂ O	NA - MilliQ H ₂ O	MP5289
2	0.040 <i>(0.040)</i>	201808D	8.22.20	MP5319
3	0.081 <i>(0.080)</i>	201807C	7.25.20	MP5290
4	0.151 <i>(0.150)</i>	201811E	11.26.20	MP5320
5	0.301 <i>(0.300)</i>	201803H	3.22.20	MP5321

3. 0.100 AC Calibration Gas for H₂O Adjustment
 - a. Lot No. 13518100A3 Cyl No. 6 Exp. Date: 8.5.20
4. Atmospheric Pressure
 - a. 957 mbar Displayed by Intoxilyzer® 8000
 - b. 958 mbar Adjusted to using barometer
 - c. 958 mbar on Auto Calibration Report printout
5. Screen displayed "Calibration Success"

6. Calibration Adjustment Printout Attached
- a. Solution 1 Avg % Abs ≤ 0.2500
 - b. Solution 2-5 REL STD DEV ≤ 3.000
 - c. Residual (g/210 L) Values for Solutions 1-5 ≤ 0.0020 for 3 μm and 9 μm channels
 - d. Dry Gas H₂O Adjustment Sum for 3 μm and 9 μm channels within ± 10

	Average		H ₂ O Adjust		
3 μm	<u>4392</u>	+	<u>369</u>	=	<u>4761</u>
9 μm	<u>4140</u>	+	<u>621</u>	=	<u>4761</u>

7. Optical Bench Calibration Verification (Level 1, S and C)

a. Wet Calibration Check

i. Low AC Known Value ≤ 0.03 AC: 0.020 AC
Sim. SN: DR3378 Lot No.: 18020 Exp. Date: 1.9.20

ii. High AC Known Value ≥ 0.25 AC: 0.250 AC
Sim. SN: DR7351 Lot No.: 2018036 Exp. Date: 3.22.20

b. Dry Calibration Check: Known Value 0.08 AC

Lot No. 34917080A3 Cyl No. 2 Exp. Date: 2/5/20

Test 1 0.080 AC Test 4 0.080 AC Test 7 0.080 AC

Test 2 0.079 AC Test 5 0.079 AC Test 8 0.080 AC

Test 3 0.080 AC Test 6 0.079 AC Test 9 0.080 AC

Average 0.080 AC

- c. Wet Calibration Check and Dry Calibration Check AC results are within ± 0.005 or $\pm 5\%$ (whichever is greater) of stated value.

D. Remarks/Maintenance: CALIBRATION ADJUSTMENT DUE TO

SOME 0.084 AC RETURNS FOR THE 0.080 AC STANDARD.

0.084 AC VALUE STILL WITHIN TOLERANCE.

REPLACED GAS REGULATOR. WIRES FRAYED, RUBBER AGAINST CASE.

EXPOSED.

Instrument is acceptable to be used in the field.



Breath Analyst Signature

NA

Reviewed by

12/16/2019

Date

NA

Date

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-007093
 12/16/2019 15:37:01

Auto Calibration

pg 1 of 2

	<<<<< 3um >>>>>		<<<<< 9um >>>>>	
	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)

Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	0.0630	(-0.0080)	0.2370	(-0.0120)
Sample #2	0.0670	(0.0280)	0.2270	(-0.0020)
Sample #3	0.0210	(0.0680)	0.2060	(0.0120)
Sample #4	0.0590	(0.0760)	0.2370	(0.0110)
Avg % Abs	0.0490	(0.0573)	0.2233	(0.0070)
STD DEV	0.0246	(0.0257)	0.0158	(0.0078)
REL STD DEV	50.156	(44.854)	7.084	(111.575)

Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	0.7780	(-0.0150)	1.6250	(-0.0220)
Sample #2	0.7680	(0.0060)	1.6000	(0.0000)
Sample #3	0.7900	(0.0160)	1.5990	(0.0090)
Sample #4	0.7460	(0.0460)	1.5630	(0.0200)
Avg % Abs	0.7680	(0.0227)	1.5873	(0.0097)
STD DEV	0.0220	(0.0208)	0.0211	(0.0100)
REL STD DEV	2.865	(91.838)	1.328	(103.621)

Solution = 0.081 g/210L or 0.3857 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	1.5070	(-0.0290)	2.9610	(-0.0250)
Sample #2	1.5300	(-0.0120)	2.9540	(0.0140)
Sample #3	1.4830	(0.0110)	2.9170	(0.0210)
Sample #4	1.5090	(0.0010)	2.9400	(0.0210)
Avg % Abs	1.5073	(0.0000)	2.9370	(0.0187)
STD DEV	0.0235	(0.0115)	0.0187	(0.0040)
REL STD DEV	1.562	(0.000)	0.636	(21.651)

Solution = 0.151 g/210L or 0.7190 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	2.6890	(-0.0050)	5.1660	(0.0070)
Sample #2	2.6560	(0.0280)	5.1240	(0.0410)
Sample #3	2.6850	(0.0310)	5.1380	(0.0420)
Sample #4	2.6850	(0.0310)	5.1430	(0.0300)
Avg % Abs	2.6753	(0.0300)	5.1350	(0.0377)
STD DEV	0.0167	(0.0017)	0.0098	(0.0067)
REL STD DEV	0.626	(5.773)	0.192	(17.677)

Solution = 0.301 g/210L or 1.4333 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	5.1790	(0.0000)	9.6520	(0.0000)
Sample #2	5.1930	(0.0150)	9.6490	(0.0280)
Sample #3	5.1640	(0.0280)	9.6420	(0.0300)
Sample #4	5.1590	(0.0300)	9.6280	(0.0440)
Avg % Abs	5.1720	(0.0243)	9.6397	(0.0340)
STD DEV	0.0184	(0.0081)	0.0107	(0.0087)
REL STD DEV	0.355	(33.471)	0.111	(25.641)

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-007093
 12/16/2019 15:37:01

Auto Calibration

pg 2 of 2

<<<<< 3um >>>>>			<<<<< 9um >>>>>		
-----			-----		
Zero Order Coef	-146.10		Zero Order Coef	-323.56	
First Order Coef	2646.15		First Order Coef	1389.33	
Second Order Coef	29.92		Second Order Coef	13.64	
-----			-----		
Act	Fit	Residual	Act	Fit	Residual
(g/210L)	(g/210L)	(g/210L)	(g/210L)	(g/210L)	(g/210L)
0.000	-0.000	0.0003	0.000	-0.000	0.0003
0.040	0.040	0.0000	0.040	0.040	-0.0002
0.081	0.082	-0.0011	0.081	0.081	-0.0004
0.151	0.150	0.0009	0.151	0.151	0.0004
0.301	0.301	-0.0001	0.301	0.301	-0.0001

<<<<< 3um >>>>>		<<<<< 9um >>>>>	
-----		-----	
Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1			
Sample			
Sample #1	4486.00	4221.00	
Sample #2	4351.00	4147.00	
Sample #3	4428.00	4159.00	
Sample #4	4397.00	4116.00	
Avg	4392.0000	4140.6665	
STD DEV	38.7427	22.1886	
REL STD DEV	0.882	0.536	
H2O adjust (mg/l*10k)	369	621	

Atmospheric Pressure = 958

*****CALIBRATION SUCCESSFUL*****

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-007093
 12/16/2019 15:37:01

Auto Calibration
 Max Power Res Value = 95
 Au+ - Power Res Value = 77



Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007093
Location = TOXL 8164.16.00 09/18
12/16/2019 16:21

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:21
02 Std. Sol.	0.020	16:22
03 Room Air	0.000	16:23
04 Std. Sol.	0.020	16:23
05 Room Air	0.000	16:24
06 Std. Sol.	0.020	16:25
07 Room Air	0.000	16:25

08 Sim Temp = 34.0°C

Simul Ser No = DR3378

Std Sol No = 18020

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

Low AC
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007093
Location = TOXL 8164.16.00 09/18
12/16/2019 16:26

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:26
02 Std. Sol.	0.246	16:27
03 Room Air	0.000	16:28
04 Std. Sol.	0.247	16:28
05 Room Air	0.000	16:29
06 Std. Sol.	0.248	16:30
07 Room Air	0.000	16:30

08 Sim Temp = 34.0°C

Simul Ser No = DR7351

Std Sol No = 201803G

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

HIGH AC
0.250 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007093
Location = TOXL 8164.16.00 09/18
12/16/2019 16:31

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:31
02 Std. Gas	0.080	16:32
03 Room Air	0.000	16:32
04 Std. Gas	0.079	16:33
05 Room Air	0.000	16:33
06 Std. Gas	0.080	16:34
07 Room Air	0.000	16:34

Lot No = 34917080A3
Cyl No = 2
Exp Date = 02/05/2020
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007093
Location = TOXL 8164.16.00 09/18
12/16/2019 16:34

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:35
02 Std. Gas	0.080	16:35
03 Room Air	0.000	16:36
04 Std. Gas	0.079	16:36
05 Room Air	0.000	16:37
06 Std. Gas	0.079	16:37
07 Room Air	0.000	16:38

Lot No = 34917080A3
Cyl No = 2
Exp Date = 02/05/2020
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007093
Location = TOXL 8164.16.00 09/18
12/16/2019 16:38

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:38
02 Std. Gas	0.080	16:39
03 Room Air	0.000	16:39
04 Std. Gas	0.080	16:40
05 Room Air	0.000	16:40
06 Std. Gas	0.080	16:40
07 Room Air	0.000	16:41

Lot No = 34917080A3
Cyl No = 2
Exp Date = 02/05/2020
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

CALIBRATION CHECK
0.080 AC

Form 106-I8000