Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

# **INTOXILYZER® 8000 INSPECTION**

Intoxil	ilyzer® 8000 Serial Number: <u>80-00 <del>7</del>0</u> 95	Location: ToxL
Α.	General Setup and Checks:  1. Diagnostic Tests Pass and Instrument in Breath tube heated  3. Date, time and location code (Level 2,E).  4. Print test (Level 1,P). Sign and attach test Tank monitor (Level 3,D,G). Display and a. Display: 352 psi Regulator: 350 psi	Re-set if necessary. et record. Regulator within 50 psi.
B.	Tests (Sign and attach test records):  Configure simulator for the following tests (Level 1. Wet Calibration Check  a.	≤ 0.03 AC: <u>0 · 0 /5</u> AC 5C Exp. Date: 5·30·20
	<ul> <li>Configure dry gas standard for the following test</li> <li>Interferent Check (Level 1,B): Known Value Sim. Ser #: DR3451 Lot #: TCS 6 Explose Explose</li></ul>	reads "RFI Detect" Note: WA  Test 7 0.083  Test 8 0.083
C.	Remarks/Maintenance: ANNUAL INTENT WAS LEFT ON MOUNTAIN TIME	
	- NO THE	E DWAING THSTECTION.
Breath	trument is acceptable to be used in the field.  Analyst Signature  Wed by	5. /2.2020 Date NA Date

02/18 Issuing Authority - CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007095
Location = TOXL 8164.16.00 09/18
05/12/2020 12:08

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 05620080A1 Standard Cyl #? 47 Standard Expiration? 04/05/2022 Oper No? 666666

Flow Cal. Date: 09/27/2018 Slope 683 Intercept -637452

IR Calibration Date: 09/27/2018
3um 9um

Oth Coef(\*100): -59629 -23097 1st Coef(\*100): 252379 137307 2nd Coef(\*100): 2489 1549 H2O adj(mg/l\*10k): 1414 594

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature CHARLES EDER

Remarks: PRINT TEST

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 05/12/2020

Alcohol Analyzer SN 80-007095 8164.16.00 09/18 12:10

### WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:10
02 Std. Sol.	0.012	12:11
03 Room Air	0.000	12:12
04 Std. Sol.	0.013	12:12
05 Room Air	0.000	12:13
06 Std. Sol.	0.013	12:14
07 Room Air	0.000	12:14

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = MP3064 Std Sol No = 201805C

County = 08

Oper No. = 666666

Operator Signature

CHARLES EDER LOW AC

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
05/12/2020

Alcohol Analyzer SN 80-007095 8164.16.00 09/18 12:16

### WET CAL CHECK

Test			AC	Time
01	Room	Air	0.000	12:16
02	Std.	Sol.	0.303	12:17
03	Room	Air	0.000	12:18
04	Std.	Sol.	0.305	12:18
05	Room	Air	0.000	12:19
06	Std.	Sol.	0.305	12:20
07	Room	Air	0.000	12:20

08 Sim Temp = 34.0°C

Simul Ser No = MP3069 Std Sol No = 201902D

County = 08

Oper No. = 666666

Operator Signature

CHARLES EDER

Remarks:

0.300 AC Form 106-18000

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007095
Location = TOXL 8164.16.00 09/18
05/12/2020 12:23

Test		AC	Time
01	Room Air	0.000	12:23
02	*Subject Test	INT*	12:24
03	Room Air	0.000	12:24

#### \*Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT

Dr. Lic. = ND/DIS921456

Lot No = 05620080A1

Cyl No = 47

Expiration Date = 04/05/2022

> Operator Signature CHARLES EDER

Remarks:

Interferent Check 0.10AC + 0.05% ACETONE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007095
Location = TOXL 8164.16.00 09/18
05/12/2020 12:25

Test	AC	Time
01 Diagnostic	OK	12:26
02 Room Air	0.000	12:26
03 *Subject Test	RFI*	12:27
04 Room Air	0.000	12:27

\*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS921456

Lot No = 05620080A1

Cyl No = 47

Expiration Date = 04/05/2022

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

RFI CHECK

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
05/12/2020

Alcohol Analyzer SN 80-007095 8164.16.00 09/18 12:28

### DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:28
02 Std. Gas	0.083	12:28
03 Room Air	0.000	12:29
04 Std. Gas	0.083	12:29
05 Room Air	0.000	12:30
06 Std. Gas	0.083	12:30
07 Room Air	0.000	12:31

Lot No = 05620080A1

Cyl No = 47

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

1 ALIBRATION CHECK -0.083 A-CEE 0.080 AC Form 106-18000

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
05/12/2020

Alcohol Analyzer SN 80-007095 8164.16.00 09/18 12:31

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:32
02 Std. Gas	0.083	12:32
03 Room Air	0.000	12:33
04 Std. Gas	0.083	12:33
05 Room Air	0.000	12:34
06 Std. Gas	0.083	12:34
07 Room Air	0.000	12:35

Lot No = 05620080A1

Cyl No = 47

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature
CHARLES EDER

ALIBRATION CHECK

Remarks:

0.080 Ac Form 106-18000

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
05/12/2020

SN 80-007095 8164.16.00 09/18 12:35

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:35
02 Std. Gas	0.083	12:36
03 Room Air	0.000	12:36
04 Std. Gas	0.083	12:37
05 Room Air	0.000	12:37
06 Std. Gas	0.082	12:38
07 Room Air	0.000	12:38

Lot No = 05620080A1

Cyl No = 47

Exp Date = 04/05/2022

County = 08

Remarks:

Oper No. = 666666

Operator Signature
CHARLES EDER

ALIBRATION CHECK

 $\sigma$