

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
SFN 59281 (06/2018)

erial Number Instrument Location 0-007095 HETT				
Reason for Install/Repair				
	er Location Change			
Other (Specify)				
Check When Done:				
☑ 2. Telephone Line Connected to Intoxilyzer® 8000.				
ヌ 3. Breath Tube Heated.				
X 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
∑ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
区 6. Run Tests:				
■ A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
☑ C. Radio Frequency Interference (RFI) Test (CMS Mode or provided to the control of the co	r Level 1, Function B or C; Key Radio Du	ring Test).		
7. Repair and/or Maintenance Performed (if any):				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.				
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.				
Field Inspector Signature 4 Collins		Date 5/13/2020		
Crime Laboratory Use Only				

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	Charla Exten	Certified Date 5/14/2020

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007095
Location = HETT 8164.16.00 09/18
05/13/2020 16:41

******* Printer Test *******

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz12345f7890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 24119080A1
Standard Cyl #? 25
Standard Expiration? 11/05/2021
Oper No? 020533

Flow Cal. Date: 09/27/2018 Slope 683 Intercept -637452

IR Calibration Date: 09/27/2018

3um
9um

Oth Coef(*100): -59629 -23097
1st Coef(*100): 252379 137307
2nd Coef(*100): 2489 1549
H2O adj(mg/l*10k): 1414 594

******* Printer Test End ********

Operator Signature TRAVIS COLLINS

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = HETT 05/13/2020

SN 80-007095 8164.16.00 09/18 17:02

DRY CAL CHECK

T	est		AC	Time
	Room		0.000	17:03
	Std. Room		0.081 0.000	17:03 17:04
	Std. Room		0.081 0.000	17:04 17:05
06	Std.	Gas	0.080	17:05
07	Room	Air	0.000	17:06

Lot No = 24119080A1

Cyl No = 25

Exp Date = 11/05/2021

County = 01

Oper No. = 020533

Operator Signature TRAVIS COLLINS

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-007095
Location = HETT	8164.16.00 09/18
05/13/2020	16:59

Test	AC	Time
01 Room Air	0.000	17:01
02 *Subject Test	RFI*	17:02
03 Room Air	0.000	17:02

*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST T
Sub DOB = 09/09/1987
Sub Sex = Female
Test = DUI

Weight = 200 Cit = 1111111

Dr. Lic. = ND/1234567 Lot No = 24119080A1

Cyl No = 25

Expiration Date = 11/05/2021

County = 01

Oper No. = 020533

Operator Signature TRAVIS COLLINS

Remarks:

Form 106-I8000