

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

80-007093	Instrument Location BOPD
eason for Install/Repair	
X Install After Receiving From Crime Laboratory Install	After Location Change
Other (Specify)	
neck When Done:	
1 Surge Protector Installed/Property Grounded	
🕱 2 Telephone Line Connected to Intoxilyzer® 8000 🕞 🔾	4.
🔀 4 Enter Preliminary Data (i.e. Date, Time, DST (Y), and Loc	ation; Level 2, Function E)
5 Scan/Enter Gas Cylinder Information (Level 1, Function S).
★ B ACA Test (Level 1, Function C)	
C Radio Frequency Interference (RFI) Test (CMS Mode	or Level 1, Function B or C; Key Radio During Test)
7 Repair and/or Maintenance Performed (if any):	nual Inspection
★ 8 Complete the Top Portion of the Intoxilyzer® Record (SFN)	50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
🕱 9 File Previous Intoxilyzer® Record (SFN504096, Form 120-	
10 Send the Following to the Crime Laboratory: Completed In 104-G). Print Test, ACA Test, and RFI Test.	ntoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form
d Inspector Signature	Date
I B	5.15.2020
ne Laboratory Use Only	
	be used for the analysis of breath to determine alcohol concentration file at the Office of Attorney General, Crime Laboratory Division, in the opy of the documents received
riewed/Certified By	Certified Date

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007093
Location = BOPD 8164.16.00 09/18
05/15/2020 20:52

********* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Enabled Data Entry Mode: Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No No Display Volume? Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry 0.080 Standard Value? Standard Lot #? Standard Cyl #? 24119080A1 33 Standard Expiration? 11/05/2021 130538 Oper No?

Flow Cal. Date: 06/07/2019 Slope 638 Intercept -611758

IR Calibration Date: 12/16/2019
3um 9um

Oth Coef(*100): -14610 -32356
1st Coef(*100): 264615 138933
2nd Coef(*100): 2992 1364
H2O adj(mg/l*10k): 369 621

******* Printer Test End ********

Operator Signature
TREVOR BERGERSON

Remarks:

Install

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007093 Location = BOPD 05/15/2020

8164.16.00 09/18 20:53

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:54
02 Std. Gas	0.078	20:55
03 Room Air	0.000	20:55
04 Std. Gas	0.078	20:55
05 Room Air	0.000	20:56
06 Std. Gas	0.077	20:56
07 Room Air	0.000	20:57

Lot No = 24119080A1

Cyl No = 33

Exp Date = 11/05/2021

County = 06

Oper No. = 130538

Operato Signature

TREVOR BERGERSON

Remarks:

Install

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BOPD
SN 80-007093
8164.16.00 09/18
05/15/2020
20:59

Test	AC	Time
01 Room Air	0.000	21:03
02 *Subject Test	RFI*	21:03
03 Room Air	0.000	21:04

*Invalid Test Inhibited - RFI

Sub Name = RFI TEST, RFI TEST RFI TEST

Sub DOB = 06/03/1988

Sub Sex = Male Weight = 155
Test = OTH Cit = NA

Dr. Lic. = ND/NA Lot No = 24119080A1

Cyl No = 33

Expiration Date = 11/05/2021

County = 06 Oper No. = 130538

Operator Signature TREVOR BERGERSON

Remarks:

Install

Form 106-I8000