Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

Intoxilyzer[®] 8000 Serial Number: 80-007094

Location: ToxL

- Α. General Setup and Checks:
 - Z Diagnostic Tests Pass and Instrument in Ready Mode 1.
 - 2. Breath tube heated
 - Z Date, time and location code (Level 2,E). Re-set if necessary. 3.
 - X Print test (Level 1,P). Sign and attach test record. 4.
 - X Tank monitor (Level 3,D,G). Display and Regulator within 50 psi. 5.
 - Display: 684 psi Regulator: 700 psi
- Tests (Sign and attach test records): Β.

X Configure simulator for the following tests (Level 1,S).

- Wet Calibration Check
 - \mathbf{X} Low AC (Level 1,C): Known Value ≤ 0.03 AC: $\mathbf{O}_{\mathbf{i}} \mathbf{O}_{\mathbf{i}} \mathbf{O}_{\mathbf{i}} \mathbf{O}_{\mathbf{i}}$ AC а. Sim. Ser #: MP 306 [Lot #: 20/810 D Exp. Date: 10.24, 20
 - X High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC b. Sim. Ser #: MP3069 Lot #: 201902D Exp. Date: 2.20.21

Configure dry gas standard for the following tests (Level 1,S).

- XInterferent Check (Level 1,B): Known Value: O- 10 AC + 0.05% Acetone 2. Sim. Ser #: <u>DR3451</u> Lot #: <u>TC96</u> Exp. Date: <u>NA</u> Display reads "Interferent Detect"
- Display reads "RFI Detect" KRFI Check (CMS Mode) 3.
- Dry Calibration Check (Level 1,C): Known Value 0.080 AC 4. $\begin{array}{c} Gas \ Cyl \ Lot \ \#: \ \underline{05620080Al} \ Cyl \ \#: \ \underline{34} \ Exp. \ Date: \ \underline{4.5.22} \\ Test \ 1 \ \underline{0.080} \ Test \ 4 \ \underline{0.080} \ Test \ 4 \ \underline{0.080} \ Test \ 7 \ \underline{0.080} \ Test \ 7 \ \underline{0.080} \ Test \ 8 \ \underline{0.080} \ Test \ 9 \ \underline{0.081} \ \end{array}$
 - Average 0.080
- Remarks/Maintenance: ANNUAL INSPECTION, INSTRUMENT LEFT ON MOUNTAIN TIME FOR FUSPECTION. C.

Instrument is acceptable to be used in the field.

Breath Analyst Signature

Reviewed by

02/18 Issuing Authority - CEE

Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-007094 Location = CARS CEE 8164.16.00 09/18 TOYL 06/08/2020 10:22 abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 05620080A1 Standard Cyl #? 34 Standard Expiration? 04/05/2022 Oper No? 666666 Flow Cal. Date: 09/26/2018 Slope 672 Intercept -607823 IR Calibration Date: 09/26/2018 3um 9um 0th Coef(*100): -28849 -15361 1st Coef(*100):2645712nd Coef(*100):3025 134954 1343 H2O adj(mq/l*10k): 851379

Operator Signature CHARLES EDER

Remarks:

PRINT TEST

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-007094Location = TOXL8164.16.00 09/1806/08/202010:25

	WET CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air 04 Std. Sol. 05 Room Air 06 Std. Sol.	0.000 0.018 0.000 0.018 0.000 0.018	10:25 10:26 10:27 10:27 10:28 10:29
07 Room Air	0.000	10:29

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = MP3061 Std Sol No = 201810D County = 08

Oper No. = 666666

rai

Operator Signature CHARLES EDER

Remarks:

Low AC 0.020AC

Form 106-18000

۱

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-007094Location = TOXL8164.16.00 09/1806/08/202010:30

	WET CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air	0.000 0.297 0.000	10:31 10:31 10:32
04 Std. Sol. 05 Room Air 06 Std. Sol. 07 Room Air	0.298 0.000 0.297 0.000	10:33 10:33 10:34 10:35

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

an

Operator Signature CHARLES EDER

HIGH AC

0,300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-007094 Location = TOXL 8164.16.00 09/18 06/08/2020 10:36 Test AC Time 01 Room Air 0.000 10:37 02 *Subject Test INT* 10:37 03 Room Air 0.000 10:38 *Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = INTERFERENT Dr. Lic. = ND/DIS921456Lot No = 05620080A1Cyl No = 34Expiration Date = 04/05/2022County = 08Oper No. = 666666

Operator Signature CHARLES EDER Remarks: Interforment Check 0,104C + 0.05%. ACETONE

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-007094 Location = TOXL 8164.16.00 09/18 06/08/2020 10:41 Test AC Time 01 Diagnostic OK 10:42 02 Room Air 0.000 10:43 03 *Subject Test RFI* 10:43 04 Room Air 0.000 10:43 *Invalid Test Inhibited - RFI Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUICit = RFI CHECK Dr. Lic. = ND/DIS921456Lot No = 05620080A1Cyl No = 34Expiration Date = 04/05/2022County = 08Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Form 106-18000

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-007094 Location = TOXL8164.16.00 09/18 06/08/2020 10:48

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.080 0.000 0.080 0.000	10:49 10:49 10:50 10:50 10:51 10:51 10:52

Lot No = 05620080A1Cyl No = 34Exp Date = 04/05/2022County = 08

Oper No. = 666666

tai

Operator Signature CHARLES EDER

ALIBRATION CHECK

0.080 AC Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-007094Location = TOXL8164.16.00 09/1806/08/202010:52

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.080 0.000 0.080 0.080	10:52 10:53 10:53 10:54 10:54 10:55 10:55

Lot No = 05620080A1 Cyl No = 34 Exp Date = 04/05/2022 County = 08

Oper No. = 666666

lu

Operator Signature CHARLES EDER

ALIBRATION CHECK

0.080 AC

Form 106-18000

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-007094 Location = TOXL 8164.16.00 09/18 06/08/2020 10:55

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.080 0.000 0.081 0.000	10:56 10:56 10:57 10:57 10:58 10:58 10:59

Lot No = 05620080A1Cyl No = 34Exp Date = 04/05/2022County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

ALIBRATION CHECK

0.080AC

Form 106-18000