

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 007093	Instrument Location UNDP	
Reason for Install/Repair		
Install After Receiving From Crime Laboratory 🔲 Ins	stall After Location Change	
Other (Specify)		
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
Telephone Line Connected to Intoxilyzer® 8000.		
☐ 3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and	Location; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Functi	on S).	
6. Run Tests:		
A. Print Test (Level1, Function P).		
B ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS M	lode or Level 1, Function B or C; Key Ra	adio During Test).
7. Repair and/or Maintenance Performed (if any):	o repair or mainten	ance performed
▼ 8. Complete the Top Portion of the Intoxilyzer® Record	(SFN50496, Form 120-G) and Place it b	by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form	120-G) at the Intoxilyzer® Location at t	he Agency.
<ol> <li>Send the Following to the Crime Laboratory: Comple 104-G), Print Test, ACA Test, and RFI Test.</li> </ol>	ted Intoxilyzer® 8000 Installation and Re	epair Checkout (SFN59281, Form
ield Inspector Signature		Date 06-09-2020
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approve the date the Field Inspector performed the installation. This recor County of Burleigh, North Dakota, is certified to be a true and corr	d on file at the Office of Attorney General	
Reviewed/Certified By	El	Certified Date 6 • 12 • 2020

FORM 104-G

Current Instrument Setup Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No No Display Volume? Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? Standard Lot #? Standard Cyl #? 0.080 24119080A1 Standard Expiration? 11/05/2021 119409 Oper No?

Flow Cal. Date: 09/26/2018 Slope 670 Intercept -596123

IR Calibration Date: 09/26/2018
3um 9um
----0th Coef(\*100): -27453 -17291

 0th Coef(\*100):
 -27453
 -17291

 1st Coef(\*100):
 248629
 134601

 2nd Coef(\*100):
 2770
 1393

 H2O adj(mg/l\*10k):
 872
 414

\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature TRACY MEIDINGER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-007092 Location = UNDP 06/09/2020

8164.16.00 09/18 10:07

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:08
02 Std. Gas	0.080	10:08 10:09
03 Room Air	0.000	10:09
04 Std. Gas	0.080 0.000	10:10
05 Room Air	0.000	10:10
06 Std. Gas	0.000	10:11
07 Room Λir	0.000	10.11

Lot No = 24119080A1

Cyl No = 3

Exp Date = 11/05/2021

County = 18

Oper No. = 119409

Operator Signature TRACY MEIDINGER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = UNDP
06/09/2020

Alcohol Analyzer
SN 80-007092
8164.16.00 09/18

Test	AC	Time
01 Room Air	0.000	10:13
02 *Subject Test	RFI*	10:14
03 Room Air	0.000	10:14

\*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST TEST

Sub DOB = 09/09/1990

Sub Sex = Male Weight = 150Test = DUI Cit = 88888888

Dr. Lic. = ND/55555555

Lot No = 24119080A1

Cyl No = 3

Expiration Date = 11/05/2021

County = 18 Oper No. = 119409

perator Signature

Remarks:

Form 106-I8000