

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROF

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number				
	surty So Or			
Reason for Install/Repair	<u> </u>			
Install After Receiving From Crime Laboratory 🔲 Install After Location Change				
Other (Specify) Took to crime Lab for Annua	al Inspection			
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function B	Ε).			
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any):				
A. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.				
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Lo	cation at the Agency.			
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installar 104-G), Print Test, ACA Test, and RFI Test.	ion and Repair Checkout (SFN59281, Form			
Field Inspector Signature	Date 06-30-2020			

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	Cha	latth	Certified Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-007094Location = CARS8164.16.00 09/1806/30/202021:32

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Inhib Printer (Y/N)? Display Volume? Disable On Memfull? # of Print Copies? Select Std (D/W/I)? Standard Value? Standard Lot #? Standard Cyl #?	Enabled DABACABA Yes Yes No No Yes 1 Dry 0.080 05620080A1 37 04/05/2022 132318 09/26/2018 672 -607823 09/26/2018 9um
0th Coef(*100): -28849	-15361
1st Coef(*100):2645712nd Coef(*100):3025H2O adj(mg/l*10k):851	L 134954 1343 379

********** Printer Test End **********

Operator Signature TERRY DILLON

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-007094Location = CARS8164.16.00 09/1806/30/202021:32

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.078 0.000 0.078 0.000 0.078 0.000	21:34 21:34 21:35 21:35 21:36 21:36 21:37

Lot No = 05620080A1 Cyl No = 37 Exp Date = 04/05/2022 County = 19

Oper No. = 132318

en

Operator Signature TERRY DILLON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-007094 Location = CARS 8164.16.00 09/18 06/30/2020 21:41 Test AC Time 01 Room Air 0.000 21:46 02 *Subject Test 21:46 RFI* 03 Room Air 21:46 0.000 *Invalid Test Inhibited - RFI Sub Name = RFITEST, INSTALL NA Sub DOB = 01/01/1900Sub Sex = Unknown Weight = NA Cit = NA Test = DUI Dr. Lic. = ND/RFI001Lot No = 05620080A1Cyl No = 37Expiration Date = 04/05/2022County = 19Oper No. = 132318

en

Operator Signature TERRY DILLON

Remarks:

Form 106-18000