

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-007097	Instrument Location Rolette County Sheriff Department				
Reason for Install/Repair					
☑ Install After Receiving From Crime Laboratory ☐ Install After					
Other (Specify)					
Check When Done:					
■ 1. Surge Protector Installed/Property Grounded.					
▼ 2. Telephone Line Connected to Intoxilyzer® 8000.					
🕱 3. Breath Tube Heated.					
▼ 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
▼ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).					
🕱 6. Run Tests:					
X A. Print Test (Level1, Function P).					
B. ACA Test (Level 1, Function C).					
☑ C. Radio Frequency Interference (RFI) Test (CMS Mode or)	Level 1, Function B or C; Key Radio During Test).				
X 7. Repair and/or Maintenance Performed (if any): N/A					
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50	0496, Form 120-G) and Place it by the Intoxilyzer® for Use.				
▼ 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.					
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.					
Field Inspector Signature HeBJLZY	Date 9/25/2020				
Crime Laboratory Use Only					
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on fill County of Burleigh, North Dakota, is certified to be a true and correct county.	le at the Office of Attorney General, Crime Laboratory Division, in the				

Reviewed/Certified By	Charles Este	Certified Date 10.2-20

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007097
Location = ROLL 8164.16.00 09/18
09/17/2020 13:45

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\%*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry 0.080 Standard Value? 05620080A1 Standard Lot #? Standard Cyl #? 04/05/2022 Standard Expiration? Oper No? 131679

Flow Cal. Date: 09/03/2020 Slope 660 Intercept -636264

 0th Coef(*100):
 -36810
 -25919

 1st Coef(*100):
 254044
 139068

 2nd Coef(*100):
 1858
 1258

 H2O adj(mg/l*10k):
 903
 584

****** Printer Test End *******

Operator Signature ALEX BREITBACH

Remarks:

Print test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = ROLL
09/17/2020

Alcohol Analyzer SN 80-007097 8164.16.00 09/18 14:09

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:10
02 Std. Gas	0.081	14:10
03 Room Air	0.000	14:11
04 Std. Gas	0.081	14:11
05 Room Air	0.000	14:12
06 Std. Gas	0.081	14:12
07 Room Air	0.000	14:13

Lot No = 05620080A1

Cyl No = 10

Exp Date = 04/05/2022

County = 40

Oper No. = 131679

Operator Signature · ALEX BREITBACH

Remarks:

ACA TEST FORM 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-007097
Location = ROLL	8164.16.00 09/18
09/17/2020	14:14

Test	AC	Time
01 Room Air	0.000	14:17
02 Reference	RFI*	14:17
03 Room Air	0.000	14:17

^{*}Invalid Test Inhibited - RFI

Sub Name = RFI, TEST INSTALL

Sub DOB = 01/01/1969

Sub Sex = Unknown

Weight = NA

Test = OTH

Cit = NA

Dr. Lic. = ND/NA

Lot No = 05620080A1

Cyl No = 10

Expiration Date = 04/05/2022

County = 40

Oper No. = 131679

Operator Signature ALEX BREITBACH

Remarks:

Form 106-I8000