

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-007097	Instrument Location Rolette County Sheriff Department			
Reason for Install/Repair				
▼ Install After Receiving From Crime Laboratory   □	Install After Location Change			
Other (Specify)				
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
▼ 2. Telephone Line Connected to Intoxilyzer® 8000.				
ヌ 3. Breath Tube Heated.				
▼ 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
▼ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
★ 6. Run Tests:				
B. ACA Test (Level 1, Function C).				
X C. Radio Frequency Interference (RFI) Test (CM	IS Mode or Level 1, Function B or C; Key Radio During Test).			
🗵 7. Repair and/or Maintenance Performed (if any):	None			
▼ 8. Complete the Top Portion of the Intoxilyzer® Record	cord (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.			
▼ 9. File Previous Intoxilyzer® Record (SFN504096, F	Form 120-G) at the Intoxilyzer® Location at the Agency.			
10. Send the Following to the Crime Laboratory: Con 104-G), Print Test, ACA Test, and RFI Test.	mpleted Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form			
Field Inspector Signature	Date			
Jeffrey Bauske	4/16/2021			
Crime Laboratory Use Only				
	proved to be used for the analysis of breath to determine alcohol concentration from			

FORM 104-G

Reviewed/Certified By

Certified Date
4.26.21

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007097
Location = ROLL 8164.16.00 09/18
04/16/2021 12:19

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? Disable On Memfull? Yes 1 # of Print Copies? Select Std (D/W/I)? Dry 0.080 Standard Value? Standard Lot #? 05620080A1 Standard Cyl #? 10 04/05/2022 Standard Expiration? 120501 Oper No?

Flow Cal. Date: 09/03/2020 Slope 660 Intercept -636264

IR Calibration Date: 09/03/2020

3um 9um

Oth Coef(\*100): -36810 -25919
1st Coef(\*100): 254044 139068
2nd Coef(\*100): 1858 1258

2nd Coef(\*100): 1858 1258 H2O adj(mg/l\*10k): 903 584

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

perator Signature JEFFREY BAUSKE

Remarks: PRINT TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = ROLL
04/16/2021

Alcohol Analyzer SN 80-007097 8164.16.00 09/18 12:23

## DRY CAL CHECK

Test	AC ′	Time
01 Room Air	0.000	12:23
02 Std. Gas	0.081	12:24
03 Room Air	0.000	12:24
04 Std. Gas	0.081	12:25
05 Room Air	0.000	12:25
06 Std. Gas	0.082	12:25
07 Room Air	0.000	12:26

Lot No = 05620080A1

Cyl No = 10

Exp Date = 04/05/2022

County = 40

Oper No. = 120501

Operator Signature JEFFREY BAUSKE

Remarks: ACA TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = ROLL 04/16/2021 SN 80-007097 8164.16.00 09/18 12:20

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:20
02 Std. Gas	0.082	12:21
03 Room Air	0.000	12:21
04 Std. Gas	0.082	12:22
05 Room Air	RFI*	12:22
06 Room Air	0.000	12:22

\*Invalid Test Inhibited - RFI

Lot No = 05620080A1

Cyl No = 10

Exp Date = 04/05/2022

County = 40

Oper No. = 120501

Operator Signature JEFFREY BAUSKE

Remarks: **RFI TEST** 

Form 106-I8000