

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-007095	Instrument Location HETT		
Reason for Install/Repair			
▼ Install After Receiving From Crime Laboratory	☐ Install After Location Change		
Other (Specify)			
Check When Done:			
∑ 2. Telephone Line Connected to Intoxilyzer® 80	00.		
☑ 3. Breath Tube Heated.			
✓ 4. Enter Preliminary Data (I.e. Date, Time, DST)	(Y), and Location; Level 2, Function E).		
	, Function S).		
X A. Print Test (Level1, Function P).			
⋈ B. ACA Test (Level 1, Function C).			
	(CMS Mode or Level 1, Function B or C; Key Radio During Test).		
7. Repair and/or Maintenance Performed (if any	):		
	Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.		
☑ 9. File Previous Intoxilyzer® Record (SFN5040)	96, Form 120-G) at the Intoxilyzer® Location at the Agency.		
10. Send the Following to the Crime Laboratory: 104-G), Print Test, ACA Test, and RFI Test.	Completed intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form		
Field Inspector Signature	Date 5/6/2021		
Crime Laboratory Use Only			
•	consequent to be used for the analysis of breath to determine alcohol concentration from		

the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the

County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

FORM 104-G

Reviewed/Certified By

Certified Date 5: 7.2021

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007095
Location = HETT 8164.16.00 09/18
05/06/2021 17:36

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\^&\*()\_+?

Current Instrument Setup Enabled Data Entry Mode: Start Test Sequence: DABACABA
Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No Yes , No Display Volume? Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 24119080A1 Standard Cyl #? 25 Standard Expiration? 11/05/2021 Oper No? 020533

Flow Cal. Date: 09/27/2018 Slope 683 Intercept -637452

IR Calibration Date: 09/27/2018
3um 9um

Oth Coef(\*100): -59629 -23097
1st Coef(\*100): 252379 137307
2nd Coef(\*100): 2489 1549
H2O adj(mg/l\*10k): 1414 594

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*

Operator Signature TRAVIS COLLINS

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = HETT 05/06/2021

Alcohol Analyzer SN 80-007095 8154.16.00 09/18 18:02

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:03
02 Std. Gas	0.082	18:03
03 Room Air	0.000	18:04
04 Std. Gas	0.082	18:04
05 Room Air	0.000	18:04
06 Std. Gas	0.081	18:05
07 Room Air	0.000	18:05

Lot No = 24119080A1

Cyl No = 25

Exp Date = 11/05/2021

County = 01

Oper No. = 020533

Operator Signature TRAVIS COLLINS

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

SN 80-007095

18:06

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = HETT 8164.16.00 09/18 05/06/2021

Test	AC	Time
01 Room Air	0.000	18:07
02 *Subject 3	Cest RFI*	18:07
03 Room Air	0.000	18:08

\*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST T Sub DOB = 06/23/1987

Sub Sex = FemaleWeight = 119Test = MICCit = 1234

Dr. Lic. = ND/98765Lot No = 24119080A1

Cyl No = 25

Expiration Date = 11/05/2021

County = 01Oper No. = 020533

> Operator Signature TRAVIS COLLINS

Remarks:

Form 106-I8000