

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 007092	Instrument Location UND-PD
Reason for Install/Repair	
Install After Receiving From Crime Laboratory	er Location Change
Other (Specify)	
Check When Done:	
1. Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
Breath Tube Heated.	
Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locat	ion; Level 2, Function E).
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
6. Run Tests:	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
Radio Frequency Interference (RFI) Test (CMS Mode of	⁻ Level 1, Function B or C; Key Radio During Test).
T. Repair and/or Maintenance Performed (if any):	
S. Complete the Top Portion of the Intoxilyzer® Record (SFN5)	
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-0	3) at the Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: Completed Int 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form
ield Inenester Signature	
ield Inspector Signeture	Date 06-04-2021

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	Cultor Ela	Certified Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-007092 Location = UNDP 8164.16.00 09/18 06/04/2021 10:32 abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 02621080A1 17 Standard Expiration? 03/05/2023 Oper No? 119409 Flow Cal. Date: 09/26/2018 Slope 670 Intercept -596123 IR Calibration Date: 09/26/2018 3um 9um _____ Oth Coef(*100): -27453 -17291
 1st Coef(*100):
 248629

 2nd Coef(*100):
 2770
134601 1393 H2O adj(mg/l*10k): 872 414

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000Alcohol Analyzer
SN 80-007092Location = UNDP8164.16.00 09/18
10:33

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.079 0.000 0.080 0.000 0.080 0.080 0.000	10:33 10:34 10:34 10:35 10:35 10:36 10:36

Lot No = 02621080A1 Cyl No = 17 Exp Date = 03/05/2023 County = 18 Ope

Oper No. = 119409

Operator Signature

TRACY MEIDINGER

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-007092 Location = UNDP 8164.16.00 09/18 06/04/2021 10:37 Test AC Time 01 Room Air 0.000 10:38 02 *Subject Test RFI* 10:39 03 Room Air 0.000 10:40 *Invalid Test Inhibited - RFI Sub Name = TEST, S TEST Sub DOB = 01/01/1999Sub Sex = Male Weight = 150Test = DUI Cit = 77777 Dr. Lic. = ND/454545Lot No = 02621080A1Cyl No = 17Expiration Date = 03/05/2023County = 18Oper No. = 119409

Operator Signature TRACY MEIDINGER

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000Alcohol Analyzer
SN 80-007092Location = UNDP8164.16.00 09/18
10:22

DIAGNOSTIC

Voltage/Current Test RAM Test	Pass Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 119409

Operator Signature TRACY MEIDINGER

Remarks:

Form 106-18000

4