



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <b>80-007097</b>	Instrument Location <b>Rolla S.O</b>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature 	Date <b>7-5-22</b>
-------------------------------	-----------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By 	Certified Date <b>07/11/2022</b>
---------------------------	-------------------------------------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-007097  
Location = ROLL                8164.16.00 09/18  
07/05/2022                      10:35

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            34921080A1  
Standard Cyl #?            61  
Standard Expiration?       02/05/2024  
Oper No?                    130918

Flow Cal. Date:            09/03/2020  
Slope                      660  
Intercept                  -636264

IR Calibration Date:       09/03/2020  
                                    3um                      9um

-----  
0th Coef(\*100):            -36810                      -25919  
1st Coef(\*100):            254044                      139068  
2nd Coef(\*100):            1858                        1258  
H2O adj(mg/l\*10k):        903                         584

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
MITCH SLATER

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-007097  
Location = ROLL                    8164.16.00 09/18  
07/05/2022                            10:41

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:42
02 Std. Gas	0.079	10:42
03 Room Air	0.000	10:43
04 Std. Gas	0.079	10:43
05 Room Air	0.000	10:43
06 Std. Gas	0.080	10:44
07 Room Air	0.000	10:44

Lot No = 34921080A1  
Cyl No = 61  
Exp Date = 02/05/2024  
County = 40                            Oper No. = 130918



Operator Signature  
MITCH SLATER

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-007097  
Location = ROLL      8164.16.00 09/18  
07/05/2022      10:45

Test	AC	Time
01 Room Air	RFI*	10:46
02 Room Air	RFI*	10:46

\*Invalid Test  
Inhibited - RFI

Sub Name = RFI TEST, RFI TEST RFI TEST  
Sub DOB = 04/24/1993  
Sub Sex = Male      Weight = 234  
Test = OTH      Cit = 123456  
Dr. Lic. = ND/SLA934567  
Lot No = 34921080A1  
Cyl No = 61  
Expiration Date = 02/05/2024  
County = 40      Oper No. = 130918



Operator Signature  
MITCH SLATER

Remarks: