



# INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION  
Toxicology Section/Breath Alcohol Program  
SFN 59281 (5-2011)

|                                                                                                                                                                                                                                                        |                                    |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|
| Serial Number<br><b>80-007091</b>                                                                                                                                                                                                                      | Instrument Location<br><b>Mott</b> |
| Reason for Install/Repair<br><input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change<br><input checked="" type="checkbox"/> Other (Specify) <u>Install After Yearly Check-Up</u> |                                    |

### Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): N/A
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

|                                                                 |                          |
|-----------------------------------------------------------------|--------------------------|
| Field Inspector Signature<br><i>J. Allen</i>                    | Date<br><b>7/1/20</b>    |
| Reviewed By (Crime Laboratory Use Only)<br><i>Charles E. Ed</i> | Date<br><b>7.15.2020</b> |

### Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By: *Charles E. Ed* Certified Date: **7.15.2020**

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-007091  
Location = MOTT                8164.16.00 09/18  
07/01/2020                      14:04

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            13518080A6  
Standard Cyl #?            18  
Standard Expiration?       08/05/2020  
Oper No?                    132666

Flow Cal. Date:            09/25/2018  
Slope                      681  
Intercept                  -651030

IR Calibration Date:       09/25/2018  
                                    3um                      9um

|                    | 3um    | 9um    |
|--------------------|--------|--------|
| 0th Coef(*100):    | -39121 | -20906 |
| 1st Coef(*100):    | 240667 | 132345 |
| 2nd Coef(*100):    | 3960   | 1516   |
| H2O adj(mg/l*10k): | 899    | 428    |

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
Operator Signature  
JEDEDIAH KOHLER

Remarks: *Print test*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-007091  
Location = MOTT      8164.16.00 09/18  
07/01/2020      14:04

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 14:05 |
| 02 Std. Gas | 0.078 | 14:05 |
| 03 Room Air | 0.000 | 14:06 |
| 04 Std. Gas | 0.077 | 14:06 |
| 05 Room Air | 0.000 | 14:07 |
| 06 Std. Gas | 0.077 | 14:07 |
| 07 Room Air | 0.000 | 14:08 |

Lot No = 13518080A6  
Cyl No = 18  
Exp Date = 08/05/2020  
County = 21      Oper No. = 132666

  
Operator Signature  
JEDEDIAH KOHLER

Remarks: *ACA Function test*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-007091  
Location = MOTT                 8164.16.00 09/18  
07/01/2020                        14:13

| Test          | AC   | Time  |
|---------------|------|-------|
| 01 Diagnostic | OK   | 14:15 |
| 02 Room Air   | RFI* | 14:15 |
| 03 Room Air   | RFI* | 14:15 |

\*Invalid Test  
  Inhibited - RFI

Sub Name = KOHLER, JEDEDIAH SCOTT  
Sub DOB = 08/20/1987  
Sub Sex = Male                     Weight = 175  
Test = DUI                         Cit = NA  
Dr. Lic. = ND/KOH872036  
Lot No = 13518080A6  
Cyl No = 18  
Expiration Date = 08/05/2020  
County = 21                        Oper No. = 132666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

  
\_\_\_\_\_  
Operator Signature  
JEDEDIAH KOHLER

Remarks: *RFI test*

Form 106-I8000