

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-006510

Location: TOXL

Check When Done

- A. Is the breath tube heated? ✓
- B. Is the diagnostic check complete? ✓
- C. Is the time, date, and year correct (re-set if necessary)? ✓
- D. Print test. (Attach test record.) ✓
- E. Low AC. Use  $\leq 0.03$  AC in ACA mode. (Attach test record.)  
Sim SN DRS114 Lot No. 201505B Rep AC 0.030 ✓
- F. High AC. Use  $\geq 0.25$  AC in ACA mode. (Attach test record.)  
Sim SN MP3068 Lot No. 14102 Rep AC 0.400 ✓
- G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.  
(Attach test record.)  
Sim SN DRS135 Lot No. ICS4 Solution 0.10AC+0.05%ACE ✓
- H. RFI check. Run CMS mode. Key radio during CMS test.  
(Attach test record.) ✓
- I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.  
(Attach test record.) 0.080 AC

Lot No. 13314080A3 Cyl No. 028 Exp. Date 6.5.2016

Test 1 <u>0.080</u>	Test 4 <u>0.081</u>	Test 7 <u>0.081</u>
Test 2 <u>0.080</u>	Test 5 <u>0.081</u>	Test 8 <u>0.080</u>
Test 3 <u>0.080</u>	Test 6 <u>0.081</u>	Test 9 <u>0.081</u>

Average 0.081

J. Tank pressure, Level 3 Func D Sub G P: 832 psi Reg ~800 psi

K. Remarks/Maintenance record: \_\_\_\_\_

4.1.2016  
Date  
01 Apr 16  
Date

[Signature]  
Field Inspector's Signature  
[Signature]  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer          Alcohol Analyzer  
North Dakota Model 8000          SN 80-006510  
Location = TOXL                  8164.13.00 06/09  
04/01/2016                                  12:17

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\* ()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\* ()\_+?

Current Instrument Setup


Data Entry Mode:                  Enabled  
Start Test Sequence:              DABACABA  
Display Prelim Rslt?              Yes  
Display Third Digit?              Yes  
Inhib Printer(Y/N)?              No  
Display Volume?                  No  
Disable On Memfull?              Yes  
# of Print Copies?                1  
Select Std (D/W/I)?              Dry  
Standard Value?                  0.080  
Standard Lot #?                  13314080A3  
Standard Cyl #?                  28  
Standard Expiration?              06/05/2016  
Oper No?                          666666

Flow Cal. Date:                  08/24/2015  
Slope                              624  
Intercept                         -539297

IR Calibration Date:              08/24/2015  
                                      3um                              9um

-----  
0th Coef(\*100):                  -23613                              -13697  
1st Coef(\*100):                  251959                              131614  
2nd Coef(\*100):                  4659                                1853  
H2O adj(mg/l\*10k):              819                                 493

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006510  
Location = TOXL      8164.13.00 06/09  
04/01/2016      12:18

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:18
02 Std. Sol.	0.028	12:19
03 Room Air	0.000	12:20
04 Std. Sol.	0.028	12:20
05 Room Air	0.000	12:21
06 Std. Sol.	0.029	12:22
07 Room Air	0.000	12:22

08 Sim Temp = 34.0°C

Simul Ser No = DR5114  
Std Sol No = 201505B  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Low AC*  
*0.030 AC*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006510  
Location = TOXL      8164.13.00 06/09  
04/01/2016      12:24

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:24
02 Std. Sol.	0.404	12:25
03 Room Air	0.000	12:26
04 Std. Sol.	0.404	12:27
05 Room Air	0.000	12:27
06 Std. Sol.	0.404	12:28
07 Room Air	0.000	12:29

08 Sim Temp = 34.0°C

Simul Ser No = MP3068  
Std Sol No = 14102  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *HIGANCE HIGH AC*  
*0.400 AC*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006510  
Location = TOXL      8164.13.00 06/09  
04/01/2016      12:30

Test	AC	Time
01 Room Air	0.000	12:31
02 *Subject Test	INT*	12:31
03 Room Air	0.000	12:32

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 13314080A3  
Cyl No = 28  
Expiration Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *INTERFERENT CHECK*  
*0.10AC + 0.05% ACETONE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006510  
Location = TOXL      8164.13.00 06/09  
04/01/2016      12:33

Test	AC	Time
01 Diagnostic	OK	12:34
02 Room Air	0.000	12:34
03 Reference	RFI*	12:34
04 Room Air	0.000	12:35

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 13314080A3  
Cyl No = 28  
Expiration Date = 06/05/2016  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *RFI CHECK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006510  
Location = TOXL      8164.13.00 06/09  
04/01/2016      12:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:36
02 Std. Gas	0.080	12:36
03 Room Air	0.000	12:37
04 Std. Gas	0.080	12:37
05 Room Air	0.000	12:38
06 Std. Gas	0.080	12:38
07 Room Air	0.000	12:38

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006510  
Location = TOXL      8164.13.00 06/09  
04/01/2016      12:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:39
02 Std. Gas	0.081	12:40
03 Room Air	0.000	12:40
04 Std. Gas	0.081	12:41
05 Room Air	0.000	12:41
06 Std. Gas	0.081	12:42
07 Room Air	0.000	12:42

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080AC



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006510  
Location = TOXL      8164.13.00 06/09  
04/01/2016      12:43

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:43
02 Std. Gas	0.081	12:43
03 Room Air	0.000	12:44
04 Std. Gas	0.080	12:44
05 Room Air	0.000	12:45
06 Std. Gas	0.081	12:45
07 Room Air	0.000	12:46

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC