



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (03-2014)

Serial Number 80-006510	Instrument Location HETT
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): NA
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>J Collins</i>	Date <i>05/13/2016</i> ✓
Reviewed By (Crime Laboratory Use Only) <i>Deb Shanaver</i>	Date <i>20 May 16</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006510
Location = HETT 8164.13.00 06/09
05/13/2016 10:07

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 29614080A2
Standard Cyl #? 34
Standard Expiration? 12/05/2016
Oper No? 020533

Flow Cal. Date: 08/24/2015
Slope 624
Intercept -539297

IR Calibration Date: 08/24/2015
 3um 9um

0th Coef(*100): -23613 -13697
1st Coef(*100): 251959 131614
2nd Coef(*100): 4659 1853
H2O adj(mg/l*10k): 819 493

***** Printer Test End *****



Operator Signature
TRAVIS COLLINS

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006510
Location = HETT 8164.13.00 06/09
05/13/2016 10:08

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:09
02 Std. Gas	0.079	10:09
03 Room Air	0.000	10:09
04 Std. Gas	0.080	10:10
05 Room Air	0.000	10:10
06 Std. Gas	0.080	10:11
07 Room Air	0.000	10:11

Lot No = 29614080A2

Cyl No = 34

Exp Date = 12/05/2016

County = 01

Oper No. = 020533



Operator Signature
TRAVIS COLLINS

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006510
Location = HETT 8164.13.00 06/09
05/13/2016 10:14

Test	AC	Time
01 Room Air	0.000	10:17
02 *Subject Test	RFI*	10:17
03 Room Air	0.000	10:17

*Invalid Test
Inhibited - RFI

Sub Name = TEST, JOE T
Sub DOB = 10/10/2014
Sub Sex = Male Weight = 188
Test = MIC Cit = 012345
Dr. Lic. = ND/012345
Lot No = 29614080A2
Cyl No = 34
Expiration Date = 12/05/2016
County = 01 Oper No. = 020533



Operator Signature
TRAVIS COLLINS

Remarks:

Form 106-I8000