



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (03-2014)

Serial Number 80-006513	Instrument Location Mott
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>Sarah Warner</i>	Date 5/16/2016 ✓
Reviewed By (Crime Laboratory Use Only) <i>Deb Shanaver</i>	Date 26 May 16

This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Certified By: *Charles E. Eln* Certified Date: **11/09/2016**

Intoxilyzer Test Record and Checklist
 NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
 North Dakota Model 8000 SN 80-006513
 Location = MOTT 8164.13.00 06/09
 05/16/2016 12:13

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-|=|
 ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|
 ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
 Start Test Sequence: DABACABA
 Display Prelim Rslt? Yes
 Display Third Digit? Yes
 Inhib Printer(Y/N)? No
 Display Volume? No
 Disable On Memfull? Yes
 # of Print Copies? 1
 Select Std (D/W/I)? Dry
 Standard Value? 0.080
 Standard Lot #? 19615080A5
 Standard Cyl #? 25
 Standard Expiration? 09/05/2017
 Oper No? 100620

 Flow Cal. Date: 08/20/2015
 Slope 680
 Intercept -696037

IR Calibration Date: 08/21/2015

	3um	9um
0th Coef(*100):	-14436	-13610
1st Coef(*100):	265490	136147
2nd Coef(*100):	4422	1710
H2O adj(mg/l*10k):	672	448

***** Printer Test End *****


 Operator Signature
 SARAH WARNER

Remarks:


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006513
Location = MOTT 8164.13.00 06/09
05/16/2016 12:14

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:15
02 Std. Gas	0.078	12:15
03 Room Air	0.000	12:16
04 Std. Gas	0.079	12:16
05 Room Air	0.000	12:17
06 Std. Gas	0.080	12:17
07 Room Air	0.000	12:18

Lot No = 19615080A5
Cyl No = 25
Exp Date = 09/05/2017
County = 21 Oper No. = 100620



Operator Signature
SARAH WARNER

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006513
Location = MOTT 8164.13.00 06/09
05/16/2016 12:19

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:20
02 Std. Gas	0.079	12:20
03 Room Air	RFI*	12:21
04 Room Air	RFI*	12:21

*Invalid Test
Inhibited - RFI

Lot No = 19615080A5
Cyl No = 25
Exp Date = 09/05/2017
County = 21 Oper No. = 100620



Operator Signature
SARAH WARNER

Remarks:

Form 106-I8000