

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-006666

A. GENERAL Setup and Checks:

1. Power On
 - 5 min. countdown during initial warm-up
 - 1 min. countdown to begin test after standby mode
2. Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:
 - Level 1 (Operator) Level 2 (Field Inspector) Level 3 (Lab Personnel)
4. Printer Setup (Level 3,S,P)
5. Forms Load (CMI, Inc., CoBrA Software)

B. LEVEL 1 Setup and Checks:

1. Internal Printer Test and Paper Feed (Level 1,P). Attach record.
2. External Printer Test (Level 1,P). Attach record.
3. Barcode Reader (Level 1,B):
 - Operator ID ND License R key: Refusal. Attach record.

C. LEVEL 2 Setup and Checks:

1. Continuous Air Blank (Level 2,A)
2. Diagnostic (Level 2,D). Attach record.
3. Prelim Data Entry (Level 2,E)
4. Ser Number/Version Info (Level 2,V):
 - Inst Ser#: 80-006666 DSP Ver: 1362.14
 - Core Ver: 8164.14 CPLD Ver: 7

D. LEVEL 3 Setup and Checks:

1. DVM Monitor (Level 3,D,D): 3: 12620 s: 0 c: 47.00
9: 13152 s: 1 B: 45 0
(RFI indicator with keyed radio): 1
2. Function (F) Keys: : F5 – Tone F9 - Fan
 F6 – Green LED F10 – Pump
 F7 – Red LED F11 – Solenoid
3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 952 CF: 1.06
4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 214 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T): C: 47.00 A: 48.0 I: 0.0
B: 45.0 V: 0.0
6. Tank Monitor (Level 3,D,G): P (display): 922 psi P (gauge): ~950 psi
7. General Setup (Level 3,S,G)
8. Configure Start Test (Level 3,S,T)
9. Network (Level 3,S,N)

E. TESTS (attach all test records):

Configure wet bath simulator standards before running tests (Level 1,S):

1. Zero AC (Level 1,C):
Sim ser #: MP3064 Lot #: MILLI QH₂O Exp date: N/A
Known value: 0.000 AC
Results: 0.000 AC 0.000 AC 0.000 AC
2. Low AC (Level 1,C):
Sim ser #: MP3060 Lot #: 201501A Exp date: 01.06.17
Known value: 0.015 AC
Results: 0.015 AC 0.015 AC 0.015 AC
3. High AC (Level 1,C):
Sim ser #: MP3063 Lot #: 16102 Exp date: 03.22.18
Known value: 0.400 AC
Results: 0.393 AC 0.394 AC 0.393 AC

Configure dry gas standard before running tests (Level 1,S):

4. Ethanol + 0.05% Acetone (Level 1,B):
Sim ser #: DR5135 Lot #: ICS4 Exp date: 9.9.17
Known value: 0.10AC + 0.05% ACETONE
 Display reads "Interferent Detect"
5. Acetone (Level 1,B):
Sim ser #: DR5143 Lot #: ACE2 Exp date: 9.29.17
Known value: 0.05% ACETONE
 Display reads "Interferent Detect"
6. Methanol (Level 1,B):
Sim ser #: DR3377 Lot #: MET2 Exp date: 9.29.17
Known value: 0.10% METHANOL
 Display reads "Interferent Detect"
7. Isopropanol (Level 1,B):
Sim ser #: DR3453 Lot #: IPA2 Exp date: 9.29.17
Known value: 0.10% ISOPROPANOL
 Display reads "Interferent Detect"

- 8. Ambient Fail (Level 1,B): Display and test record read correct.
- 9. Invalid Sample (Level 1,B): Display and test record read correct.
- 10. Sequence Aborted (Level 1,B): Display and test record read correct.
- 11. ABA (Level 1,B): Test record reads correct.
- 12. Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- 13. RFI Detected (B,C, or CMS mode): Display and test record read correct.
- 14. Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 15. Deficient Samples (CMS mode): Display and test record read correct.
- 16. Improper Sample (CMS mode): Display and test record read correct.
- 17. Difference Too Great (CMS mode): Display and test record read correct.
- 18. CMS (CMS mode): Test record reads correct.

19. Stability Test (Level 3,T,S):
Ave: 0.0813 Std Dev: 0.0006 Rel Std Dev: 0.68%

Maintenance/Repairs/Notes: _____

Instrument is acceptable to be used in the field.

Charles E. Eda
Field Inspector's Signature

27 DEC. 2016
Date

Deb Shanaver
Reviewed by

27 Dec 16
Date

Intoxilyzer Test Record and Checklist
NDODAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/20/2016 16:33

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTUUVWX

32 characters per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTUUVWXYZ123456

42 characters per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTUUVWXYZ1234567890


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Menfull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 8765432
Standard Cyl #? 12345
Standard Expiration? 11/11/2016
Oper No? 666666

Flow Cal. Date: 08/15/2016
Slope 698
Intercept -608944

IR Calibration Date: 08/15/2016

	3um	9um
0th Coef(*100):	-14698	-21487
1st Coef(*100):	257953	133392
2nd Coef(*100):	2905	1310
H2O adj(mg/l*10k):	486	492

***** Test End *****



Operator Signature
CHARLES EDER

Remarks: INTERNAL
PRINT
TEST

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/20/2016 16:36

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=
|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=
|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 8765432
Standard Cyl #? 12345
Standard Expiration? 11/11/2016
Oper No? 666666

Flow Cal. Date: 08/15/2016
Slope 698
Intercept -608944

IR Calibration Date: 08/15/2016
 3um 9um

0th Coef(*100): -14698 -21487
1st Coef(*100): 257953 133392
2nd Coef(*100): 2905 1310
H2O adj(mg/l*10k): 486 492

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

EXTERNAL PRINT
TEST

Intoxilyzer Test Record and Checklist
 NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
 North Dakota Model 8000 SN 80-006666
 Location = TOXL 8164.14.00 09/16
 12/20/2016 16:38

Test	AC	Time
01 Room Air	0.000	16:39
02 *Subject Test	REF*	16:39
03 Room Air	0.000	16:40

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
 Sub DOB = 01/01/1982
 Sub Sex = Female Weight = 150
 Test = DUI Cit = R KEY REFUSAL
 Dr. Lic. = ND/DIS821456
 Lot No = 8765432
 Cyl No = 12345
 Expiration Date = 11/11/2016
 County = 08 Oper No. = 666666



Operator Signature
 CHARLES EDER

Remarks: *R KEY*
REFUSAL

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/20/2016 16:41

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: **DIAGNOSTIC**

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/21/2016 10:07

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:07
02 Std. Sol.	0.000	10:08
03 Room Air	0.000	10:08
04 Std. Sol.	0.000	10:09
05 Room Air	0.000	10:10
06 Std. Sol.	0.000	10:10
07 Room Air	0.000	10:11

08 Sim Temp = 34.0°C

Simul Ser No = MP3064
Std Sol No = MILLIQ H2O
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

0.000 AC

ZERO AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/21/2016 10:13

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:13
02 Std. Sol.	0.015	10:14
03 Room Air	0.000	10:15
04 Std. Sol.	0.015	10:15
05 Room Air	0.000	10:16
06 Std. Sol.	0.015	10:17
07 Room Air	0.000	10:17

08 Sim Temp = 34.0°C

Simul Ser No = MP3060

Std Sol No = 201501A

County = 08

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Low AC
0.015 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/21/2016 10:26

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:26
02 Std. Sol.	0.393	10:27
03 Room Air	0.000	10:28
04 Std. Sol.	0.394	10:29
05 Room Air	0.000	10:29
06 Std. Sol.	0.393	10:30
07 Room Air	0.000	10:30

08 Sim Temp = 34.0°C

Simul Ser No = MP3063

Std Sol No = 16102

County = 08

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

HIGH AC

0.400 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/21/2016 10:40

Test	AC	Time
01 Room Air	0.000	10:41
02 *Subject Test	INT*	10:41
03 Room Air	0.000	10:42

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = 0.10AC+0.05ACE
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent
0.10AC + 0.05% ACETONE

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/21/2016 10:43

Test	AC	Time
01 Room Air	0.000	10:44
02 *Subject Test	INT*	10:46
03 Room Air	0.000	10:46

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = 0.05 ACETONE
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent
0.05% ACETONE

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/21/2016 11:23

Test	AC	Time
01 Room Air	0.000	11:24
02 *Subject Test	INT*	11:24
03 Room Air	0.000	11:25

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = 0.10 METHANOL
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent
0.10% METHANOL

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/21/2016 11:26

Test	AC	Time
01 Room Air	0.000	11:26
02 *Subject Test	INT*	11:27
03 Room Air	0.000	11:28

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = 0.10 IPA
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent
0.10% ISOPROPANOL

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/21/2016 12:37

Test	AC	Time
01 Room Air	AMB*	12:38
02 Room Air	0.000	12:39

*Invalid Test
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = AMBIENT FAIL
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

*AMBIENT
FAIL*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/21/2016 12:40

Test	AC	Time
01 Room Air	0.000	12:40
02 Invalid Sample	X.XXX	12:41
03 Room Air	0.000	12:41

*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INVALID SAMPLE
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

INVALID SAMPLE

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/21/2016 12:42

Test	AC	Time
01 Room Air	0.000	12:43
02 *Subject Test	ABT*	12:43
03 Room Air	0.000	12:43

*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = SEQ ABORTED
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

**SEQUENCE
ABORTED**

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/21/2016 12:48

Test	AC	Time
01 Diagnostic	OK	12:49
02 Room Air	0.000	12:49
03 Subject Test 1	0.217	12:50
04 Room Air	0.000	12:52
05 Std. Gas	0.086*	12:53
06 Room Air	0.000	12:54

*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = CAL OUT TOL
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *CALIBRATION CHECK
OUT OF TOLERANCE*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/27/2016 13:27

Test	AC	Time
01 Diagnostic	OK	13:28
02 Room Air	0.000	13:28
03 *Subject Test	0.000*	13:31
04 Room Air	0.000	13:33
05 Std. Gas	0.082	13:34
06 Room Air	0.000	13:36
07 *Subject Test	0.000*	13:39
08 Room Air	0.000	13:40
No RFI Detected		

*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = DEFICIENT SAMP
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
CHARLES EDER

Remarks: **DEFICIENT
SAMPLES**

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/27/2016 13:41

Test	AC	Time
01 Diagnostic	OK	13:43
02 Room Air	0.000	13:43
03 Subject Test 1	0.214	13:44
04 Room Air	0.000	13:46
05 Std. Gas	0.082	13:47
06 Room Air	IPS*	13:47
07 Room Air	0.000	13:48

*Improper Sample

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = IMPROPER SAMPLE
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *Improper Sample*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/27/2016 13:51

Test	AC	Time
01 Diagnostic	OK	13:52
02 Room Air	0.000	13:53
03 Subject Test 1	0.282*	13:53
04 Room Air	0.000	13:55
05 Std. Gas	0.080	13:56
06 Room Air	0.000	13:58
07 Subject Test 2	0.211*	13:58
08 Room Air	0.000	13:59

*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = DIFF TOO GREAT
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

**DIFFERENCE
TOO GREAT**

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/27/2016 14:01

Test	AC	Time
01 Diagnostic	OK	14:01
02 Room Air	0.000	14:02
03 Subject Test 1	0.212	14:02
04 Room Air	0.000	14:05
05 Std. Gas	0.080	14:06
06 Room Air	0.000	14:07
07 Subject Test 2	0.212	14:08
08 Room Air	0.000	14:09
09 Reported AC	0.212	14:08

Difference OK
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = CMS
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

C.M.S.

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
12/27/2016 14:10

Test	AC	Time
01 Room Air	0.000	14:10
02 Std. Gas	0.080	14:11
03 Room Air	0.000	14:11
04 Std. Gas	0.081	14:12
05 Room Air	0.000	14:12
06 Std. Gas	0.081	14:12
07 Room Air	0.000	14:13
08 Std. Gas	0.081	14:13
09 Room Air	0.000	14:14
10 Std. Gas	0.081	14:14
11 Room Air	0.000	14:15
12 Std. Gas	0.082	14:15
13 Room Air	0.000	14:16
14 Std. Gas	0.082	14:16
15 Room Air	0.000	14:16
16 Std. Gas	0.082	14:17
17 Room Air	0.000	14:17
18 Std. Gas	0.081	14:18
19 Room Air	0.000	14:18
20 Std. Gas	0.081	14:18
21 Room Air	0.000	14:19
22 Std. Gas	0.082	14:19
23 Room Air	0.000	14:20
24 Std. Gas	0.081	14:20
25 Room Air	0.000	14:21
26 Std. Gas	0.082	14:21
27 Room Air	0.000	14:22
28 Std. Gas	0.081	14:22
29 Room Air	0.000	14:22
30 Std. Gas	0.081	14:23
31 Room Air	0.000	14:23
32 Std. Gas	0.081	14:24
33 Room Air	0.000	14:24
34 Std. Gas	0.081	14:25
35 Room Air	0.000	14:25
36 Std. Gas	0.082	14:25
37 Room Air	0.000	14:26
38 Std. Gas	0.081	14:26
39 Room Air	0.000	14:27
40 Std. Gas	0.081	14:27
41 Room Air	0.000	14:28
42 Cal Check Stats:		
43 Average	0.0813	
44 Std Dev	0.0006	
45 Rel Std Dev	0.68%	



Operator Signature
CHARLES EDER

Remarks:

STABILITY TEST
0.080 AC F.G.S.

Form 106-I8000