

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-006670

A. GENERAL Setup and Checks:

1. Power On
 5 min. countdown during initial warm-up
 1 min. countdown to begin test after standby mode
2. Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:
 Level 1 (Operator) Level 2 (Field Inspector) Level 3 (Lab Personnel)
4. Printer Setup (Level 3,S,P)
5. Forms Load (CMI, Inc., CoBrA Software)

B. LEVEL 1 Setup and Checks:

1. Internal Printer Test and Paper Feed (Level 1,P). Attach record.
2. External Printer Test (Level 1,P). Attach record.
3. Barcode Reader (Level 1,B):
 Operator ID ND License R key: Refusal. Attach record.

C. LEVEL 2 Setup and Checks:

1. Continuous Air Blank (Level 2,A)
2. Diagnostic (Level 2,D). Attach record.
3. Prelim Data Entry (Level 2,E)
4. Ser Number/Version Info (Level 2,V):
Inst Ser#: 80-006670 DSP Ver: 1362.14
Core Ver: 8164.14 CPLD Ver: 7

D. LEVEL 3 Setup and Checks:

1. DVM Monitor (Level 3,D,D): 3: 12842 s: 1 c: 47.00
9: 13148 s: 1 B: 45 0
(RFI indicator with keyed radio): 1
2. Function (F) Keys: : F5 – Tone F9 - Fan
 F6 – Green LED F10 – Pump
 F7 – Red LED F11 – Solenoid
3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 945 CF: 1.07
4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 238 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T): C: 47.0 A: 49.5 I: 0.0
B: 45.0 V: 0.0
6. Tank Monitor (Level 3,D,G): P (display): 96 psi P (gauge): ~100 psi
7. General Setup (Level 3,S,G)
8. Configure Start Test (Level 3,S,T)
9. Network (Level 3,S,N)

E. TESTS (attach all test records):

Configure wet bath simulator standards before running tests (Level 1,S):

1. Zero AC (Level 1,C):
Sim ser #: MP3064 Lot #: MILLI QH₂O Exp date: N/A
Known value: 0.000 AC
Results: 0.000 AC 0.000 AC 0.000 AC
2. Low AC (Level 1,C):
Sim ser #: MP3003 Lot #: 201603H Exp date: 3-16-18
Known value: 0.010 AC
Results: 0.012 AC 0.011 AC 0.012 AC
3. High AC (Level 1,C):
Sim ser #: MP3063 Lot #: 16102 Exp date: 3-22-18
Known value: 0.400 AC
Results: 0.394 AC 0.397 AC 0.397 AC

Configure dry gas standard before running tests (Level 1,S):

4. Ethanol + 0.05% Acetone (Level 1,B):
Sim ser #: DR5135 Lot #: ICS4 Exp date: 9.9.17
Known value: 0.10 AC + 0.05% ACETONE
 Display reads "Interferent Detect"
5. Acetone (Level 1,B):
Sim ser #: DR5143 Lot #: ACE2 Exp date: 9.29.17
Known value: 0.05% ACETONE
 Display reads "Interferent Detect"
6. Methanol (Level 1,B):
Sim ser #: DR3377 Lot #: MET2 Exp date: 9.29.17
Known value: 0.10% METHANOL
 Display reads "Interferent Detect"
7. Isopropanol (Level 1,B):
Sim ser #: DR3453 Lot #: IPA2 Exp date: 9.29.17
Known value: 0.10% ISOPROPANOL
 Display reads "Interferent Detect"

- 8. Ambient Fail (Level 1,B): Display and test record read correct.
- 9. Invalid Sample (Level 1,B): Display and test record read correct.
- 10. Sequence Aborted (Level 1,B): Display and test record read correct.
- 11. LABA (Level 1,B): Test record reads correct.
- 12. Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- 13. RFI Detected (B,C, or CMS mode): Display and test record read correct.
- 14. Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 15. Deficient Samples (CMS mode): Display and test record read correct.
- 16. Improper Sample (CMS mode): Display and test record read correct.
- 17. Difference Too Great (CMS mode): Display and test record read correct.
- 18. CMS (CMS mode): Test record reads correct.
- 19. Stability Test (Level 3,T,S):

Ave: 0.0823 Std Dev: 0.0006 Rel Std Dev: 0.67%

Maintenance/Repairs/Notes: _____

Instrument is acceptable to be used in the field.

Charles E. Ed
Field Inspector's Signature

1/10/2017
Date

Deb Shanaver
Reviewed by

13 Jan '17
Date

Intoxilyzer Test Record and Checklist
NDDAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 60-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 13:48

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstuvwx
ABCDEFGHIJKLMNQRSTUWX

32 characters per line
abcdefghijklmnopqrstuvwxy123456
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line
abcdefghijklmnopqrstuvwxy1234567890abcdef
ABCDEFGHIJKLMNQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 19615080A5
Standard Cyl #? 26
Standard Expiration? 09/05/2017
Oper No? 666666

Flow Cal. Date: 08/15/2016
Slope 691
Intercept -681667

IR Calibration Date: 08/15/2016

	3um	9um
0th Coef(*100):	-36090	-26909
1st Coef(*100):	270161	131597
2nd Coef(*100):	3492	1526
H2O adj(mg/l*10k):	790	617

***** Test End *****



Operator Signature
CHARLES EDER

Remarks:

Internal Print
Test

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 13:50

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMN OPQRSTUVWXYZ!@#%&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMN OPQRSTUVWXYZ!@#%&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 80
Standard Cyl #? 1123456
Standard Expiration? 11/11/2016
Oper No? 666666

Flow Cal. Date: 08/15/2016
Slope 691
Intercept -681667

IR Calibration Date: 08/15/2016
 3um 9um

0th Coef(*100): -36090 -26909
1st Coef(*100): 270161 131597
2nd Coef(*100): 3492 1526
H2O adj(mg/l*10k): 790 617

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks: EXTERNAL PRINT
TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 13:52

Test	AC	Time
01 Room Air	0.000	13:53
02 *Subject Test	REF*	13:53
03 Room Air	0.000	13:54

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = R KEY REFUSAL
Dr. Lic. = ND/DIS821456
Lot No = 80
Cyl No = 1123456
Expiration Date = 11/11/2016
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: "R" KEY REFUSAL

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 13:55

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

DIAGNOSTIC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 14:03

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:04
02 Std. Sol.	0.000	14:04
03 Room Air	0.000	14:05
04 Std. Sol.	0.000	14:06
05 Room Air	0.000	14:06
06 Std. Sol.	0.000	14:07
07 Room Air	0.000	14:07

08 Sim Temp = 34.0°C

Simul Ser No = MP3064
Std Sol No = MILLIQ H2O
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

ZERO AC
0.000 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 14:11

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:11
02 Std. Sol.	0.012	14:12
03 Room Air	0.000	14:12
04 Std. Sol.	0.011	14:13
05 Room Air	0.000	14:14
06 Std. Sol.	0.012	14:14
07 Room Air	0.000	14:15

08 Sim Temp = 34.0°C

Simul Ser No = MP3003

Std Sol No = 201603H

County = 08

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Low AC*
0.010 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 14:21

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:21
02 Std. Sol.	0.394	14:22
03 Room Air	0.000	14:23
04 Std. Sol.	0.397	14:24
05 Room Air	0.000	14:24
06 Std. Sol.	0.397	14:25
07 Room Air	0.000	14:25

08 Sim Temp = 34.0°C

Simul Ser No = MP3063
Std Sol No = 16102
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

HIGH AC
0.400 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 14:30

Test	AC	Time
01 Room Air	0.000	14:31
02 *Subject Test	INT*	14:31
03 Room Air	0.000	14:32

*Invalid Test
Interferent Detected

Sub Name = TEST, BOBBY NONE
Sub DOB = 01/01/1985
Sub Sex = Male Weight = NA
Test = DUI Cit = 0.10AC+0.05ACE
Dr. Lic. = ND/TES851012
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent

CEE
0.10 AC + ~~0.05~~ 0.05% ACETONE

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 14:33

Test	AC	Time
01 Room Air	0.000	14:34
02 *Subject Test	INT*	14:34
03 Room Air	0.000	14:35

*Invalid Test
Interferent Detected

Sub Name = TEST, BOBBY NONE
Sub DOB = 01/01/1985
Sub Sex = Male Weight = 300
Test = DUI Cit = 0.05 ACETONE
Dr. Lic. = ND/TES851012
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent
0.05% ACETONE

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 14:35

Test	AC	Time
01 Room Air	0.000	14:36
02 *Subject Test	INT*	14:36
03 Room Air	0.000	14:37

*Invalid Test
Interferent Detected

Sub Name = TEST, BOBBY NONE
Sub DOB = 01/01/1985
Sub Sex = Male Weight = 456
Test = DUI Cit = 0.10 METHANOL
Dr. Lic. = ND/TES851012
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

0.10 % METHANOL

Interferent

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 14:37

Test	AC	Time
01 Room Air	0.000	14:38
02 *Subject Test	INT*	14:38
03 Room Air	0.000	14:39

*Invalid Test
Interferent Detected

Sub Name = TEST, BOBBY NONE
Sub DOB = 01/01/1985
Sub Sex = Male Weight = 123
Test = DUI Cit = 0.10 IPA
Dr. Lic. = ND/TES851012
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent

0.10% ISOPROPA

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 14:40

Test	AC	Time
01 Room Air	AMB*	14:41
02 Room Air	0.000	14:41

*Invalid Test
Check Ambient Conditions

Sub Name = TEST, BOBBY NONE
Sub DOB = 01/01/1985
Sub Sex = Male Weight = 234
Test = DUI Cit = AMBIENT FAIL
Dr. Lic. = ND/TES851012
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

AMBIENT FAIL

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 14:43

Test	AC	Time
01 Room Air	0.000	14:44
02 Invalid Sample	X.XXX	14:44
03 Room Air	0.000	14:45

*Invalid Test - Mouth Alcohol

Sub Name = TEST, BOBBY NONE
Sub DOB = 01/01/1985
Sub Sex = Male Weight = 111
Test = DUI Cit = INVALID SAMPLE
Dr. Lic. = ND/TES851012
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *INVALID SAMPLE*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 14:45

Test	AC	Time
01 Room Air	0.000	14:46
02 *Subject Test	ABT*	14:46
03 Room Air	0.000	14:47

*Sequence Aborted

Sub Name = TEST, BOBBY NONE
Sub DOB = 01/01/1985
Sub Sex = Male Weight = 67
Test = DUI Cit = SEQUENCE ABORT
Dr. Lic. = ND/TES851012
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: SEQUENCE ABORTED

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 14:47

Test	AC	Time
01 Room Air	0.000	14:48
02 Subject Test 1	0.212	14:48
03 Room Air	0.000	14:49
04 Reported AC	0.212	14:48

No RFI Detected

Sub Name = TEST, BOBBY NONE
Sub DOB = 01/01/1985
Sub Sex = Male Weight = 345
Test = DUI Cit = ABA MODE
Dr. Lic. = ND/TES851012
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: ABA MODE
 SINGLE BREATH TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 14:49

Test	AC	Time
01 Diagnostic	OK	14:50
02 Room Air	0.000	14:51
03 Subject Test 1	0.014	14:52
04 Room Air	0.000	14:54
05 Std. Gas	0.090*	14:55
06 Room Air	0.000	14:55

*Cal Check Out of Tolerance

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = 444
Test = DUI Cit = CAL OUT OF TOL
Dr. Lic. = ND/TES989643
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *CALIBRATION CHECK
OUT OF TOLERANCE*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 14:58

Test	AC	Time
01 Diagnostic	OK	14:59
02 Room Air	RFI*	14:59
03 Room Air	0.000	15:00

*Invalid Test
Inhibited - RFI

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = 666
Test = DUI Cit = RFI DETECTED
Dr. Lic. = ND/TES989643
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *RFI DETECTED*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 15:01

Test	AC	Time
01 Diagnostic	OK	15:02
02 Room Air	0.000	15:02
03 *Subject Test	0.000*	15:05
04 Room Air	0.000	15:12
05 Std. Gas	0.082	15:13
06 Room Air	0.000	15:14
07 *Subject Test	0.000*	15:17
08 Room Air	0.000	15:18

No RFI Detected

*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = 987
Test = DUI Cit = DEFICIENT SAMP
Dr. Lic. = ND/TES989643
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

DEFICIENT
SAMPLES

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 15:19

Test	AC	Time
01 Diagnostic	OK	15:20
02 Room Air	0.000	15:20
03 Reference	IPS*	15:20
04 Room Air	0.000	15:21

*Improper Sample

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = 234
Test = DUI Cit = IMPROPER SAMP
Dr. Lic. = ND/TES989643
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

*IMPROPER
SAMPLE*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 15:21

Test	AC	Time
01 Diagnostic	OK	15:23
02 Room Air	0.000	15:23
03 Subject Test 1	0.273*	15:23
04 Room Air	0.000	15:26
05 Std. Gas	0.082	15:27
06 Room Air	0.000	15:28
07 Subject Test 2	0.210*	15:29
08 Room Air	0.000	15:29

*Difference Too Great

Sub Name = TEST, AAA FFF
Sub DOB = 01/01/1994
Sub Sex = Male Weight = 765
Test = DUI Cit = DIFF TOO GREAT
Dr. Lic. = ND/TES940351
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: DIFFERENCE TOO
GREAT

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 15:31

Test	AC	Time
01 Diagnostic	OK	15:32
02 Room Air	0.000	15:32
03 Subject Test 1	0.272	15:32
04 Room Air	0.000	15:35
05 Std. Gas	0.081	15:36
06 Room Air	0.000	15:37
07 Subject Test 2	0.272	15:38
08 Room Air	0.000	15:39
09 Reported AC	0.272	15:32

Difference OK
No RFI Detected

Sub Name = TEST, AAA FFF
Sub DOB = 01/01/1994
Sub Sex = Male Weight = 234
Test = DUI Cit = CMS MODE
Dr. Lic. = ND/TES940351
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
CHARLES EDER

Remarks: CMS MODE

EVIDENTIAL BREATH TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006670
Location = TOXL 8164.14.00 09/16
01/10/2017 15:40

Test	AC	Time
01 Room Air	0.000	15:41
02 Std. Gas	0.081	15:41
03 Room Air	0.000	15:41
04 Std. Gas	0.082	15:42
05 Room Air	0.000	15:42
06 Std. Gas	0.082	15:43
07 Room Air	0.000	15:43
08 Std. Gas	0.082	15:44
09 Room Air	0.000	15:44
10 Std. Gas	0.082	15:44
11 Room Air	0.000	15:45
12 Std. Gas	0.082	15:45
13 Room Air	0.000	15:46
14 Std. Gas	0.082	15:46
15 Room Air	0.000	15:47
16 Std. Gas	0.083	15:47
17 Room Air	0.000	15:48
18 Std. Gas	0.083	15:48
19 Room Air	0.000	15:48
20 Std. Gas	0.082	15:49
21 Room Air	0.000	15:49
22 Std. Gas	0.083	15:50
23 Room Air	0.000	15:50
24 Std. Gas	0.083	15:50
25 Room Air	0.000	15:51
26 Std. Gas	0.082	15:51
27 Room Air	0.000	15:52
28 Std. Gas	0.082	15:52
29 Room Air	0.000	15:53
30 Std. Gas	0.083	15:53
31 Room Air	0.000	15:54
32 Std. Gas	0.082	15:54
33 Room Air	0.000	15:54
34 Std. Gas	0.082	15:55
35 Room Air	0.000	15:55
36 Std. Gas	0.083	15:56
37 Room Air	0.000	15:56
38 Std. Gas	0.082	15:56
39 Room Air	0.000	15:57
40 Std. Gas	0.082	15:57
41 Room Air	0.000	15:58
42 Cal Check Stats:		
43 Average	0.0823	
44 Std Dev	0.0006	
45 Rel Std Dev	0.67%	



Operator Signature
CHARLES EDER

Remarks:

STABILITY TEST

Form 106-I8000

LOT: 19615080A5

CYL: 026

EXP: 9/5/2017