

**INTOXILYZER® 8000 INITIAL INSPECTION**

Intoxilyzer® 8000 Serial Number: 80-006689

**A. GENERAL Setup and Checks:**

1.  Power On
  - 5 min. countdown during initial warm-up
  - 1 min. countdown to begin test after standby mode
2.  Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:
  - Level 1 (Operator)    Level 2 (Field Inspector)    Level 3 (Lab Personnel)
4.  Printer Setup (Level 3,S,P)
5.  Forms Load (CMI, Inc., CoBrA Software)

**B. LEVEL 1 Setup and Checks:**

1.  Internal Printer Test and Paper Feed (Level 1,P). Attach record.
2.  External Printer Test (Level 1,P). Attach record.
3. Barcode Reader (Level 1,B):
  - Operator ID    ND License    R key: Refusal. Attach record.

**C. LEVEL 2 Setup and Checks:**

1.  Continuous Air Blank (Level 2,A)
2.  Diagnostic (Level 2,D). Attach record.
3.  Prelim Data Entry (Level 2,E)
4. Ser Number/Version Info (Level 2,V):
  - Inst Ser#: 80-006689   DSP Ver: 1362.14
  - Core Ver: 8164.14   CPLD Ver: 7

**D. LEVEL 3 Setup and Checks:**

1. DVM Monitor (Level 3,D,D): 3: 12934   s: 0   c: 47.00  
9: 13321   s: 1   B: 45.   0  
(RFI indicator with keyed radio): 1
2. Function (F) Keys :  F5 – Tone                       F9 - Fan  
 F6 – Green LED                       F10 – Pump  
 F7 – Red LED                           F11 – Solenoid
3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 947   CF: 1.06
4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 233   F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T): C: 47.0 A: 50.0 I: 0.0  
B: 45.0 V: 0.0
6. Tank Monitor (Level 3,D,G): P (display): 708 psi P (gauge): -725 psi
7.  General Setup (Level 3,S,G)
8.  Configure Start Test (Level 3,S,T)
9.  Network (Level 3,S,N)

**E. TESTS (attach all test records):**

**Configure wet bath simulator standards before running tests (Level 1,S):**

1. Zero AC (Level 1,C): <sup>4 eff</sup>  
Sim ser #: MP3068 Lot #: MILLI Q H<sub>2</sub>O Exp date: NA  
Known value: 0.000 AC  
Results: 0.000 AC 0.000 AC 0.000 AC
2. Low AC (Level 1,C):  
Sim ser #: MP3060 Lot #: 201507F Exp date: 7/28/17  
Known value: 0.020 AC  
Results: 0.018 AC 0.019 AC 0.019 AC
3. High AC (Level 1,C):  
Sim ser #: MP3070 Lot #: 16060 Exp date: 02.10.18  
Known value: 0.300 AC  
Results: 0.297 AC 0.297 AC 0.297 AC

**Configure dry gas standard before running tests (Level 1,S):**

4. Ethanol + 0.05% Acetone (Level 1,B):  
Sim ser #: DR5135 Lot #: ICS4 Exp date: 9/9/17  
Known value: 0.10AC + 0.05% ACETONE  
 Display reads "Interferent Detect"
5. Acetone (Level 1,B):  
Sim ser #: DR5143 Lot #: ACE2 Exp date: 9/29/17  
Known value: 0.05% ACETONE  
 Display reads "Interferent Detect"
6. Methanol (Level 1,B):  
Sim ser #: DR3377 Lot #: MET2 Exp date: 9/29/17  
Known value: 0.10% METHANOL  
 Display reads "Interferent Detect"
7. Isopropanol (Level 1,B):  
Sim ser #: DR3453 Lot #: IPA2 Exp date: 9/29/17  
Known value: 0.10% ISOPROPANOL  
 Display reads "Interferent Detect"

- 8.  Ambient Fail (Level 1,B): Display and test record read correct.
- 9.  Invalid Sample (Level 1,B): Display and test record read correct.
- 10.  Sequence Aborted (Level 1,B): Display and test record read correct.
- 11.  ABA (Level 1,B): Test record reads correct.
- 12.  Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- 13.  RFI Detected (B,C, or CMS mode): Display and test record read correct.
- 14.  Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 15.  Deficient Samples (CMS mode): Display and test record read correct.
- 16.  Improper Sample (CMS mode): Display and test record read correct.
- 17.  Difference Too Great (CMS mode): Display and test record read correct.
- 18.  CMS (CMS mode): Test record reads correct.
- 19. Stability Test (Level 3,T,S):

Ave: 0.0811 Std Dev: 0.0006 Rel Std Dev: 0.69%

Maintenance/Repairs/Notes: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

Instrument is acceptable to be used in the field.

Charles E. Eda  
Field Inspector's Signature

3/23/2017  
Date

Deb Shearer  
Reviewed by

24 Mar 17  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-006689  
Location = TOXL 8164.14.00 09/16  
03/23/2017 13:17

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
abcdefghijkl  
ABCDEFGHIJKL

16 char per line  
abcdefghijklmnop  
ABCDEFGHIJKLMN

21 chars per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNQRSTU

24 characters per line  
abcdefghijklmnopqrstuwx  
ABCDEFGHIJKLMNQRSTUWX

32 characters per line  
abcdefghijklmnopqrstuwx123456  
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line  
abcdefghijklmnopqrstuwx1234567890abcdef  
ABCDEFGHIJKLMNQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup  
Data Entry Mode: Enabled  
Start Test Sequence: DABACABA  
Display Prelim Rslt? Yes  
Display Third Digit? Yes  
Inhib Printer(Y/N)? No  
Display Volume? No  
Disable On Menfull? Yes  
# of Print Copies? 1  
Select Std (D/W/I)? Dry  
Standard Value? 0.080  
Standard Lot #? 17316080A2  
Standard Cyl #? 26  
Standard Expiration? 08/05/2018  
Oper No? 666666

Flow Cal. Date: 09/07/2016  
Slope 688  
Intercept -608758

IR Calibration Date: 09/07/2016  
3um 9um

0th Coef(\*100): -26009 -27854

1st Coef(\*100): 275168 133238

2nd Coef(\*100): 2128 1198

H2O adj(ng/l\*10K): 643 602

\*\*\*\*\* Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks: INTERNAL  
PRINTER  
TEST  
Form 106-18000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006689  
Location = TOXL                    8164.14.00 09/16  
03/23/2017                            13:19

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            17316080A2  
Standard Cyl #?            26  
Standard Expiration?        08/05/2018  
Oper No?                    666666

Flow Cal. Date:            09/07/2016  
Slope                        688  
Intercept                   -608758

IR Calibration Date:        09/07/2016  
                                    3um                            9um

-----  
0th Coef(\*100):            -26009                        -27854  
1st Coef(\*100):            275168                        133238  
2nd Coef(\*100):            2128                            1198  
H2O adj(mg/l\*10k):        643                            602

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:    EXTERNAL PRINTER  
                                    TEST



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006689  
Location = TOXL      8164.14.00 09/16  
03/23/2017      13:22

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *DIAGNOSTIC*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006689  
Location = TOXL      8164.14.00 09/16  
03/23/2017      13:35

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:36
02 Std. Sol.	0.000	13:37
03 Room Air	0.000	13:37
04 Std. Sol.	0.000	13:38
05 Room Air	0.000	13:38
06 Std. Sol.	0.000	13:39
07 Room Air	0.000	13:39

08 Sim Temp = 34.0°C

Simul Ser No = MP3064  
Std Sol No = MILLIQ H2O  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: ZERO AC  
0.000 AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006689  
Location = TOXL      8164.14.00 09/16  
03/23/2017      13:43

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:44
02 Std. Sol.	0.018	13:45
03 Room Air	0.000	13:45
04 Std. Sol.	0.019	13:46
05 Room Air	0.000	13:47
06 Std. Sol.	0.019	13:47
07 Room Air	0.000	13:48


08 Sim Temp = 34.0°C

Simul Ser No = MP3060

Std Sol No = 201507F

County = 08

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Low AC*

*0.020 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006689  
Location = TOXL      8164.14.00 09/16  
03/23/2017      14:11

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:11
02 Std. Sol.	0.297	14:12
03 Room Air	0.000	14:13
04 Std. Sol.	0.297	14:14
05 Room Air	0.000	14:14
06 Std. Sol.	0.297	14:15
07 Room Air	0.000	14:15

08 Sim Temp = 34.0°C

Simul Ser No = MP3070  
Std Sol No = 16060  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: HIGH AC  
0.300 AC


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006689  
Location = TOXL      8164.14.00 09/16  
03/23/2017      14:17

Test	AC	Time
01 Room Air	0.000	14:18
02 *Subject Test	INT*	14:18
03 Room Air	0.000	14:19

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent*

*0.10 AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006689  
Location = TOXL      8164.14.00 09/16  
03/23/2017      14:20

Test	AC	Time
01 Room Air	0.000	14:21
02 *Subject Test	INT*	14:21
03 Room Air	0.000	14:22

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent* 0.05% ACETONE  
~~0.10% METHANOL~~ CEE

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006689  
Location = TOXL                    8164.14.00 09/16  
03/23/2017                            14:26

Test	AC	Time
01 Room Air	0.000	14:27
02 *Subject Test	INT*	14:27
03 Room Air	0.000	14:28

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent*  
*0.10% METHANOL*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006689  
Location = TOXL                    8164.14.00 09/16  
03/23/2017                            15:04

Test	AC	Time
01 Room Air	0.000	15:05
02 *Subject Test	INT*	15:05
03 Room Air	0.000	15:06

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent*  
*0.10% ISOPROPAOL*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006689  
Location = TOXL                    8164.14.00 09/16  
03/23/2017                            15:07

Test	AC	Time
01 Room Air	AMB*	15:08
02 Room Air	0.000	15:09

\*Invalid Test  
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = AMBIENT FAIL  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:                    *AMBIENT  
FAIL*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-006689  
Location = TOXL                                      8164.14.00 09/16  
03/23/2017    15:09

Test	AC	Time
01 Room Air	0.000	15:10
02 Invalid Sample	X.XXX	15:10
03 Room Air	0.000	15:11

\*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                                      Weight = 150  
Test = DUI                                      Cit = INVALID SAMPLE  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08                                      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *INVALID SAMPLE*




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006689  
Location = TOXL                8164.14.00 09/16  
03/23/2017                        15:13

Test	AC	Time
01 Room Air	0.000	15:14
02 *Subject Test	ABT*	15:14
03 Room Air	0.000	15:14

\*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                Weight = 150  
Test = DUI                        Cit = SEQUENCE ABORT  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08                        Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: SEQUENCE  
              ABORTED

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006689  
Location = TOXL                    8164.14.00 09/16  
03/23/2017                            15:22

Test	AC	Time
01 Room Air	0.000	15:23
02 Subject Test 1	0.000	15:23
03 Room Air	0.000	15:23
04 Reported AC	0.000	15:23

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = ABA MODE  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

ABA MORE

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006689  
Location = TOXL                    8164.14.00 09/16  
03/23/2017                            15:41

Test	AC	Time
01 Diagnostic	OK	15:43
02 Room Air	0.000	15:43
03 Subject Test 1	0.143	15:44
04 Room Air	0.000	15:46
05 Std. Gas	0.089*	15:47
06 Room Air	0.000	15:48

\*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = CAL CHK OUT TOL  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
OUT OF TOLERANCE

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006689  
Location = TOXL                 8164.14.00 09/16  
03/23/2017                        15:50

Test	AC	Time
01 Diagnostic	OK	15:51
02 Room Air	RFI*	15:51
03 Room Air	0.000	15:51

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                 Weight = 150  
Test = DUI                        Cit = RFI DETECTED  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08                        Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *RFI DETECTED*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006689  
Location = TOXL                    8164.14.00 09/16  
03/23/2017                            15:52

Test	AC	Time
01 Diagnostic	OK	15:53
02 Room Air	0.000	15:53
03 *Subject Test	0.000*	15:56
04 Room Air	0.000	15:59
05 Std. Gas	0.081	16:00
06 Room Air	0.000	16:01
07 *Subject Test	0.000*	16:04
08 Room Air	0.000	16:05


No RFI Detected

\*Deficient Sample - Value Printed was  
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = DEFICIENT SAMP  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08                            Oper No. = 666666

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the  
instructions displayed by the Intoxilyzer  
in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: DEFICIENT  
SAMPLES  
Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006689  
Location = TOXL      8164.14.00 09/16  
03/23/2017      16:06

Test	AC	Time
01 Diagnostic	OK	16:07
02 Room Air	IPS*	16:07
03 Room Air	0.000	16:07

\*Improper Sample

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = IMPROPER SAMPLE  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*Improper*

*Sample*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006689  
Location = TOXL                    8164.14.00 09/16  
03/23/2017                            16:09

Test	AC	Time
01 Diagnostic	OK	16:09
02 Room Air	0.000	16:10
03 Subject Test 1	0.099*	16:10
04 Room Air	0.000	16:13
05 Std. Gas	0.081	16:14
06 Room Air	0.000	16:15
07 Subject Test 2	0.142*	16:16
08 Room Air	0.000	16:17

\*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = DIFF TOO GREAT  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: DIFFERENCE  
            TOO GREAT

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006689  
Location = TOXL                    8164.14.00 09/16  
03/23/2017                            16:19

Test	AC	Time
01 Diagnostic	OK	16:21
02 Room Air	0.000	16:21
03 Subject Test 1	0.141	16:22
04 Room Air	0.000	16:24
05 Std. Gas	0.080	16:25
06 Room Air	0.000	16:27
07 Subject Test 2	0.141	16:27
08 Room Air	0.000	16:28
09 Reported AC	0.141	16:22

Difference OK  
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = CMS MODE  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08                            Oper No. = 666666

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *CMS MODE*



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006689  
Location = TOXL      8164.14.00 09/16  
03/23/2017      16:29

Test	AC	Time
01 Room Air	0.000	16:30
02 Std. Gas	0.080	16:30
03 Room Air	0.000	16:31
04 Std. Gas	0.081	16:31
05 Room Air	0.000	16:32
06 Std. Gas	0.081	16:32
07 Room Air	0.000	16:32
08 Std. Gas	0.081	16:33
09 Room Air	0.000	16:33
10 Std. Gas	0.081	16:34
11 Room Air	0.000	16:34
12 Std. Gas	0.082	16:34
13 Room Air	0.000	16:35
14 Std. Gas	0.081	16:35
15 Room Air	0.000	16:36
16 Std. Gas	0.080	16:36
17 Room Air	0.000	16:37
18 Std. Gas	0.081	16:37
19 Room Air	0.000	16:38
20 Std. Gas	0.081	16:38
21 Room Air	0.000	16:38
22 Std. Gas	0.081	16:39
23 Room Air	0.000	16:39
24 Std. Gas	0.081	16:40
25 Room Air	0.000	16:40
26 Std. Gas	0.082	16:41
27 Room Air	0.000	16:41
28 Std. Gas	0.081	16:41
29 Room Air	0.000	16:42
30 Std. Gas	0.081	16:42
31 Room Air	0.000	16:43
32 Std. Gas	0.082	16:43
33 Room Air	0.000	16:44
34 Std. Gas	0.081	16:44
35 Room Air	0.000	16:45
36 Std. Gas	0.081	16:45
37 Room Air	0.000	16:45
38 Std. Gas	0.082	16:46
39 Room Air	0.000	16:46
40 Std. Gas	0.081	16:47
41 Room Air	0.000	16:47
42 Cal Check Stats:		
43 Average	0.0811	
44 Std Dev	0.0006	
45 Rel Std Dev	0.69%	



Operator Signature  
CHARLES EDER

Remarks: STABILITY TEST  
0.080 AC

Form 106-I8000

Lot: 17316080 A2

Cyl: 026

Exp: 8/5/2018