

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-006671

A. GENERAL Setup and Checks:

1. Power On
 - 5 min. countdown during initial warm-up
 - 1 min. countdown to begin test after standby mode
2. Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:
 - Level 1 (Operator) Level 2 (Field Inspector) Level 3 (Lab Personnel)
4. Printer Setup (Level 3,S,P)
5. Forms Load (CMI, Inc., CoBrA Software)

B. LEVEL 1 Setup and Checks:

1. Internal Printer Test and Paper Feed (Level 1,P). Attach record.
2. External Printer Test (Level 1,P). Attach record.
3. Barcode Reader (Level 1,B):
 - Operator ID ND License R key: Refusal. Attach record.

C. LEVEL 2 Setup and Checks:

1. Continuous Air Blank (Level 2,A)
2. Diagnostic (Level 2,D). Attach record.
3. Prelim Data Entry (Level 2,E)
4. Ser Number/Version Info (Level 2,V):
 - Inst Ser#: 80-006671 DSP Ver: 1362.14
 - Core Ver: 8164.14 CPLD Ver: 7

D. LEVEL 3 Setup and Checks:

1. DVM Monitor (Level 3,D,D): 3: 12872 s: 1 c: 47.06
9: 13181 s: 1 B: 45 0
(RFI indicator with keyed radio): 1
2. Function (F) Keys: : F5 – Tone F9 - Fan
 F6 – Green LED F10 – Pump
 F7 – Red LED F11 – Solenoid
3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 967 CF: 1.04
4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 252 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T): C: 47.0 A: 49.5 I: 0.0
B: 95.0 V: 0.0
6. Tank Monitor (Level 3,D,G): P (display): 917 psi P (gauge): ~900 psi
7. General Setup (Level 3,S,G)
8. Configure Start Test (Level 3,S,T)
9. Network (Level 3,S,N)

E. TESTS (attach all test records):

Configure wet bath simulator standards before running tests (Level 1,S):

1. Zero AC (Level 1,C):
Sim ser #: MP3064 Lot #: MILL Q H₂O Exp date: NA
Known value: 0.000 AC
Results: 0.000 AC 0.000 AC 0.000 AC
2. Low AC (Level 1,C):
Sim ser #: MP3003 Lot #: 201603H Exp date: 3.16.18
Known value: 0.010 AC
Results: 0.011 AC 0.011 AC 0.012 AC
3. High AC (Level 1,C):
Sim ser #: MP3063 Lot #: 16102 Exp date: 3.22.18
Known value: 0.400 AC
Results: 0.394 AC 0.396 AC 0.396 AC

Configure dry gas standard before running tests (Level 1,S):

4. Ethanol + 0.05% Acetone (Level 1,B):
Sim ser #: DR5135 Lot #: ICS4 Exp date: 9.9.2017
Known value: 0.10 AC + 0.05% ACETONE
 Display reads "Interferent Detect"
5. Acetone (Level 1,B):
Sim ser #: DR5143 Lot #: ACE2 Exp date: 9.29.17
Known value: 0.05% ACETONE
 Display reads "Interferent Detect"
6. Methanol (Level 1,B):
Sim ser #: DR3377 Lot #: MET2 Exp date: 9.29.17
Known value: 0.10% METHANOL
 Display reads "Interferent Detect"
7. Isopropanol (Level 1,B):
Sim ser #: DR3453 Lot #: IPA2 Exp date: 9.29.17
Known value: 0.10% ISOPROPANOL
 Display reads "Interferent Detect"

- 8. Ambient Fail (Level 1,B): Display and test record read correct.
- 9. Invalid Sample (Level 1,B): Display and test record read correct.
- 10. Sequence Aborted (Level 1,B): Display and test record read correct.
- 11. ABA (Level 1,B): Test record reads correct.
- 12. Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- 13. RFI Detected (B,C, or CMS mode): Display and test record read correct.
- 14. Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 15. Deficient Samples (CMS mode): Display and test record read correct.
- 16. Improper Sample (CMS mode): Display and test record read correct.
- 17. Difference Too Great (CMS mode): Display and test record read correct.
- 18. CMS (CMS mode): Test record reads correct.
- 19. Stability Test (Level 3,T,S):

Ave: 0.0831 Std Dev: 0.0008 Rel Std Dev: 1.00%

Maintenance/Repairs/Notes: INTERNAL PRINT TEST WAS RE-RUN DUE
TO TAPE OBSCURING (FADING) CONTENT OF ORIGINAL
PRINT TEST. 1/24/2017 Charles E. Edr

Instrument is acceptable to be used in the field.

Charles E. Edr
Charles E. Edr
Field Inspector's Signature

1/24/2017
1/13/17
Date

Deb Shanever
Reviewed by

13 Jan '17^{DS}
Date 27 Mar 17

Intoxilyzer Test Record and Checklist
NDDAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 13:59

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstuwx
ABCDEFGHIJKLMNQRSTUWX

32 characters per line
abcdefghijklmnopqrstuwx123456
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line

Data Entry Mode: enauiet
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 80
Standard Cyl #? 1234567
Standard Expiration? 11/11/2017
Oper No? 666666

Flow Cal. Date: 08/22/2016
Slope 703
Intercept -691868

IR Calibration Date: 08/22/2016
3um 9um

0th Coef(*100):	-17135	-24345
1st Coef(*100):	254798	132856
2nd Coef(*100):	3386	1354
H2O adj(mg/l*10K):	610	539

***** Test End *****



Operator Signature
CHARLES EDER

Remarks: *Internal*
PRINT TEST
Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/24/2017 12:30

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstuvw
ABCDEFGHIJKLMNQRSTUW

32 characters per line
abcdefghijklmnopqrstuvwxy123456
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line
abcdefghijklmnopqrstuvwxy1234567890abcdef
ABCDEFGHIJKLMNQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On MemFull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 09216080A3
Standard Cyl #? 38
Standard Expiration? 05/05/2018
Oper No? 666666

Flow Cal. Date: 08/22/2016
Slope 703
Intercept -691868

IR Calibration Date: 08/22/2016

	3um	9um
0th Coef(*100):	-17135	-24345
1st Coef(*100):	254798	132856
2nd Coef(*100):	3386	1354
H2O adj(ng/l)*10K):	610	539

***** Test End *****



Operator Signature
CHARLES EDER

Remarks: INTERNAL
PRINT TEST

Form 106-18000

ORIGINAL PRINT OUT
HAD SOME CONTENT
OBSCURED BY TAPE
RE-PRINT

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 14:01

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 80
Standard Cyl #? 1234567
Standard Expiration? 11/11/2017
Oper No? 666666

Flow Cal. Date: 08/22/2016
Slope 703
Intercept -691868

IR Calibration Date: 08/22/2016
 3um 9um

0th Coef(*100): -17135 -24345
1st Coef(*100): 254798 132856
2nd Coef(*100): 3386 1354
H2O adj(mg/l*10k): 610 539

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

EXTERNAL PRINT
TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 14:10

Test	AC	Time
01 Room Air	0.000	14:11
02 *Subject Test	REF*	14:11
03 Room Air	0.000	14:11

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 03/01/1994
Sub Sex = Male Weight = 98
Test = DUI Cit = R KEY REFUSAL
Dr. Lic. = ND/DIS941449
Lot No = 80
Cyl No = 1234567
Expiration Date = 11/11/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: "R" KEY REFUSAL

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 14:12

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

DIAGNOSTIC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 14:19

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:19
02 Std. Sol.	0.000	14:20
03 Room Air	0.000	14:21
04 Std. Sol.	0.000	14:21
05 Room Air	0.000	14:22
06 Std. Sol.	0.000	14:23
07 Room Air	0.000	14:23

08 Sim Temp = 34.0°C

Simul Ser No = MP3064
Std Sol No = MILLIQ H2O
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: ZERO AC
0.000AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 14:28

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:28
02 Std. Sol.	0.011	14:29
03 Room Air	0.000	14:30
04 Std. Sol.	0.011	14:30
05 Room Air	0.000	14:31
06 Std. Sol.	0.012	14:31
07 Room Air	0.000	14:32

08 Sim Temp = 34.0°C

Simul Ser No = MP3003

Std Sol No = 201603H

County = 08

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Low AC*

0.010 AC
Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 14:34

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:35
02 Std. Sol.	0.394	14:36
03 Room Air	0.000	14:36
04 Std. Sol.	0.396	14:37
05 Room Air	0.000	14:38
06 Std. Sol.	0.396	14:38
07 Room Air	0.000	14:39

08 Sim Temp = 34.0°C

Simul Ser No = MP3063
Std Sol No = 16102
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *HIGH AC*
0.400 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 14:41

Test	AC	Time
01 Room Air	0.000	14:42
02 *Subject Test	INT*	14:42
03 Room Air	0.000	14:43

*Invalid Test
Interferent Detected

Sub Name = BILLS, THE BUFFALO
Sub DOB = 06/14/1993
Sub Sex = Male Weight = 78
Test = DUI Cit = 0.10AC+0.05ACE
Dr. Lic. = ND/BIL933277
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent
0.10 AC + 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 14:45

Test	AC	Time
01 Room Air	0.000	14:46
02 *Subject Test	INT*	14:47
03 Room Air	0.000	14:47

*Invalid Test
Interferent Detected

Sub Name = BILLS, THE BUFFALO
Sub DOB = 06/14/1993
Sub Sex = Male Weight = 78
Test = DUI Cit = 0.05 ACETONE
Dr. Lic. = ND/BIL933277
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent
0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 14:48

Test	AC	Time
01 Room Air	0.000	14:50
02 *Subject Test	INT*	14:51
03 Room Air	0.000	14:51

*Invalid Test
Interferent Detected

Sub Name = AREENOUGHALREADY, APRIL SHOWERS
Sub DOB = 02/28/1994
Sub Sex = Female Weight = 67
Test = DUI Cit = 0.10 METHANOL
Dr. Lic. = ND/ARE942282
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent
0.10% METHANOL

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 14:53

Test	AC	Time
01 Room Air	0.000	14:53
02 *Subject Test	INT*	14:54
03 Room Air	0.000	14:54

*Invalid Test
Interferent Detected

Sub Name = AREENOUGHALREADY, APRIL SHOWERS
Sub DOB = 02/28/1994
Sub Sex = Female Weight = NA
Test = DUI Cit = 0.10 IPA
Dr. Lic. = ND/ARE942282
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent

0.10% ISOPROPAVOL
Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 14:56

Test	AC	Time
01 Room Air	AMB*	14:56
02 Room Air	0.000	14:57

*Invalid Test
Check Ambient Conditions

Sub Name = AREENOUGHALREADY, APRIL SHOWERS
Sub DOB = 02/28/1994
Sub Sex = Female Weight = 09
Test = DUI Cit = AMBIENT FAIL
Dr. Lic. = ND/ARE942282
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *AMBIENT FAIL*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 14:58

Test	AC	Time
01 Room Air	0.000	14:59
02 Invalid Sample	X.XXX	14:59
03 Room Air	0.000	15:00

*Invalid Test - Mouth Alcohol

Sub Name = BILLS, THE BUFFALO
Sub DOB = 06/14/1993
Sub Sex = Male Weight = 789
Test = DUI Cit = INVALID SAMPLE
Dr. Lic. = ND/BIL933277
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *INVALID SAMPLE*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 15:02

Test	AC	Time
01 Room Air	0.000	15:03
02 *Subject Test	ABT*	15:03
03 Room Air	0.000	15:03

*Sequence Aborted

Sub Name = AREENOUGHALREADY, APRIL SHOWERS
Sub DOB = 02/28/1994
Sub Sex = Female Weight = 567
Test = DUI Cit = SEQUENCE ABORT
Dr. Lic. = ND/ARE942282
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: SEQUENCE ABORTED

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 15:04

Test	AC	Time
01 Room Air	0.000	15:05
02 Subject Test 1	0.012	15:05
03 Room Air	0.000	15:06
04 Reported AC	0.012	15:05

No RFI Detected

Sub Name = BILLS, THE BUFFALO
Sub DOB = 06/14/1993
Sub Sex = Male Weight = 234
Test = DUI Cit = ABA MODE
Dr. Lic. = ND/BIL933277
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *ABA MODE*
SINGLE BREATH TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/12/2017 15:07

Test	AC	Time
01 Diagnostic	OK	15:08
02 Room Air	0.000	15:09
03 Subject Test 1	0.273	15:09
04 Room Air	0.000	15:11
05 Std. Gas	0.091*	15:12
06 Room Air	0.000	15:13

*Cal Check Out of Tolerance

Sub Name = TEST, FIVE DRIVER
Sub DOB = 06/15/1997
Sub Sex = Male Weight = 123
Test = DUI Cit = CAL OUT TOL
Dr. Lic. = ND/TES972018
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK OF CEC

OUT OF TOLERANCE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/13/2017 09:08

Test	AC	Time
01 Diagnostic	OK	09:10
02 Room Air	RFI*	09:10
03 Room Air	0.000	09:10

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 09/08/1995
Sub Sex = Female Weight = 90
Test = DUI Cit = RFI DETECTED
Dr. Lic. = ND/DIS950910
Lot No = 19615080A5
Cyl No = 26
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: RFI DETECTED

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/13/2017 09:15

Test	AC	Time
01 Diagnostic	OK	09:16
02 Room Air	0.000	09:16
03 *Subject Test	0.000*	09:19
04 Room Air	0.000	09:21
05 Std. Gas	0.083	09:22
06 Room Air	0.000	09:24
07 *Subject Test	0.000*	09:27
08 Room Air	0.000	09:28

No RFI Detected

*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 09/08/1995
Sub Sex = Female Weight = 120
Test = DUI Cit = DEFICIENT SAMP
Dr. Lic. = ND/DIS950910
Lot No = 09216080A3
Cyl No = 38
Expiration Date = 05/05/2018
County = 08 Oper No. = 666666

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
CHARLES EDER

Remarks: DEFICIENT SAMPLES

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/13/2017 10:12

Test	AC	Time
01 Diagnostic	OK	10:13
02 Room Air	IPS*	10:13
03 Room Air	0.000	10:14

*Improper Sample

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 03/01/1982
Sub Sex = Male Weight = 180
Test = DUI Cit = IMPROPER SAMPLE
Dr. Lic. = ND/DIS821449
Lot No = 09216080A3
Cyl No = 38
Expiration Date = 05/05/2018
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

Improper Sample

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/13/2017 10:14

Test	AC	Time
01 Diagnostic	OK	10:15
02 Room Air	0.000	10:15
03 Subject Test 1	0.272*	10:16
04 Room Air	0.000	10:18
05 Std. Gas	0.083	10:19
06 Room Air	0.000	10:21
07 Subject Test 2	0.211*	10:21
08 Room Air	0.000	10:22

*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1994
Sub Sex = Female Weight = 150
Test = DUI Cit = DIFF TOO GREAT
Dr. Lic. = ND/DIS941456
Lot No = 09216080A3
Cyl No = 38
Expiration Date = 05/05/2018
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: DIFFERENCE TOO
GREAT

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/13/2017 10:43

Test	AC	Time
01 Diagnostic	OK	10:44
02 Room Air	0.000	10:45
03 Subject Test 1	0.272	10:45
04 Room Air	0.000	10:47
05 Std. Gas	0.082	10:49
06 Room Air	0.000	10:50
07 Subject Test 2	0.270	10:51
08 Room Air	0.000	10:51
09 Reported AC	0.270	10:51

Difference OK
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 09/08/1995
Sub Sex = Female Weight = 120
Test = DUI Cit = CMS MODE
Dr. Lic. = ND/DIS950910
Lot No = 09216080A3
Cyl No = 38
Expiration Date = 05/05/2018
County = 08 Oper No. = 666666

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: CMS MODE

EVIDENTIAL BREATH TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = TOXL 8164.14.00 09/16
01/13/2017 10:57

Test	AC	Time
01 Room Air	0.000	10:57
02 Std. Gas	0.081	10:58
03 Room Air	0.000	10:58
04 Std. Gas	0.082	10:58
05 Room Air	0.000	10:59
06 Std. Gas	0.082	10:59
07 Room Air	0.000	11:00
08 Std. Gas	0.084	11:00
09 Room Air	0.000	11:01
10 Std. Gas	0.083	11:01
11 Room Air	0.000	11:02
12 Std. Gas	0.083	11:02
13 Room Air	0.000	11:02
14 Std. Gas	0.082	11:03
15 Room Air	0.000	11:03
16 Std. Gas	0.083	11:04
17 Room Air	0.000	11:04
18 Std. Gas	0.084	11:05
19 Room Air	0.000	11:05
20 Std. Gas	0.084	11:05
21 Room Air	0.000	11:06
22 Std. Gas	0.084	11:06
23 Room Air	0.000	11:07
24 Std. Gas	0.084	11:07
25 Room Air	0.000	11:08
26 Std. Gas	0.083	11:08
27 Room Air	0.000	11:09
28 Std. Gas	0.083	11:09
29 Room Air	0.000	11:09
30 Std. Gas	0.083	11:10
31 Room Air	0.000	11:10
32 Std. Gas	0.083	11:11
33 Room Air	0.000	11:11
34 Std. Gas	0.083	11:11
35 Room Air	0.000	11:12
36 Std. Gas	0.083	11:12
37 Room Air	0.000	11:13
38 Std. Gas	0.083	11:13
39 Room Air	0.000	11:14
40 Std. Gas	0.084	11:14
41 Room Air	0.000	11:15
42 Cal Check Stats:		
43 Average	0.0831	
44 Std Dev	0.0008	
45 Rel Std Dev	1.00%	



Operator signature
CHARLES EDER

Remarks: STABILITY TEST

Form 106-I8000

LOT: 09216080A3

Exp: 5/5/2018

CYL: 038