

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-006666

Location: TOXL

- Check When Done
- A. Is the breath tube heated? ✓
  - B. Is the diagnostic check complete? ✓
  - C. Is the time, date, and year correct (re-set if necessary)? ✓
  - D. Print test. (Attach test record.) ✓
  - E. Low AC. Use  $\leq 0.03$  AC in ACA mode. (Attach test record.)  
Sim SN MP3003 Lot No. 201610A Rep AC 0.020 ✓
  - F. High AC. Use  $\geq 0.25$  AC in ACA mode. (Attach test record.)  
Sim SN MP3063 Lot No. 16102 Rep AC 0.400 ✓
  - G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.  
(Attach test record.)  
Sim SN DR5135 Lot No. ICS4 Solution 0.10 AC + 0.05% ACETONE ✓
  - H. RFI check. Run CMS mode. Key radio during CMS test.  
(Attach test record.) ✓
  - I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.  
(Attach test record.) 0.080 AC

Lot No. 17316080A2 Cyl No. 026 Exp. Date 8/5/2018

Test 1 <u>0.080</u>	Test 4 <u>0.082</u>	Test 7 <u>0.081</u>
Test 2 <u>0.081</u>	Test 5 <u>0.081</u>	Test 8 <u>0.081</u>
Test 3 <u>0.081</u>	Test 6 <u>0.080</u>	Test 9 <u>0.081</u>

Average 0.081

- J. Tank pressure, Level 3 Func D Sub G P: 394 psi Reg ~400 psi
- K. Remarks/Maintenance record: ANNUAL INSPECTION

4.6.2017  
Date

11 Apr 17  
Date

[Signature]  
Field Inspector's Signature

Doh Shanaver  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006666  
Location = TOXL                    8164.14.00 09/16  
04/06/2017                            09:44

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:                    Enabled  
Start Test Sequence:                DABACABA  
Display Prelim Rslt?                Yes  
Display Third Digit?                Yes  
Inhib Printer(Y/N)?                No  
Display Volume?                    No  
Disable On Memfull?                Yes  
# of Print Copies?                  1  
Select Std (D/W/I)?                Dry  
Standard Value?                    0.080  
Standard Lot #?                    17316080A2  
Standard Cyl #?                    26  
Standard Expiration?                08/05/2018  
Oper No?                            666666

Flow Cal. Date:                    08/15/2016  
Slope                                698  
Intercept                            -608944

IR Calibration Date:                08/15/2016  
                                          3um                                9um

-----  
0th Coef(\*100):                    -14698                                -21487  
1st Coef(\*100):                    257953                                133392  
2nd Coef(\*100):                    2905                                    1310  
H2O adj(mg/l\*10k):                486                                    492

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006666  
Location = TOXL      8164.14.00 09/16  
04/06/2017      09:45

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:46
02 Std. Sol.	0.017	09:46
03 Room Air	0.000	09:47
04 Std. Sol.	0.017	09:48
05 Room Air	0.000	09:48
06 Std. Sol.	0.018	09:49
07 Room Air	0.000	09:49

08 Sim Temp = 34.0°C

Simul Ser No = MP3003  
Std Sol No = 201610A  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Low AC*  
*0.020 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006666  
Location = TOXL      8164.14.00 09/16  
04/06/2017      09:52

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:52
02 Std. Sol.	0.393	09:53
03 Room Air	0.000	09:54
04 Std. Sol.	0.395	09:54
05 Room Air	0.000	09:55
06 Std. Sol.	0.395	09:56
07 Room Air	0.000	09:56

08 Sim Temp = 34.0°C

Simul Ser No = MP3063

Std Sol No = 16102

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

HIGH AC  
0.400 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006666  
Location = TOXL      8164.14.00 09/16  
04/06/2017      09:58

Test	AC	Time
01 Room Air	0.000	09:59
02 *Subject Test	INT*	09:59
03 Room Air	0.000	10:00

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

~~RFE CHE 09E~~  
INTERFERENT 0.10AC + 0.05%  
ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006666  
Location = TOXL      8164.14.00 09/16  
04/06/2017      10:00

Test	AC	Time
01 Diagnostic	OK	10:01
02 Room Air	0.000	10:01
03 *Subject Test	RFI*	10:01
04 Room Air	0.000	10:02

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 26  
Expiration Date = 08/05/2018  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

RFI CHECK


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006666  
Location = TOXL      8164.14.00 09/16  
04/06/2017      10:02

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:03
02 Std. Gas	0.080	10:03
03 Room Air	0.000	10:03
04 Std. Gas	0.081	10:04
05 Room Air	0.000	10:04
06 Std. Gas	0.081	10:05
07 Room Air	0.000	10:05

Lot No = 17316080A2  
Cyl No = 26  
Exp Date = 08/05/2018  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

CALIBRATION CHECK  
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006666  
Location = TOXL      8164.14.00 09/16  
04/06/2017      10:06

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:06
02 Std. Gas	0.082	10:06
03 Room Air	0.000	10:07
04 Std. Gas	0.081	10:07
05 Room Air	0.000	10:08
06 Std. Gas	0.080	10:08
07 Room Air	0.000	10:09

Lot No = 17316080A2  
Cyl No = 26  
Exp Date = 08/05/2018  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006666  
Location = TOXL      8164.14.00 09/16  
04/06/2017      10:09

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:09
02 Std. Gas	0.081	10:10
03 Room Air	0.000	10:10
04 Std. Gas	0.081	10:11
05 Room Air	0.000	10:11
06 Std. Gas	0.081	10:12
07 Room Air	0.000	10:12

Lot No = 17316080A2  
Cyl No = 26  
Exp Date = 08/05/2018  
County = 08      Oper No. = 666666

  
Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC