

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-006688

Location: TOXL

- |  | Check When Done |
|--|-----------------|
| A. Is the breath tube heated?  | <u>X</u>        |
| B. Is the diagnostic check complete?   | <u>X</u>        |
| C. Is the time, date, and year correct (re-set if necessary)?  | <u>X</u>        |
| D. Print test. (Attach test record.)   | <u>X</u>        |
| E. Low AC. Use $\leq 0.03$ AC in ACA mode. (Attach test record.)<br>Sim SN <u>0R3846</u> Lot No. <u>16200</u> Rep AC <u>0.020</u>  | <u>X</u>        |
| F. High AC. Use $\geq 0.25$ AC in ACA mode. (Attach test record.)<br>Sim SN <u>0R7111</u> Lot No. <u>16060</u> Rep AC <u>0.300</u>   | <u>X</u>        |
| G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.<br>(Attach test record.)<br>Sim SN <u>G10329</u> Lot No. <u>1054</u> Solution <u>0.10AC + 0.05%</u><br><span style="margin-left: 150px;">Acetone</span> | <u>X</u>        |
| H. RFI check. Run CMS mode. Key radio during CMS test.<br>(Attach test record.)  | <u>X</u>        |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.) <u>0.080</u> AC  |                 |

Lot No. 17316080A2 Cyl No. 27 Exp. Date 5 Aug 18

Test 1 0.081 Test 4 0.082 Test 7 0.083

Test 2 0.082 Test 5 0.082 Test 8 0.083

Test 3 0.082 Test 6 0.082 Test 9 0.082

Average 0.082

J. Tank pressure, Level 3 Func D Sub G

P: <sup>897 ~~RAW~~</sup> 950 psi Reg <sup>900 ~~RAW~~</sup> 980 psi

K. Remarks/Maintenance record: Annual Inspection.

17 Apr 17  
Date

17 Apr 17  
Date

Roberta Mueger-Nimmo  
Field Inspector's Signature

Djb Shanaver  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000       SN 80-006688  
Location = TOXL                8164.14.00 09/16  
04/17/2017                      10:58

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?         No  
Display Volume?             No  
Disable On Memfull?         Yes  
# of Print Copies?          1  
Select Std (D/W/I)?         Dry  
Standard Value?             0.080  
Standard Lot #?             17316080A2  
Standard Cyl #?             27  
Standard Expiration?        08/05/2018  
Oper No?                    888888

Flow Cal. Date:             09/07/2016  
Slope                        698  
Intercept                   -642832

IR Calibration Date:        09/07/2016

	3um	9um
0th Coef(*100):	-10822	-23152
1st Coef(*100):	262468	134492
2nd Coef(*100):	3796	1436
H2O adj(mg/l*10k):	665	568

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: *Print Test*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006688  
Location = TOXL      8164.14.00 09/16  
04/17/2017      10:58

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:59
02 Std. Sol.	0.021	11:00
03 Room Air	0.000	11:00
04 Std. Sol.	0.020	11:01
05 Room Air	0.000	11:02
06 Std. Sol.	0.021	11:02
07 Room Air	0.000	11:03

08 Sim Temp = 34.0°C

Simul Ser No = DR3846  
Std Sol No = 16200  
County = 08      Oper No. = 888888

  
Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: *Low AC*  
*0.020 AC*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006688  
Location = TOXL                    8164.14.00 09/16  
04/17/2017                            11:04

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:04
02 Std. Sol.	0.304	11:05
03 Room Air	0.000	11:06
04 Std. Sol.	0.303	11:06
05 Room Air	0.000	11:07
06 Std. Sol.	0.304	11:08
07 Room Air	0.000	11:08

08 Sim Temp = 34.0°C

Simul Ser No = DR7111  
Std Sol No = 16060  
County = 08                            Oper No. = 888888

  
Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: High AC  
          0.300 AC


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006688  
Location = TOXL      8164.14.00 09/16  
04/17/2017      11:12

Test	AC	Time
01 Room Air	0.000	11:13
02 *Subject Test	INT*	11:14
03 Room Air	0.000	11:14

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 27  
Expiration Date = 08/05/2018  
County = 08      Oper No. = 888888

  
Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: Interferent Check

0.10 AC + 0.05% Acetone

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006688  
Location = TOXL                8164.14.00 09/16  
04/17/2017                        11:14

Test	AC	Time
01 Diagnostic	OK	11:15
02 Room Air	RFI*	11:15
03 Room Air	0.000	11:16

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                Weight = 150  
Test = DUI                        Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 27  
Expiration Date = 08/05/2018  
County = 08                        Oper No. = 888888

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: *RFI Check*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006688  
Location = TOXL      8164.14.00 09/16  
04/17/2017      11:17

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:18
02 Std. Gas	0.081	11:18
03 Room Air	0.000	11:19
04 Std. Gas	0.082	11:19
05 Room Air	0.000	11:19
06 Std. Gas	0.082	11:20
07 Room Air	0.000	11:20

Lot No = 17316080A2  
Cyl No = 27  
Exp Date = 08/05/2018  
County = 08      Oper No. = 888888

Roberta Grieger-Nimmo  
Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: Calibration Check  
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006688  
Location = TOXL      8164.14.00 09/16  
04/17/2017      11:21

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:21
02 Std. Gas	0.082	11:21
03 Room Air	0.000	11:22
04 Std. Gas	0.082	11:22
05 Room Air	0.000	11:23
06 Std. Gas	0.082	11:23
07 Room Air	0.000	11:24

Lot No = 17316080A2  
Cyl No = 27  
Exp Date = 08/05/2018  
County = 08      Oper No. = 888888

  
Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: *Calibration Check*  
*0.080 AC*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006688  
Location = TOXL      8164.14.00 09/16  
04/17/2017      11:24

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:25
02 Std. Gas	0.083	11:25
03 Room Air	0.000	11:26
04 Std. Gas	0.083	11:26
05 Room Air	0.000	11:26
06 Std. Gas	0.082	11:27
07 Room Air	0.000	11:27

Lot No = 17316080A2  
Cyl No = 27  
Exp Date = 08/05/2018  
County = 08      Oper No. = 888888

  
Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: Calibration Check  
0.080 AC