

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-006510

Location: HETT Adams CO SO

- | | Check When Done |
|---|-----------------|
| A. Is the breath tube heated? | <u>X</u> |
| B. Is the diagnostic check complete? | <u>X</u> |
| C. Is the time, date, and year correct (re-set if necessary)? | <u>X</u> |
| D. Print test. (Attach test record.) | <u>X</u> |
| E. Low AC. Use ≤ 0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5113</u> Lot No. <u>201507F</u> Rep AC <u>0.020 AC</u> | <u>X</u> |
| F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR3378</u> Lot No. <u>16060</u> Rep AC <u>0.300 AC</u> | <u>X</u> |
| G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.
(Attach test record.)
Sim SN <u>DR3803</u> Lot No. <u>1CS5</u> Solution <u>0.100AC + 0.05% Acetone</u> | <u>X</u> |
| H. RFI check. Run CMS mode. Key radio during CMS test.
(Attach test record.) | <u>X</u> |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) <u>0.080 AC</u> | |

Lot No. <u>09216080A3</u>	Cyl No. <u>40</u>	Exp. Date <u>5 May 18</u>
Test 1 <u>0.081</u>	Test 4 <u>0.081</u>	Test 7 <u>0.081</u>
Test 2 <u>0.081</u>	Test 5 <u>0.081</u>	Test 8 <u>0.081</u>
Test 3 <u>0.081</u>	Test 6 <u>0.081</u>	Test 9 <u>0.081</u>
Average <u>0.081</u>		

J. Tank pressure, Level 3 Func D Sub G P: 993 psi Reg 1000 psi

K. Remarks/Maintenance record: Annual Inspection,

30 May 17
Date

Rebecca Suega-Timm
Field Inspector's Signature

31 May 17
Date

Deb Shanaver
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006510
Location = HETT 8164.14.00 09/16
05/30/2017 10:05

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:05
02 Std. Sol.	0.019	10:06
03 Room Air	0.000	10:07
04 Std. Sol.	0.019	10:07
05 Room Air	0.000	10:08
06 Std. Sol.	0.019	10:09
07 Room Air	0.000	10:09


08 Sim Temp = 34.0°C

Simul Ser No = DR5133

Std Sol No = 201507F

County = 01

Oper No. = 888888



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: low AC
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006510
Location = HETT 8164.14.00 09/16
05/30/2017 10:23

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:24
02 Std. Sol.	0.299	10:26
03 Room Air	0.000	10:27
04 Std. Sol.	0.299	10:28
05 Room Air	0.000	10:28
06 Std. Sol.	0.300	10:29
07 Room Air	0.000	10:30


08 Sim Temp = 34.0°C

Simul Ser No = DR3378

Std Sol No = 16060

County = 01

Oper No. = 888888



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: High AC
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006510
Location = HETT 8164.14.00 09/16
05/30/2017 10:36

Test	AC	Time
01 Diagnostic	OK	10:37
02 Room Air	0.000	10:38
03 *Subject Test	RFI*	10:38
04 Room Air	0.000	10:38

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS821456
Lot No = 09216080A3
Cyl No = 40
Expiration Date = 05/05/2018
County = 01 Oper No. = 888888

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.


Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: *RFI Check*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006510
Location = HETT 8164.14.00 09/16
05/30/2017 10:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:40
02 Std. Gas	0.081	10:40
03 Room Air	0.000	10:41
04 Std. Gas	0.081	10:41
05 Room Air	0.000	10:41
06 Std. Gas	0.081	10:42
07 Room Air	0.000	10:42

Lot No = 09216080A3
Cyl No = 40
Exp Date = 05/05/2018
County = 01 Oper No. = 888888



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: Calibration Check
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006510
Location = HETT 8164.14.00 09/16
05/30/2017 10:43

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:43
02 Std. Gas	0.081	10:44
03 Room Air	0.000	10:44
04 Std. Gas	0.081	10:44
05 Room Air	0.000	10:45
06 Std. Gas	0.081	10:45
07 Room Air	0.000	10:46

Lot No = 09216080A3
Cyl No = 40
Exp Date = 05/05/2018
County = 01 Oper No. = 888888


Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: Calibration Check
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006510
Location = HETT 8164.14.00 09/16
05/30/2017 10:46

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:47
02 Std. Gas	0.081	10:47
03 Room Air	0.000	10:48
04 Std. Gas	0.081	10:48
05 Room Air	0.000	10:49
06 Std. Gas	0.081	10:49
07 Room Air	0.000	10:49

Lot No = 09216080A3
Cyl No = 40
Exp Date = 05/05/2018
County = 01 Oper No. = 888888


Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: Calibration Check
0.080AC

Form 106-I8000