



ETHANOL GAS STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <u>STEVEN NAGEL #071109</u>		
Location <u>CARRINGTON PD</u>	<u>1010 MAIN ST CARRINGTON, NC 28421</u>	Intoxilyzer® Serial Number <u>80-006512</u>
Gas Lot Number <u>20016080A2</u>	Gas Cylinder Number <u>44</u>	Gas Expiration Date <u>09/25/2018</u>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 A. 0.080 AC
 B. 0.081 AC
 C. 0.081 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <u>Steven Nagel</u>	Date <u>09/27/2017</u>
Reviewed By (Crime Laboratory Use Only) <u>Roberta Yancy - Thomas</u>	Date <u>10/2/2017</u>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006512
Location = CARR 8164.14.00 09/16
09/27/2017 19:13

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:14
02 Std. Gas	0.080	19:15
03 Room Air	0.000	19:15
04 Std. Gas	0.081	19:16
05 Room Air	0.000	19:16
06 Std. Gas	0.081	19:17
07 Room Air	0.000	19:17

Lot No = 20016080A2
Cyl No = 44
Exp Date = 09/05/2018
County = 16 Oper No. = 071109



Operator Signature
STEVEN NAGEL

Remarks: