

**INTOXILYZER® 8000 INSPECTION**

Intoxilyzer® 8000 Serial Number: 80-006681

Location: TOXL

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 426 psi Regulator: 450 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check

- a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. Ser #: MP3002 Lot #: 17160 Exp. Date: 6/8/2019
- b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.400 AC  
Sim. Ser #: MP3061 Lot #: 201605F Exp. Date: 5.19.2018

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR7089 Lot #: ICS 5 Exp. Date: 5/15/2019  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 03417080A3 Cyl #: 27 Exp. Date: 3/5/2019

Test 1	<u>0.081</u>	Test 4	<u>0.081</u>	Test 7	<u>0.081</u>
Test 2	<u>0.081</u>	Test 5	<u>0.081</u>	Test 8	<u>0.080</u>
Test 3	<u>0.081</u>	Test 6	<u>0.081</u>	Test 9	<u>0.081</u>
Average	<u>0.081</u>				

C. Remarks/Maintenance: ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

Charles E. Edin  
Breath Analyst Signature

4.2.2018  
Date

Roberta Mungo - Nunno  
Reviewed by

3 Apr 18  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006681  
Location = TOXL                    8164.14.00 09/16  
04/02/2018                            15:06

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&\*( )\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&\*( )\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            03417080A3  
Standard Cyl #?            27  
Standard Expiration?       03/05/2019  
Oper No?                    666666

Flow Cal. Date:            09/09/2016  
Slope                      694  
Intercept                  -677594

IR Calibration Date:       09/09/2016  
                                  3um                            9um

-----  
0th Coef(\*100):        -10880                    -18848  
1st Coef(\*100):        255177                   129254  
2nd Coef(\*100):        4336                     1756  
H2O adj(mg/l\*10k):    468                      501

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:            PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006681  
Location = TOXL      8164.14.00 09/16  
04/02/2018      15:18

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:18
02 Std. Sol.	0.019	15:19
03 Room Air	0.000	15:20
04 Std. Sol.	0.019	15:20
05 Room Air	0.000	15:21
06 Std. Sol.	0.020	15:22
07 Room Air	0.000	15:22

08 Sim Temp = 34.0°C

Simul Ser No = MP3002

Std Sol No = 17160

County = 08

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Low AC*

*0.020 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006681  
Location = TOXL      8164.14.00 09/16  
04/02/2018      15:23

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:23
02 Std. Sol.	0.409	15:24
03 Room Air	0.000	15:24
04 Std. Sol.	0.408	15:25
05 Room Air	0.000	15:26
06 Std. Sol.	0.407	15:26
07 Room Air	0.000	15:27

08 Sim Temp = 34.0°C

Simul Ser No = MP3061  
Std Sol No = 201605F  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.400AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006681  
Location = TOXL                8164.14.00 09/16  
04/02/2018                        15:12

Test	AC	Time
01 Room Air	0.000	15:13
02 *Subject Test	INT*	15:13
03 Room Air	0.000	15:14

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                Weight = 150  
Test = DUI                        Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 03417080A3  
Cyl No = 27  
Expiration Date = 03/05/2019  
County = 08                        Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006681  
Location = TOXL      8164.14.00 09/16  
04/02/2018      15:03

Test	AC	Time
01 Diagnostic	OK	15:04
02 Room Air	RFI*	15:04
03 Room Air	0.000	15:05

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 03417080A3  
Cyl No = 27  
Expiration Date = 03/05/2019  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*RFI CHECK*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006681  
Location = TOXL      8164.14.00 09/16  
04/02/2018      15:07

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:08
02 Std. Gas	0.081	15:08
03 Room Air	0.000	15:09
04 Std. Gas	0.081	15:09
05 Room Air	0.000	15:09
06 Std. Gas	0.081	15:10
07 Room Air	0.000	15:10

Lot No = 03417080A3  
Cyl No = 27  
Exp Date = 03/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006681  
Location = TOXL      8164.14.00 09/16  
04/02/2018      15:14

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:14
02 Std. Gas	0.081	15:15
03 Room Air	0.000	15:15
04 Std. Gas	0.081	15:16
05 Room Air	0.000	15:16
06 Std. Gas	0.081	15:17
07 Room Air	0.000	15:17

Lot No = 03417080A3  
Cyl No = 27  
Exp Date = 03/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006681  
Location = TOXL      8164.14.00 09/16  
04/02/2018      15:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:28
02 Std. Gas	0.081	15:28
03 Room Air	0.000	15:29
04 Std. Gas	0.080	15:29
05 Room Air	0.000	15:30
06 Std. Gas	0.081	15:30
07 Room Air	0.000	15:31

Lot No = 03417080A3  
Cyl No = 27  
Exp Date = 03/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC