

**INTOXILYZER® 8000 INSPECTION**

Intoxilyzer® 8000 Serial Number: 80-006673 Location: TOXL

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 879 psi Regulator: 900 psi

B. Tests (Sign and attach test records):

- Configure simulator for the following tests (Level 1,S).
1. Wet Calibration Check
    - a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. Ser #: MP3002 Lot #: 17160 Exp. Date: 6.8.19
    - b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.400 AC  
Sim. Ser #: MP3061 Lot #: 201605F Exp. Date: 5.19.18

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR7089 Lot #: ICS #5 Exp. Date: 5.15.19  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 19817080A4 Cyl #: 21 Exp. Date: 9.5.19

Test 1	<u>0.083</u>	Test 4	<u>0.082</u>	Test 7	<u>0.082</u>
Test 2	<u>0.082</u>	Test 5	<u>0.082</u>	Test 8	<u>0.082</u>
Test 3	<u>0.082</u>	Test 6	<u>0.082</u>	Test 9	<u>0.082</u>
Average <u>0.082</u>					

C. Remarks/Maintenance: ANNUAL INSPECTION, REPLACED PRESSURE GAUGE ASSEMBLY, TANK SENSOR CALIBRATION, REPLACED SIMULATOR RETURN O-RING.

Instrument is acceptable to be used in the field.

Chalala E. Eder  
Breath Analyst Signature

4.9.2018  
Date

Roberta Yunque Nimmo  
Reviewed by

10 Apr 18  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000           SN 80-006673  
Location = TOXL                   8164.14.00 09/16  
04/09/2018                           16:22

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:           Enabled  
Start Test Sequence:       DABACABA  
Display Prelim Rslt?       Yes  
Display Third Digit?       Yes  
Inhib Printer(Y/N)?       No  
Display Volume?           No  
Disable On Memfull?       Yes  
# of Print Copies?         1  
Select Std (D/W/I)?       Dry  
Standard Value?           0.080  
Standard Lot #?           19817080A4  
Standard Cyl #?           21  
Standard Expiration?      09/05/2019  
Oper No?                   666666

Flow Cal. Date:           03/30/2017  
Slope                      697  
Intercept                 -667725

IR Calibration Date:      03/30/2017  
                            3um                   9um

0th Coef(*100):	-9136	-17613
1st Coef(*100):	267126	133333
2nd Coef(*100):	5438	1698
H2O adj(mg/l*10k):	684	589

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006673  
Location = TOXL      8164.14.00 09/16  
04/09/2018      16:23

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:23
02 Std. Sol.	0.019	16:24
03 Room Air	0.000	16:25
04 Std. Sol.	0.019	16:25
05 Room Air	0.000	16:26
06 Std. Sol.	0.019	16:27
07 Room Air	0.000	16:27

08 Sim Temp = 34.0°C

Simul Ser No = MP3002  
Std Sol No = 17160  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC  
0.020 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006673  
Location = TOXL      8164.14.00 09/16  
04/09/2018      16:28

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:29
02 Std. Sol.	0.396	16:29
03 Room Air	0.000	16:30
04 Std. Sol.	0.400	16:31
05 Room Air	0.000	16:31
06 Std. Sol.	0.399	16:32
07 Room Air	0.000	16:32

08 Sim Temp = 34.0°C

Simul Ser No = MP3061

Std Sol No = 201605F

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

HIGH AC  
0.400 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006673  
Location = TOXL      8164.14.00 09/16  
04/09/2018      16:34

Test	AC	Time
01 Room Air	0.000	16:34
02 *Subject Test	INT*	16:35
03 Room Air	0.000	16:35

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 19817080A4  
Cyl No = 21  
Expiration Date = 09/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10 AE + 0.05% ACETONE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006673  
Location = TOXL      8164.14.00 09/16  
04/09/2018      16:38

Test	AC	Time
01 Diagnostic	OK	16:39
02 Room Air	0.000	16:40
03 *Subject Test	RFI*	16:40
04 Room Air	0.000	16:40

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 19817080A4  
Cyl No = 21  
Expiration Date = 09/05/2019  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*RFI CHECK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006673  
Location = TOXL      8164.14.00 09/16  
04/09/2018      16:43

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:44
02 Std. Gas	0.083	16:44
03 Room Air	0.000	16:45
04 Std. Gas	0.082	16:45
05 Room Air	0.000	16:46
06 Std. Gas	0.082	16:46
07 Room Air	0.000	16:47

Lot No = 19817080A4  
Cyl No = 21  
Exp Date = 09/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006673  
Location = TOXL      8164.14.00 09/16  
04/09/2018      16:47

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:48
02 Std. Gas	0.082	16:48
03 Room Air	0.000	16:49
04 Std. Gas	0.082	16:49
05 Room Air	0.000	16:49
06 Std. Gas	0.082	16:50
07 Room Air	0.000	16:50

Lot No = 19817080A4  
Cyl No = 21  
Exp Date = 09/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006673  
Location = TOXL      8164.14.00 09/16  
04/09/2018      16:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:51
02 Std. Gas	0.082	16:51
03 Room Air	0.000	16:52
04 Std. Gas	0.082	16:52
05 Room Air	0.000	16:53
06 Std. Gas	0.082	16:53
07 Room Air	0.000	16:54

Lot No = 19817080A4  
Cyl No = 21  
Exp Date = 09/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000