

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-006507

Location: WATF  
(McKENZIE COUNTY S.O.)

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 959 psi Regulator: 950 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check
  - a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. Ser #: DR3377 Lot #: 201610A Exp. Date: 10-4-18
  - b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.300 AC  
Sim. Ser #: DR5191 Lot #: 17350 Exp. Date: 10-11-19

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR3451 Lot #: ICS#5 Exp. Date: 5-15-19  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 19817080A4 Cyl #: 18 Exp. Date: 9-5-2019

Test 1	<u>0.079</u>	Test 4	<u>0.079</u>	Test 7	<u>0.080</u>
Test 2	<u>0.079</u>	Test 5	<u>0.079</u>	Test 8	<u>0.080</u>
Test 3	<u>0.079</u>	Test 6	<u>0.079</u>	Test 9	<u>0.079</u>
Average	<u>0.079</u>				

C. Remarks/Maintenance: ANNUAL INSPECTION, REPLACED EXTERNAL PRINTER

Instrument is acceptable to be used in the field.

Charles E. Ed  
Breath Analyst Signature

5-7-2018  
Date

Roberta Muegge-Nimmer  
Reviewed by

14 May 18  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000       SN 80-006507  
Location = WATF                8164.13.00 06/09  
05/07/2018                     11:53

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?         No  
Display Volume?             No  
Disable On Memfull?         Yes  
# of Print Copies?          1  
Select Std (D/W/I)?         Dry  
Standard Value?             0.080  
Standard Lot #?             19817080A4  
Standard Cyl #?             18  
Standard Expiration?        09/05/2019  
Oper No?                    666666

Flow Cal. Date:             08/20/2015  
Slope                        659  
Intercept                   -658315

IR Calibration Date:        08/20/2015

	3um	9um
0th Coef(*100):	-9452	-11722
1st Coef(*100):	258808	139588
2nd Coef(*100):	4476	1674
H2O adj(mg/l*10k):	499	376

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006507  
Location = WATF      8164.13.00 06/09  
05/07/2018      12:13

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:13
02 Std. Sol.	0.018	12:14
03 Room Air	0.000	12:14
04 Std. Sol.	0.018	12:15
05 Room Air	0.000	12:16
06 Std. Sol.	0.018	12:16
07 Room Air	0.000	12:17

08 Sim Temp = 34.0°C

Simul Ser No = DR3377  
Std Sol No = 201610A  
County = 27      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC  
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006507  
Location = WATF      8164.13.00 06/09  
05/07/2018      12:22

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:23
02 Std. Sol.	0.306	12:23
03 Room Air	0.000	12:24
04 Std. Sol.	0.305	12:25
05 Room Air	0.000	12:25
06 Std. Sol.	0.305	12:26
07 Room Air	0.000	12:27

08 Sim Temp = 34.0°C

Simul Ser No = DR5191

Std Sol No = 17350

County = 27

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

HIGH AC

0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006507  
Location = WATF                8164.13.00 06/09  
05/07/2018                        12:02

Test	AC	Time
01 Room Air	0.000	12:03
02 *Subject Test	INT*	12:04
03 Room Air	0.000	12:04

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 19817080A4  
Cyl No = 18  
Expiration Date = 09/05/2019  
County = 27                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Interferent Check*  
*0.10 AC + 0.05% ACETONE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006507  
Location = WATF                8164.13.00 06/09  
05/07/2018                        12:00

Test	AC	Time
01 Diagnostic	OK	12:01
02 Room Air	0.000	12:02
03 *Subject Test	RFI*	12:02
04 Room Air	0.000	12:02

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 19817080A4  
Cyl No = 18  
Expiration Date = 09/05/2019  
County = 27                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:            *RFI CHECK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006507  
Location = WATF      8164.13.00 06/09  
05/07/2018      12:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:05
02 Std. Gas	0.079	12:05
03 Room Air	0.000	12:06
04 Std. Gas	0.079	12:06
05 Room Air	0.000	12:07
06 Std. Gas	0.079	12:07
07 Room Air	0.000	12:08

Lot No = 19817080A4  
Cyl No = 18  
Exp Date = 09/05/2019  
County = 27      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006507  
Location = WATF      8164.13.00 06/09  
05/07/2018      12:08

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:08
02 Std. Gas	0.079	12:09
03 Room Air	0.000	12:09
04 Std. Gas	0.079	12:10
05 Room Air	0.000	12:10
06 Std. Gas	0.079	12:11
07 Room Air	0.000	12:11

Lot No = 19817080A4  
Cyl No = 18  
Exp Date = 09/05/2019  
County = 27      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC  
Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006507  
Location = WATF      8164.13.00 06/09  
05/07/2018      12:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:28
02 Std. Gas	0.080	12:28
03 Room Air	0.000	12:29
04 Std. Gas	0.080	12:29
05 Room Air	0.000	12:29
06 Std. Gas	0.079	12:30
07 Room Air	0.000	12:30

Lot No = 19817080A4  
Cyl No = 18  
Exp Date = 09/05/2019  
County = 27      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 Ac

Form 106-I8000