



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Thompson, Chad A</i>		
Location <i>Cass Jail</i>	Intoxilyzer® Serial Number <i>80-006668</i>	
Gas Lot Number <i>17316080A2</i>	Gas Cylinder Number <i>34</i>	Gas Expiration Date <i>8-5-18</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0.082* AC
 - B. *0.083* AC
 - C. *0.083* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i> <i>3561</i>	Date <i>5-25-18</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Yungue-Nimmo</i>	Date <i>5-29-18</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006668
Location = CASS 8164.14.00 09/16
05/25/2018 01:13

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	01:14
02 Std. Gas	0.082	01:14
03 Room Air	0.000	01:14
04 Std. Gas	0.083	01:15
05 Room Air	0.000	01:15
06 Std. Gas	0.083	01:16
07 Room Air	0.000	01:16

Lot No = 17316080A2
Cyl No = 34
Exp Date = 08/05/2018
County = 09 Oper No. = 050910

 3561

Operator Signature
CHAD THOMPSON

Remarks: Installed low PSI Tank, Ran AEA
Test. No subjects were Ran. Then
Form 106-I8000 Installed new cylinder