



# INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION  
Toxicology Section/Breath Alcohol Program  
SFN 59281 (10-2016)

|  |                                    |
|--|------------------------------------|
| Serial Number<br><i>8000 6681</i>  | Instrument Location<br><i>FINL</i> |
| Reason for Install/Repair  |                                    |
| <input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change<br><input checked="" type="checkbox"/> Other (Specify) <u><i>Reinstall after re-certification of field inspector</i></u> |                                    |

### Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- NA* 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): *NA*
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

|   |                       |
|---|-----------------------|
| Field Inspector Signature<br><i>[Signature]</i> | Date<br><i>7/5/18</i> |
|---|-----------------------|

### Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

|   |                                 |
|---|---------------------------------|
| Reviewed/Certified by<br><i>Beata Hugel-Nimmo</i> | Certified Date<br><i>7/6/18</i> |
|---|---------------------------------|

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006681  
Location = FINL      8164.14.00 09/16  
07/05/2018      12:59

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%\$^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%\$^&\*()\_+?

Current Instrument Setup


Data Entry Mode:      Enabled  
Start Test Sequence:      DABACABA  
Display Prelim Rslt?      Yes  
Display Third Digit?      Yes  
Inhib Printer(Y/N)?      No  
Display Volume?      No  
Disable On Memfull?      Yes  
# of Print Copies?      1  
Select Std (D/W/I)?      Dry  
Standard Value?      0.080  
Standard Lot #?      34917080A3  
Standard Cyl #?      1  
Standard Expiration?      02/05/2020  
Oper No?      061123

Flow Cal. Date:      09/09/2016  
Slope      694  
Intercept      -677594

IR Calibration Date:      09/09/2016

|                    | 3um    | 9um    |
|--------------------|--------|--------|
| 0th Coef(*100):    | -10880 | -18848 |
| 1st Coef(*100):    | 255177 | 129254 |
| 2nd Coef(*100):    | 4336   | 1756   |
| H2O adj(mg/l*10k): | 468    | 501    |

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
\_\_\_\_\_  
Operator Signature  
WESLEY STRAIGHT

Remarks:


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006681  
Location = FINL      8164.14.00 09/16  
07/05/2018      13:00

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 13:01 |
| 02 Std. Gas | 0.078 | 13:01 |
| 03 Room Air | 0.000 | 13:02 |
| 04 Std. Gas | 0.079 | 13:02 |
| 05 Room Air | 0.000 | 13:02 |
| 06 Std. Gas | 0.078 | 13:03 |
| 07 Room Air | 0.000 | 13:03 |

Lot No = 34917080A3  
Cyl No = 1  
Exp Date = 02/05/2020  
County = 46      Oper No. = 061123



Operator Signature  
WESLEY STRAIGHT

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006681  
Location = FINL      8164.14.00 09/16  
07/05/2018      13:04

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | RFI*  | 13:04 |
| 02 Room Air | 0.000 | 13:05 |

\*Invalid Test  
Inhibited - RFI

Lot No = 34917080A3  
Cyl No = 1  
Exp Date = 02/05/2020  
County = 46      Oper No. = 061123

  
\_\_\_\_\_  
Operator Signature  
WESLEY STRAIGHT

Remarks: