

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-006687

A. GENERAL Setup and Checks:

- Power On
 - 5 min. countdown during initial warm-up
 - 1 min. countdown to begin test after standby mode
- Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:
 - Level 1 (Operator) Level 2 (Field Inspector) Level 3 (Lab Personnel)
- Printer Setup (Level 3,S,P)
- Forms Load (CMI, Inc., CoBrA Software)

B. LEVEL 1 Setup and Checks:

- Internal Printer Test and Paper Feed (Level 1,P). Attach record.
- External Printer Test (Level 1,P). Attach record.
- Barcode Reader (Level 1,B):
 - Operator ID ND License R key: Refusal. Attach record.

C. LEVEL 2 Setup and Checks:

- Continuous Air Blank (Level 2,A)
- Diagnostic (Level 2,D). Attach record.
- Prelim Data Entry (Level 2,E)
- Ser Number/Version Info (Level 2,V):
 - Inst Ser#: 80-006687 DSP Ver: 1362.14
 - Core Ver: 8164.14 CPLD Ver: 7

D. LEVEL 3 Setup and Checks:

- DVM Monitor (Level 3,D,D): 3: 12475 s: 0 C: 47.00
9: 12856 s: 1 B: 445 0
(RFI indicator with keyed radio): 1
- Function (F) Keys: : F5 – Tone F9 - Fan
 F6 – Green LED F10 – Pump
 F7 – Red LED F11 – Solenoid
- Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 964 CF: 1.05
- Flow Monitor (Level 3,D,F): Flow Sensor, R: 224 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T): C: 47.0 A: 50.0 I: 0.0
B: 45.0 V: 0.0
6. Tank Monitor (Level 3,D,G): P (display): 566 psi P (gauge): 550 psi
7. General Setup (Level 3,S,G)
8. Configure Start Test (Level 3,S,T)
9. Network (Level 3,S,N)

E. TESTS (attach all test records):

Configure wet bath simulator standards before running tests (Level 1,S):

1. Zero AC (Level 1,C):
Sim ser #: DR3377 Lot #: milli Q Exp date: NA
Known value: 0.000 AC
Results: 0.000 AC 0.000 AC 0.000 AC
2. Low AC (Level 1,C):
Sim ser #: DR7111 Lot #: 17160 Exp date: 8 Jun 19
Known value: 0.020 AC
Results: 0.020 AC 0.020 AC 0.020 AC
3. High AC (Level 1,C):
Sim ser #: DR3375 Lot #: 17350 Exp date: 11 Oct 19
Known value: 0.300 AC
Results: 0.296 AC 0.298 AC 0.297 AC

Configure dry gas standard before running tests (Level 1,S):

4. Ethanol + 0.05% Acetone (Level 1,B):
Sim ser #: DR7089 Lot #: ICSS Exp date: 15 May 19
Known value: 0.10 AC + 0.05% Acetone
 Display reads "Interferent Detect"
5. Acetone (Level 1,B):
Sim ser #: DR3846 Lot #: ACE1 Exp date: NA
Known value: 0.05% Acetone
 Display reads "Interferent Detect"
6. Methanol (Level 1,B):
Sim ser #: DR3455 Lot #: MET2 Exp date: NA
Known value: 0.10% Methanol
 Display reads "Interferent Detect"
7. Isopropanol (Level 1,B):
Sim ser #: DR5189 Lot #: IPA2 Exp date: NA
Known value: 0.10% Isopropanol
 Display reads "Interferent Detect"

- 8. Ambient Fail (Level 1,B): Display and test record read correct.
- 9. Invalid Sample (Level 1,B): Display and test record read correct.
- 10. Sequence Aborted (Level 1,B): Display and test record read correct.
- 11. ABA (Level 1,B): Test record reads correct.
- 12. Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- 13. RFI Detected (B,C, or CMS mode): Display and test record read correct.
- 14. Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 15. Deficient Samples (CMS mode): Display and test record read correct.
- 16. Improper Sample (CMS mode): Display and test record read correct.
- 17. Difference Too Great (CMS mode): Display and test record read correct.
- 18. CMS (CMS mode): Test record reads correct.
- 19. Stability Test (Level 3,T,S):

Ave: 0.0821 Std Dev: 0.0004 Rel Std Dev: 0.49%

Maintenance/Repairs/Notes: Initial inspection after receiving
instrument from CRT, ^{Inc.} warranty repair. Completed
flash update of software as instrument locked up
during initial inspection tests.

Instrument is acceptable to be used in the field.

Roberta Hueque Nimmis
Field Inspector's Signature

21 NOV 18
Date

Charles E. Eder
Reviewed by

21 Nov. 2018
Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/19/2018 14:55

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMN

24 characters per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMN

32 characters per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMN

42 characters per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMN

Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/L)? Wet
Standard Value? 0.100
Standard Solution #?
Sim Serial #?
Oper No? 888888

Flow Cal. Date: 10/05/2018
Slope 696
Intercept -609977

IR Calibration Date: 10/05/2018
3um 9um

0th Coef(*100): -11478 -18771
1st Coef(*100): 255256 132135
2nd Coef(*100): 4432 1676
H2O adj(ng/l*10k): 604 575

***** Test End *****

Roberta Grieger-Nimmo
Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks:

Internal Print Test

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/19/2018 14:55

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#%\$^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#%\$^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/I)? Wet
Standard Value? 0.100
Standard Solution #?
Sim Serial #?
Oper No? 888888

Flow Cal. Date: 10/05/2018
Slope 696
Intercept -609977

IR Calibration Date: 10/05/2018
 3um 9um

0th Coef(*100): -11478 -18771
1st Coef(*100): 255256 132135
2nd Coef(*100): 4432 1676
H2O adj(mg/l*10k): 604 575

***** Printer Test End *****



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: External Printer Test


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/19/2018 14:56

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 14:57 |
| 02 *Subject Test | REF* | 14:57 |
| 03 Room Air | 0.000 | 14:58 |

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = BARCODE READER
Dr. Lic. = ND/DIS821456
Simul Ser No =
Std Sol No =
County = 08 Oper No. = 888888



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: *Barcode Reader*

R Key Refusal

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/19/2018 14:58

DIAGNOSTIC

| | |
|---------------------------|------|
| Voltage/Current Test | Pass |
| RAM Test | Pass |
| EEPROM Checksum Test | Pass |
| Real Time Clock Test | Pass |
| DSP Test | Pass |
| Analytical Stability Test | Pass |
| Modem Test | Pass |
| Temp Regulation Test | Pass |

Oper No. = 888888

Roberta Grieger-Nimmo

Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: Diagnostic

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/19/2018 15:07

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 15:08 |
| 02 Std. Sol. | 0.000 | 15:09 |
| 03 Room Air | 0.000 | 15:09 |
| 04 Std. Sol. | 0.000 | 15:10 |
| 05 Room Air | 0.000 | 15:10 |
| 06 Std. Sol. | 0.000 | 15:11 |
| 07 Room Air | 0.000 | 15:12 |

08 Sim Temp = 34.0°C

Simul Ser No = DR3377
Std Sol No = MILLIQ
County = 08 Oper No. = 888888



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: *Zero AC*
0.000 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/19/2018 15:13

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 15:14 |
| 02 Std. Sol. | 0.020 | 15:14 |
| 03 Room Air | 0.000 | 15:15 |
| 04 Std. Sol. | 0.020 | 15:16 |
| 05 Room Air | 0.000 | 15:16 |
| 06 Std. Sol. | 0.020 | 15:17 |
| 07 Room Air | 0.000 | 15:17 |

08 Sim Temp = 34.0°C

Simul Ser No = DR7111
Std Sol No = 17160
County = 08 Oper No. = 888888


Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: Low AC
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/19/2018 15:18

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 15:19 |
| 02 Std. Sol. | 0.296 | 15:20 |
| 03 Room Air | 0.000 | 15:20 |
| 04 Std. Sol. | 0.298 | 15:21 |
| 05 Room Air | 0.000 | 15:22 |
| 06 Std. Sol. | 0.297 | 15:22 |
| 07 Room Air | 0.000 | 15:23 |


08 Sim Temp = 34.0°C

Simul Ser No = DR3375

Std Sol No = 17350

County = 08

Oper No. = 888888



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: High AC
0.300AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/19/2018 15:24

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 15:25 |
| 02 *Subject Test | INT* | 15:25 |
| 03 Room Air | 0.000 | 15:26 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 17618080A1
Cyl No = 28
Expiration Date = 08/05/2020
County = 08 Oper No. = 888888



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: *Interferent Detect*
0.100 AC + 0.05% Acetone

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/19/2018 15:26

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 15:27 |
| 02 *Subject Test | INT* | 15:28 |
| 03 Room Air | 0.000 | 15:28 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 17618080A1
Cyl No = 28
Expiration Date = 08/05/2020
County = 08 Oper No. = 888888



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: Interferent Detect
0.059% Acetone

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/19/2018 15:29

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 15:30 |
| 02 *Subject Test | INT* | 15:30 |
| 03 Room Air | 0.000 | 15:31 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 17618080A1
Cyl No = 28
Expiration Date = 08/05/2020
County = 08 Oper No. = 888888

Roberta Grieger Nimmo
Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: Interferent Detect
0.10% methanol

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/19/2018 15:31

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 15:32 |
| 02 *Subject Test | INT* | 15:33 |
| 03 Room Air | 0.000 | 15:33 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 17618080A1
Cyl No = 28
Expiration Date = 08/05/2020
County = 08 Oper No. = 888888


Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: Interferent Detect
0.10% Isopropanol

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDDAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/21/2018 08:14

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | AMB* | 08:15 |
| 02 Room Air | 0.000 | 08:16 |

*Invalid Test
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 03417080A3
Cyl No = 7
Expiration Date = 03/05/2019
County = 08 Oper No. = 888888


Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks:

Ambient fail

Form 106-18000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/21/2018 10:48

| Test | AC | Time |
|-------------------|-------|-------|
| 01 Room Air | 0.000 | 10:48 |
| 02 Invalid Sample | X.XXX | 10:49 |
| 03 Room Air | 0.000 | 10:49 |

*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 03417080A3
Cyl No = 7
Expiration Date = 03/05/2019
County = 08 Oper No. = 888888



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: Invalid Sample

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/21/2018 08:39

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | ABT* | 08:40 |
| 02 Room Air | 0.000 | 08:40 |

*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 03417080A3
Cyl No = 7
Expiration Date = 03/05/2019
County = 08 Oper No. = 888888



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: *Sequence Aborted*

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/21/2018 08:40

| Test | AC | Time |
|-------------------|-------|-------|
| 01 Room Air | 0.000 | 08:41 |
| 02 Subject Test 1 | 0.000 | 08:41 |
| 03 Room Air | 0.000 | 08:42 |
| 04 Reported AC | 0.000 | 08:41 |

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 03417080A3
Cyl No = 7
Expiration Date = 03/05/2019
County = 08 Oper No. = 888888


Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: ABA

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/21/2018 08:53

| Test | AC | Time |
|-------------------|--------|-------|
| 01 Diagnostic | OK | 08:54 |
| 02 Room Air | 0.000 | 08:54 |
| 03 Subject Test 1 | 0.000 | 08:55 |
| 04 Room Air | 0.000 | 08:57 |
| 05 Std. Gas | 0.090* | 08:58 |
| 06 Room Air | 0.000 | 08:59 |

*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 03417080A3
Cyl No = 7
Expiration Date = 03/05/2019
County = 08 Oper No. = 888888

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.


Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: Calibration Check out of tolerance

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc: Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/21/2018 09:00

| Test | AC | Time |
|------------------|-------|-------|
| 01 Diagnostic | OK | 09:01 |
| 02 Room Air | 0.000 | 09:01 |
| 03 *Subject Test | RFI* | 09:04 |
| 04 Room Air | 0.000 | 09:04 |

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 03417080A3
Cyl No = 7
Expiration Date = 03/05/2019
County = 08 Oper No. = 888888

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.


Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: RFI detected

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/21/2018 09:05

| Test | AC | Time |
|------------------|--------|-------|
| 01 Diagnostic | OK | 09:06 |
| 02 Room Air | 0.000 | 09:06 |
| 03 *Subject Test | 0.000* | 09:09 |
| 04 Room Air | 0.000 | 09:13 |
| 05 Std. Gas | 0.082 | 09:14 |
| 06 Room Air | 0.000 | 09:16 |
| 07 *Subject Test | 0.000* | 09:19 |
| 08 Room Air | 0.000 | 09:20 |


No RFI Detected

*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 03417080A3
Cyl No = 7
Expiration Date = 03/05/2019
County = 08 Oper No. = 888888

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: Deficient Samples

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/21/2018 09:25

| Test | AC | Time |
|---------------|-------|-------|
| 01 Diagnostic | OK | 09:26 |
| 02 Room Air | IPS* | 09:26 |
| 03 Room Air | 0.000 | 09:26 |

*Improper Sample

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 03417080A3
Cyl No = 7
Expiration Date = 03/05/2019
County = 08 Oper No. = 888888

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.


Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: Improper Sample

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/21/2018 10:37

| Test | AC | Time |
|-------------------|--------|-------|
| 01 Diagnostic | OK | 10:37 |
| 02 Room Air | 0.000 | 10:38 |
| 03 Subject Test 1 | 0.000* | 10:38 |
| 04 Room Air | 0.000 | 10:40 |
| 05 Std. Gas | 0.082 | 10:42 |
| 06 Room Air | 0.000 | 10:43 |
| 07 Subject Test 2 | 0.244* | 10:44 |
| 08 Room Air | 0.000 | 10:45 |

*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 03417080A3
Cyl No = 7
Expiration Date = 03/05/2019
County = 08 Oper No. = 888888

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: *Difference too great*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/21/2018 11:01


| Test | AC | Time |
|-------------------|-------|-------|
| 01 Diagnostic | OK | 11:02 |
| 02 Room Air | 0.000 | 11:03 |
| 03 Subject Test 1 | 0.000 | 11:03 |
| 04 Room Air | 0.000 | 11:08 |
| 05 Std. Gas | 0.081 | 11:09 |
| 06 Room Air | 0.000 | 11:11 |
| 07 Subject Test 2 | 0.000 | 11:11 |
| 08 Room Air | 0.000 | 11:12 |
| 09 Reported AC | 0.000 | 11:03 |

Difference OK
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 03417080A3
Cyl No = 7
Expiration Date = 03/05/2019
County = 08 Oper No. = 888888

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: CMS test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006687
Location = TOXL 8164.14.00 09/16
11/21/2018 09:27

| Test | AC | Time |
|---------------------|--------|-------|
| 01 Room Air | 0.000 | 09:27 |
| 02 Std. Gas | 0.081 | 09:28 |
| 03 Room Air | 0.000 | 09:28 |
| 04 Std. Gas | 0.082 | 09:29 |
| 05 Room Air | 0.000 | 09:29 |
| 06 Std. Gas | 0.082 | 09:30 |
| 07 Room Air | 0.000 | 09:30 |
| 08 Std. Gas | 0.082 | 09:30 |
| 09 Room Air | 0.000 | 09:31 |
| 10 Std. Gas | 0.082 | 09:31 |
| 11 Room Air | 0.000 | 09:32 |
| 12 Std. Gas | 0.082 | 09:32 |
| 13 Room Air | 0.000 | 09:33 |
| 14 Std. Gas | 0.082 | 09:33 |
| 15 Room Air | 0.000 | 09:33 |
| 16 Std. Gas | 0.082 | 09:34 |
| 17 Room Air | 0.000 | 09:34 |
| 18 Std. Gas | 0.082 | 09:35 |
| 19 Room Air | 0.000 | 09:35 |
| 20 Std. Gas | 0.082 | 09:36 |
| 21 Room Air | 0.000 | 09:36 |
| 22 Std. Gas | 0.082 | 09:36 |
| 23 Room Air | 0.000 | 09:37 |
| 24 Std. Gas | 0.082 | 09:37 |
| 25 Room Air | 0.000 | 09:38 |
| 26 Std. Gas | 0.083 | 09:38 |
| 27 Room Air | 0.000 | 09:39 |
| 28 Std. Gas | 0.082 | 09:39 |
| 29 Room Air | 0.000 | 09:40 |
| 30 Std. Gas | 0.082 | 09:40 |
| 31 Room Air | 0.000 | 09:40 |
| 32 Std. Gas | 0.082 | 09:41 |
| 33 Room Air | 0.000 | 09:41 |
| 34 Std. Gas | 0.082 | 09:42 |
| 35 Room Air | 0.000 | 09:42 |
| 36 Std. Gas | 0.082 | 09:43 |
| 37 Room Air | 0.000 | 09:43 |
| 38 Std. Gas | 0.083 | 09:43 |
| 39 Room Air | 0.000 | 09:44 |
| 40 Std. Gas | 0.082 | 09:44 |
| 41 Room Air | 0.000 | 09:45 |
| 42 Cal Check Stats: | | |
| 43 Average | 0.0821 | |
| 44 Std Dev | 0.0004 | |
| 45 Rel Std Dev | 0.49% | |

Roberta Suegra-Nimm

Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: Gas Cyl Lot # 03417080A3
Cyl # 007
Exp Date 5 Mar 19

Form 106-I8000