

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location	
80-006676	135PD	
Reason for Install/Repair		
☐ Install After Receiving From Crime Laboratory ☐ Install Afte	er Location Change	
Other (Specify) New Field Inspector change o	ver from Pathenz to Sof. 1	Tion Bleth
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
☑ 3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location	on; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
✓ 6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	uring Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50)496, Form 120-G) and Place it by the I	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Age	ency.
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	xilyzer® 8000 Installation and Repair C	Checkout (SFN59281, Form
Field Inspector Signature		Date
1 - Blett		1-7-19
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on file County of Burleigh, North Dakota, is certified to be a true and correct cop	e at the Offiœ of Attorney General, Crin	
Reviewed/Certified By		Certified Date

FORM 104-G

1-9-19

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006675
Location = BSPD 8164.14.00 09/16
01/07/2019 13:04

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34917080A3 Standard Cyl #? 16 Standard Expiration? 02/05/2020 Oper No? 069823

Flow Cal. Date: 08/25/2016 Slope 701 Intercept -701324

IR Calibration Date: 08/25/2016

Oth Coef(*100): -8608 -21568
1st Coef(*100): 258710 128604
2nd Coef(*100): 2478 1279
H2O adj(mg/l*10k): 381 528

****** Printer Test End *******

Operator Signature TIMOTHY BLETH

Remarks: New Field Inspector

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BSPD
01/07/2019

Alcohol Analyzer SN 80-006675 8164.14.00 09/16 13:05

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:06
02 Std. Gas	0.080	13:06
03 Room Air	0.000	13:07
04 Std. Gas	0.080	13:07
05 Room Air	0.000	13:08
06 Std. Gas	0.081	13:08
07 Room Air	0.000	13:08

Lot No = 34917080A3

Cyl No = 16

Exp Date = 02/05/2020

County = 08

Oper No. = 069823

Operator Signature TIMOTHY BLETH

Remarks: ACA Test - New Field Inspector

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006675
Location = BSPD 8164.14.00 09/16
01/07/2019 13:09

Test	AC	Time
01 Room Air	RFI*	13:12
02 Room Air	0.000	13:12

*Invalid Test Inhibited - RFI

Sub Name = RFI, TEST NA
Sub DOB = 01/01/2001
Sub Sex = Unknown
Test = OTH

Weight = NA Cit = NA

Dr. Lic. = ND/NA Lot No = 34917080A3 Cyl No = 16

Expiration Date = 02/05/2020

County = 08

Oper No. = 069823

Operator Signature TIMOTHY BLETH

Remarks: RFI Test - New Field Inspector

Form 106-I8000