

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location	lina Me : An invest
80-00 66 76 Reason for Install/Repair	Meci Jajem PU	400 main Ave NPW Sale
	ter Location Change	
Other (Specify)	to Loodion ondingo	
Content (Openity)		
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxifyzer® 8000.		
☑ 3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Loca	tion; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
6. Run Tests: CEE 5.23.19 TESTS PERF	ORMED	
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode of	r Level 1, Function B or C; K	ey Radio During Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SFNs	50496, Form 120-G) and Plac	ce it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-	G) at the Intoxilyzer® Locatio	on at the Agency.
10. Send the Following to the Crime Laboratory: Completed In 104-G), Print Test, ACA Test, and RFI Test.	toxilyzer® 8000 Installation a	and Repair Checkout (SFN59281, Form
Field Inspector Signature		Date
		<b>9</b> 3-7-19
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to the date the Field Inspector performed the installation. This record on		
County of Burleigh, North Dakota, is certified to be a true and correct of		
Reviewed/Certified By	2	Certified Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006676
Location = NS 8164.14.00 09/16
03/07/2019 15:14

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
# of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 34917080A3
Standard Cyl #? 35
Standard Expiration? 02/05/2020
Oper No? 020611

Flow Cal. Date: 08/22/2016 Slope 710 Intercept -713681

TR Calibration Date: 08/22/2016
3um 9um

Oth Coef(\*100): -16667 -22327
1st Coef(\*100): 259190 132291
2nd Coef(\*100): 3393 1374

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

546

Operator Signature
TOM SHARP

H20 adj(mg/1\*10k): 611

Remarks:

Form 196-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = NS 03/07/2019

SN 80-006676 8164.14.00 09/16 15:15

## DRY CAL CHECK

Te	est		AC	Time
01	Room	Air	0.000	15:16
02	Std.	Gas	0.083	15:16
03	Room	Air	0.000	15:17
04	Std.	Gas	0.082	15:17
05	Room	Air	0.000	15:18
06	Std.	Gas	0.083	15:18
07	Room	Air	0.000	15:18

Lot No = 34917080A3

Cyl No = 35

Exp Date = 02/05/2020

County = 30

Oper No. = 020611

Operator Signature TOM SHARP

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-006676
Location = NS	8164.14.00 09/16
03/07/2019	15:19

Test	AC	Time
01 Room Air	0.000	15:21
02 *Subject Test	RFI*	15:21
03 Room Air	0.000	15:22

\*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST TEST

Sub DOB = 01/01/2001

Sub Sex = Male Weight = 200Test = DUL Cit = 000000

Dr. Lic. = ND/TES010000

Lot No = 34917080A3

Cyl No = 35

Expiration Date = 02/05/2020

County = 30 Oper No. = 020611

Operator Signature
TOM SHARP

Remarks:

Form 106-I8000