### OFFICE OF ATTORNEY GENERAL **CRIME LABORATORY DIVISION**

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXII	$\_YZER^{@}$	8000	INSF	ECTI	ON
---------	--------------	------	------	------	----

Intoxilyzer® 8000 Serial Number: 80-006675 Location: Fort CEE  A. General Setup and Checks:  1. □ Diagnostic Tests Pass and Instrument in Ready Mode 2. □ Breath tube heated 3. □ Date, time and location code (Level 2,E). Re-set if necessary. 4. □ Print test (Level 1,P). Sign and attach test record. 5. □ Tank monitor (Level 3,D,G). Display and Regulator within 50 psi. a. Display: □ psi Regulator: □ psi  B. Tests (Sign and attach test records): □ Configure simulator for the following tests (Level 1,S). 1. □ Wet Calibration Check a. □ Low AC (Level 1,C): Known Value ≤ 0.03 AC: □ □ 2 □ AC Sign. Ser # □ PR 3575 Lot # □ 7 1 □ □ Exp. Date: □ 0.01 □ 1.01  □ Configure dry gas standard for the following tests (Level 1,S). 2. □ Interferent Check (Level 1,B): Known Value: □ 0.00 AC + 0.05% Acetone Sign. Ser # □ PR 361 Lot # □ 1 S □ Exp. Date: □ 0.01 □ 1.01  □ Display reads "Interferent Detect" 3. □ PR FI Check (CMS Mode) □ Display reads "RFI Detect" 4. □ Dry Calibration Check (Level 1,C): Known Value □ 0.00 AC Gas Cyl Lot # 34 □ 1 7 08 AC Cyl # □ Exp. Date: □ 0.00 Test 1 0.00 Test 2 0.00 Test 3 0.00 Test 5 0.00 Test 8 0.00 Test 9 0.00 Test 3 0.00 Test 5 0.00 Test 6 0.00 Test 9 0.00 Test		logy Section/Breath Alcohol Program zer® 8000 Inspection
1.		INTOXILYZER® 8000 INSPECTION
1.	Intoxil	yzer® 8000 Serial Number: 80-006675 Location: 70x P
Configure simulator for the following tests (Level 1,S).  1. Wet Calibration Check a. ✓ Low AC (Level 1,C): Known Value ≤ 0.03 AC: ○ ○ ○ ○ ○ ○ AC Sim. Ser #: ○ ○ ○ ○ ○ ○ ○ ○ ○ ○ ○ ○ ○ ○ ○ ○ ○ ○	Α.	<ol> <li>Diagnostic Tests Pass and Instrument in Ready Mode</li> <li>Breath tube heated</li> <li>Date, time and location code (Level 2,E). Re-set if necessary.</li> <li>Print test (Level 1,P). Sign and attach test record.</li> <li>Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.</li> </ol>
2. Interferent Check (Level 1,B): Known Value: O.Ook Acetone Sim. Ser #: DR7089 Lot #: ICS S Exp. Date: NA Display reads "Interferent Detect"  3. Interferent Detect"  4. Dry Calibration Check (Level 1,C): Known Value O.OOO AC Gas Cyl Lot #: 3491708043 Cyl #: Lexp. Date: 2.5-20 Test 1 0.000 Test 5 0.000 Test 7 0.000 Test 2 0.000 Test 5 0.000 Test 8 0.000 Test 3 0.000 Test 6 0.000 Test 9 0.000 Average O.OOO  C. Remarks/Maintenance: Annual Tuspecaron  Minstrument is acceptable to be used in the field.  Minstrument Signature  Date	В.	Configure simulator for the following tests (Level 1,S).  1. Wet Calibration Check  a. ✓ Low AC (Level 1,C): Known Value ≤ 0.03 AC: OCO AC Sim. Ser #: DR3375 Lot #: 17160 Exp. Date: OCO AC Sim. Ser #: DR7111 Lot #: 17350 Exp. Date: OCO AC Sim. Ser #: DR7111 Lot #: 17350 Exp. Date: OCO AC Sim. Ser #: DR7111 Lot #: 17350 Exp. Date: OCO AC Sim. Ser #: DR7111 Lot #: OCO AC AC Sim. Ser #: DR7111 Lot #: OCO AC
Sim. Ser #: <u>DR7081</u> Lot #: <u>TCS S Exp. Date: DR Display reads "Interferent Detect"</u> 3. PRFI Check (CMS Mode) Display reads "RFI Detect"  4. Dry Calibration Check (Level 1,C): Known Value O.080 AC Gas Cyl Lot #: <u>3411708043 Cyl #: LC Exp. Date: 25520 Test 1 0.080 Test 2 0.080 Test 5 0.080 Test 8 0.080 Test 8 0.080 Test 3 0.081 Test 9 0.081 Test 9 0.081 Average O.080  C. Remarks/Maintenance: ANNUAL TASPECTION  Winstrument is acceptable to be used in the field.  Date  NA</u>		V Company of the Comp
Instrument is acceptable to be used in the field.  Breath Analyst Signature  Date		Sim. Ser #: DR7089 Lot #: 2CS S Exp. Date: DR Display reads "Interferent Detect"  3. PRFI Check (CMS Mode) Display reads "RFI Detect"  4. Dry Calibration Check (Level 1,C): Known Value 0.080 AC Gas Cyl Lot #: 3491708043 Cyl #: Exp. Date: 2.5720  Test 1 0.081 Test 4 0.080 Test 7 0.080 Test 2 0.080 Test 5 0.080 Test 8 0.080 Test 3 0.081 Test 6 0.081 Test 9 0.081 Average 0.080
Breath Analyst Signature  NA  MA  4-//-2019  Date	C.	Remarks/Maintenance: <u>ANNUAL TUS PECTION</u>
		Analyst Signature Date
•	Revie	

02/18

Issuing Authority - CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006675
Location = BSPD 8164.14.00 09/16
04/11/2019 13:45

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34917080A3 Standard Cyl #? 16 Standard Expiration? 02/05/2020 Oper No? 666666

Flow Cal. Date: 08/25/2016 Slope 701 Intercept -701324

IR Calibration Date: 08/25/2016
3um 9um

Oth Coef(\*100): -8608 -21568

1st Coef(\*100): 258710 128604

2nd Coef(\*100): 2478 1279

H2O adj(mg/l\*10k): 381 528

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature CHARLES EDER

Remarks:

PRINT TEST

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = BSPD 04/11/2019 Alcohol Analyzer SN 80-006675 8164.14.00 09/16 14:00

# WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:01
02 Std. Sol.	0.019	14:02
03 Room Air	0.000	14:02
04 Std. Sol.	0.020	14:03
05 Room Air	0.000	14:04
06 Std. Sol.	0.020	14:04
07 Room Air	0.000	14:05

08 Sim Temp = 34.0°C

Simul Ser No = DR3375 Std Sol No = 17160 County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Low AC

0,020 AC

Form 106-I8000

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = BSPD 04/11/2019 Alcohol Analyzer SN 80-006675 8164.14.00 09/16 14:07

### WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:07
02 Std. Sol.	0.296	14:08
03 Room Air	0.000	14:09
04 Std. Sol.	0.297	14:09
05 Room Air	0.000	14:10
06 Std. Sol.	0.297	14:11
07 Room Air	0.000	14:11

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = DR7111 Std Sol No = 17350 County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

0.300 AC

HIGH AC

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006675 Location = BSPD 8164.14.00 09/16 04/11/2019 13:53

Test	AC	Time
01 Room Air	0.000	13:54
02 *Subject Test	INT*	13:54
03 Room Air	0.000	13:55

\*Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 01/01/1982

Sub Sex = Female Weight = 150 Test = DUI Cit = INTERFERENT

Dr. Lic. = ND/DIS821456

Lot No = 34917080A3

Cyl No = 16

Expiration Date = 02/05/2020

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Interferent Check 0.10AC + 0.05% ACETONE Form 106-I8000

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006675
Location = BSPD 8164.14.00 09/16
04/11/2019 13:56

Test	AC	Time
01 Diagnostic	OK	13:57
02 Room Air	0.000	13:58
03 *Subject Test	RFI*	13:58
04 Room Air	0.000	13:58

\*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 01/01/1982

Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS821456

Lot No = 34917080A3

Cyl No = 16

Expiration Date = 02/05/2020

County = 08

Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BSPD
04/11/2019

Alcohol Analyzer SN 80-006675 8164.14.00 09/16 13:47

#### DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:48
02 Std. Gas	0.081	13:48
03 Room Air	0.000	13:49
04 Std. Gas	0.080	13:49
05 Room Air	0.000	13:50
06 Std. Gas	0.081	13:50
07 Room Air	0.000	13:50

Lot No = 34917080A3

Cyl No = 16

Exp Date = 02/05/2020

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

0.080 AC

Form 106-I8000

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = BSPD 04/11/2019 Alcohol Analyzer SN 80-006675 8164.14.00 09/16 14:14

# DRY CAL CHECK

Te	est		AC	Time
01	Room	Air	0.000	14:14
02	Std.	Gas	0.080	14:15
03	Room	Air	0.000	14:15
04	Std.	Gas	0.080	14:16
05	Room	Air	0.000	14:16
06	Std.	Gas	0.081	14:17
07	Room	Air	0.000	14:17

Lot No = 34917080A3

Cyl No = 16

Exp Date = 02/05/2020

County = 08

Oper No. = 666666

Operator Signature
CHARLES EDER
CHECK

Remarks:

Form 106-I8000

D.080 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BSPD
04/11/2019

Alcohol Analyzer SN 80-006675 8164.14.00 09/16 14:17

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:18
02 Std. Gas	0.080	14:18
03 Room Air	0.000	14:19
04 Std. Gas	0.080	14:19
05 Room Air	0.000	14:20
06 Std. Gas	0.081	14:20
07 Room Air	0.000	14:20

Lot No = 34917080A3

Cyl No = 16

Exp Date = 02/05/2020

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

CALIBRATION CHECK

0.080 AC Form 106-18000