# OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

# **INTOXILYZER® 8000 INSPECTION**

Intoxil	lyzer <sup>®</sup> 8000 Serial Number: <u>80-006691</u>	Location: 78XL
A.	General Setup and Checks:  1. Diagnostic Tests Pass and Instrument in Recognition 2. Breath tube heated 3. Date, time and location code (Level 2,E). Recognition 2. Print test (Level 1,P). Sign and attach test in Tank monitor (Level 3,D,G). Display and Recognition 2. Display: 162 psi Regulator: 170 psi	e-set if necessary. record. egulator within 50 psi.
В.	Tests (Sign and attach test records):  Configure simulator for the following tests (Level 1):  Wet Calibration Check  a.	0.03 AC:O.20AC Exp. Date:G.8.19 2.0.25 AC:O.300_AC Exp. Date:/0.11.19
	Configure dry gas standard for the following tests  2. Interferent Check (Level 1,B): Known Value Sim. Ser #: DR 7089 Lot #: TCS 5 Exp. Display reads "Interferent Detect"  3. Interferent Check (CMS Mode)	e: 0.10 AC + 0.05% Acetone Date: \( \sum A \)  eads "RFI Detect" Value \( \frac{0.060}{0.060} \) AC  Exp. Date: \( \frac{2.5.20}{0.082} \)  Test 7 \( \frac{0.082}{0.082} \)
C.	Remarks/Maintenance: ANNUAL TUSPEC	non
$\mathcal{L}$	trument is acceptable to be used in the field.  Analyst Signature	4.//.2019 Date
Revie	wed by	Date

02/18

Issuing Authority - CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006691
Location = TOXL 8164.14.00 09/16
04/10/2019 16:24

\*\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34917080A3 Standard Cyl #? Standard Expiration? 02/05/2020 Oper No? 666666

Flow Cal. Date: 08/28/2018 Slope 699 Intercept -660610

IR Calibration Date: 08/28/2018

3um 9um

 0th Coef(\*100):
 -9426
 -20764

 1st Coef(\*100):
 258193
 134855

 2nd Coef(\*100):
 3532
 1465

 H2O adj(mg/l\*10k):
 408
 547

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature CHARLES EDER

Remarks: PRINT TEST

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/10/2019

Alcohol Analyzer SN 80-006691 8164.14.00 09/16 16:26

# WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:26
02 Std. Sol	0.020	16:27
03 Room Air	0.000	16:28
04 Std. Sol	0.020	16:28
05 Room Air	0.000	16:29
06 Std. Sol	0.020	16:30
07 Room Air	0.000	16:30

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = DR3375 Std Sol No = 17160

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Low AC

Remarks:

0,020 AC Form 106-18000

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/10/2019

Alcohol Analyzer SN 80-006691 8164.14.00 09/16 16:31

## WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:32
02 Std. Sol.	0.300	16:33
03 Room Air	0.000	16:33
04 Std. Sol.	0.300	16:34
05 Room Air	0.000	16:34
06 Std. Sol.	0.301	16:35
07 Room Air	0.000	16:36

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = DR7111 Std Sol No = 17350

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

0.300 AC

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006691 Location = TOXL 8164.14.00 09/16 04/10/2019

Test			AC	Time
01	Room Air		0.000	16:57
02	*Subject '	Test	INT*	16:57
03	Room Air		0.000	16:58

16:56

#### \*Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 01/01/1982

Sub Sex = Female Weight = 150

Test = DUI Cit = INTERFERENT

Dr. Lic. = ND/DIS821456

Lot No = 34917080A3

Cyl No = 28

Remarks:

Expiration Date = 02/05/2020

County = 08Oper No. = 666666

Operator Signature

Form 106-I8000

CHARLES EDER

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006691
Location = TOXL 8164.14.00 09/16
04/10/2019 16:58

Test	AC	Time	
01 Diagnostic	OK	16:59	
02 Room Air	0.000	17:00	
03 *Subject Test	RFI*	17:00	
04 Room Air	0.000	17:00	

\*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 01/01/1982

Sub Sex = Female Weight = 150 Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS821456

Lot No = 34917080A3

Cyl No = 28

Expiration Date = 02/05/2020

County = 08

Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/10/2019

Alcohol Analyzer SN 80-006691 8164.14.00 09/16 17:01

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:01
02 Std. Gas	0.082	17:02
03 Room Air	0.000	17:02
04 Std. Gas	0.082	17:02
05 Room Air	0.000	17:03
06 Std. Gas	0.082	17:03
07 Room Air	0.000	17:04

Lot No = 34917080A3

Cyl No = 28

Exp Date = 02/05/2020

County = 08

Oper No. = 666666

Operator Signature
CHARLES EDER

CACIBRATION CHECK

Remarks:

0.080AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/11/2019

Alcohol Analyzer SN 80-006691 8164.14.00 09/16 10:08

## DRY CAL CHECK

Test			AC		Time	
	01	Room	Air	(	0.000	10:09
	02	Std.	Gas	(	0.082	10:09
	03	Room	Air	(	0.000	10:10
	04	Std.	Gas	(	0.082	10:10
	05	Room	Air	(	0.000	10:11
	06	Std.	Gas	(	0.082	10:11
	07	Room	Air	(	0.000	10:11

Lot No = 34917080A3

Cyl No = 28

Exp Date = 02/05/2020

County = 08

Oper No. = 666666

Operator Signature
CHARLES EDER

LAUBRATION CHECK

Remarks:

0.080 AC Form 106-18000

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/11/2019

Alcohol Analyzer SN 80-006691 8164.14.00 09/16 10:12

### DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	10:12
02	Std.	Gas	0.082	10:12
03	Room	Air	0.000	10:13
04	Std.	Gas	0.082	10:13
05	Room	Air	0.000	10:14
06	Std.	Gas	0.082	10:14
07	Room	Air	0.000	10:15

Lot No = 34917080A3

Cyl No = 28

Exp Date = 02/05/2020

County = 08

Oper No. = 666666

Operator Signature

PACIBRATION CHECK

Remarks:

Form 106-I8000

0.080 AC