Intoxilyzer® 8000 Inspection

## **INTOXILYZER® 8000 INSPECTION**

Intoxilyzer<sup>®</sup> 8000 Serial Number: <u>80-006687</u>

Location: ToxL

- A. General Setup and Checks:
  - 1. Diagnostic Tests Pass and Instrument in Ready Mode
  - 2. Breath tube heated
  - 3. A Date, time and location code (Level 2,E). Re-set if necessary.
  - 4. XPrint test (Level 1,P). Sign and attach test record.
  - 5. XTank monitor (Level 3,D,G). Display and Regulator within 50 psi.
    - a. Display: <u>789</u> psi Regulator: <u>775</u> psi
- B. Tests (Sign and attach test records):
  - $\square$  Configure simulator for the following tests (Level 1,S).
  - 1. Wet Calibration Check
    - a. X Low AC (Level 1,C): Known Value ≤ 0.03 AC: \_\_\_\_\_AC Sim. Ser #: <u>DR3375</u> Lot #: <u>/7/66</u> \_\_ Exp. Date: <u>6.8.19</u>
    - b.
       Image: AC
       Imag

Configure dry gas standard for the following tests (Level 1,S).

- Interferent Check (Level 1,B): Known Value: O·IO AC + 0.05% Acetone Sim. Ser #: OR 7089 Lot #: ICS 5 Exp. Date: NA
  ▲ Display reads "Interferent Detect"
- 3. ARFI Check (CMS Mode) Display reads "RFI Detect"
- 4. Dry Calibration Check (Level 1,C): Known Value <u>0.080</u> AC Gas Cyl Lot #:<u>34419080AR</u> Cyl #: <u>45</u> Exp. Date: <u>2.5.21</u> Test 1 <u>0.082</u> Test 4 <u>0.081</u> Test 7 <u>0.081</u> Test 2 <u>0.082</u> Test 5 <u>0.081</u> Test 8 <u>0.082</u> Test 3 <u>0.082</u> Test 6 <u>0.082</u> Test 9 <u>0.082</u> Average <u>0.082</u>
- C. Remarks/Maintenance: <u>ANNUAL INSPECTION</u>

Instrument is acceptable to be used in the field.

colle.

Breath Analyst Signature

Reviewed by

.17.2019 JA

Date

Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006687 Location = TOXL8164.14.00 09/16 04/17/2019 09:38 \*\*\*\*\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*\*\*\*\* abcdefqhijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34418080A2 Standard Cyl #? 45 Standard Expiration? 02/05/2021 Oper No? 666666 Flow Cal. Date: 10/05/2018 Slope 696 Intercept -609977 IR Calibration Date: 10/05/2018 3um 9um \_\_\_\_\_ Oth Coef(\*100):, -11478 1st Coef(\*100): 255256 2nd Coef(\*100): 4432 -18771132135 1676 H2O adj(mg/l\*10k): 604575 

Operator Signature CHARLES EDER

Remarks:

INT TES

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006687Location = TOXL8164.14.00 09/1604/17/201909:39

|    |           | WET CAL CHECK |       |
|----|-----------|---------------|-------|
|    | İ         |               |       |
| Т  | est       | AC            | Time  |
| 01 | Room Air  | 0.000         | 09:39 |
| 02 | Std. Sol. | 0.020         | 09:40 |
| 03 | Room Air  | 0.000         | 09:40 |
| 04 | Std. Sol. | 0.020         | 09:41 |
| 05 | Room Air  | 0.000         | 09:42 |
| 06 | Std. Sol. | 0.020         | 09:42 |
| 07 | Room Air  | 0.000         | 09:43 |
| 00 | Cim Elomo | 24 089        |       |

08 Sim Temp = 34.0°C

Simul Ser No = DR3375 Std Sol No = 17160 County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Low A . eC O. C

Form 106-18000

0.020 AC

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006687Location = TOXL8164.14.00 09/1604/17/201909:44

|   | WET CAL CHECK                                      |  |
|---|--|--|
| Test  | AC   | Time   |
| 01 Room Air<br>02 Std. Sol.<br>03 Room Air<br>04 Std. Sol.<br>05 Room Air<br>06 Std. Sol. | 0.000<br>0.294<br>0.000<br>0.294<br>0.000<br>0.295 | 09:44<br>09:45<br>09:45<br>09:46<br>09:47<br>09:47 |
| 07 Room Air   | 0.295  | 09:48  |

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = DR7111 Std Sol No = 17350 County = 08

Oper No. = 666666

States

Operator Signature CHARLES EDER

HIGH AC 0.300 AC

Form 106-I8000

NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006687 Location = TOXL8164.14.00 09/16 04/17/2019 09:51 Test AC Time 01 Room Air 0.000 09:52 02 \*Subject Test INT\* 09:52 03 Room Air 0.000 09:53 \*Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 01/01/1982Sub Sex = Female Weight = 150Test = DUI Cit = INTERFERENT Dr. Lic. = ND/DIS821456Lot No = 34418080A2Cyl No = 45Expiration Date = 02/05/2021County = 08Oper No. = 666666

Intoxilyzer Test Record and Checklist

Operator Signature CHARLES EDER

.as%. ACETONE er er

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006687 Location = TOXL8164.14.00 09/16 04/17/2019 09:53 Test AC Time 01 Diagnostic 09:54 OK 02 Room Air 0.000 09:55 03 \*Subject Test RFI\* 09:55 04 Room Air 0.000 09:55 \*Invalid Test Inhibited - RFI Sub Name = DISCOVER, THE SPIRIT Sub DOB = 01/01/1982Sub Sex = Female Weight = 150Test = DUI Cit = RFI CHECK Dr. Lic. = ND/DIS821456Lot No = 34418080A2Cyl No = 45Expiration Date = 02/05/2021County = 08Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006687Location = TOXL8164.14.00 09/1604/17/201909:55

| г           |               |       |
|-------------|---------------|-------|
|             | DRY CAL CHECK |       |
| L           |               |       |
| Test        | AC            | Time  |
| 01 Room Air | 0.000         | 09:56 |
| 02 Std. Gas | 0.082         | 09:56 |
| 03 Room Air | 0.000         | 09:57 |
| 04 Std. Gas | 0.082         | 09:57 |
| 05 Room Air | 0.000         | 09:58 |
| 06 Std. Gas | 0.082         | 09:58 |
| 07 Room Air | 0.000         | 09:59 |
|             |               |       |

Lot No = 34418080A2 Cyl No = 45 Exp Date = 02/05/2021 County = 08

Oper No. = 666666

Contine C Ċ

Operator Signature CHARLES EDER

IECE 7*0*a 0.09

Form 106-18000

| CMI, Inc. Intoxilyzer   | Alcohol Analyzer |
|-------------------------|------------------|
| North Dakota Model 8000 | SN 80-006687     |
| Location = TOXL         | 8164.14.00 09/16 |
| 04/17/2019              | 09:59            |

|    | Ņ    |     | DRY CAL CHECK |       |
|----|------|-----|---------------|-------|
| Te | est  |     | AC            | Time  |
| 01 | Room | Air | 0.000         | 09:59 |
| 02 | Std. | Gas | 0.081         | 10:00 |
| 03 | Room | Air | 0.000         | 10:00 |
| 04 | Std. | Gas | 0.081         | 10:01 |
| 05 | Room | Air | 0.000         | 10:01 |
| 06 | Std. | Gas | 0.082         | 10:01 |
| 07 | Room | Air | 0.000         | 10:02 |

Operator Signature CHARLES EDER

IBRATION CHECK 0.080 AC

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006687Location = TOXL8164.14.00 09/1604/17/201910:02

|             | DRY CAL CHECK |       |
|-------------|---------------|-------|
| Test        | AC            | Time  |
| 01 Room Air | 0.000         | 10:03 |
| 02 Std. Gas | 0.081         | 10:03 |
| 03 Room Air | 0.000         | 10:04 |
| 04 Std. Gas | 0.082         | 10:04 |
| 05 Room Air | 0.000         | 10:04 |
| 06 Std. Gas | 0.082         | 10:05 |
| 07 Room Air | 0.000         | 10:05 |

Lot No = 34418080A2 Cyl No = 45 Exp Date = 02/05/2021 County = 08

Oper No. = 666666

. CCCC

Operator Signature CHARLES EDER

Remarks:

A CIBRATION CHECK 0,080 AC

Form 106-18000